ceraMotion®
Me

Clické de : © Christian Ferrario

Gebrauchsanweisung | Instructions for use | Mode d'emploi |
Modo de empleo | Modalità d'uso
CAD/CAM Primer

D
DENTAURUM

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Classification

ceraMotion® Me primer CAD/CAM belongs to the group of veneering ceramics Class 1 according to EN ISO 6872:2015. When layering with ceraMotion® Me, the primer is used as the first layer on a framework made with a dental alloy instead of a first layer of opaque.

Indication

ceraMotion® Me primer CAD/CAM is suitable for all dental alloys with a coefficient of thermal expansion (CTE) from $13.9 - 15,1 \cdot 10^{-6} \text{ K}^{-1}$ (25–500°C / 77 – 932 °F), in particular for CAD/CAM generated frameworks made of non-precious alloys (milling, laser melting ...).

ceraMotion® Me primer CAD/CAM may not be used on frameworks with a CTE outside the specified values.

The product may not be used if there is an intolerance towards specific contents.

Primer CAD/CAM REF 252-500-03

Use

Framework preparation:

- Manufacture framework, avoiding sharp edges and keeping to an anatomical tooth shape, reduced in size. The ceramic layering should not exceed 2 mm maximum.
- Adhere to the manufacturer's instructions for finishing, blasting and firing.
- Blast the remanium® alloys from Dentaureum with Al_2O_3 (110 μm) at 2 bar; clean the framework using a steam cleaner and dry.
- No oxide firing is required. Should oxide firing take place to ensure the alloys have been correctly processed, blast the oxide layer, then clean the framework.
- Do not touch the alloy afterwards.

Primer application:

- Apply the amount of paste needed onto a glass plate.
- Mix the paste before use in the pot; the paste should have a creamy consistency.
- If necessary, use a small amount of ceraMotion® Paste Liquid (REF 254-006-02) to adapt the consistency. Avoid adding liquid in the pot.

- Use a clean and dry brush. Clean the brush with ceraMotion® Paste Liquid (REF 254-006-02) or steam, then dry it.
- Spread the creamy primer evenly. Apply the primer similarly to a layer of opaque in appropriate thickness.
- Firing according to Firing Chart 1.
- The ceraMotion® Me primer has a shiny appearance after firing.

Opaque application:

- When the framework has been given a layer of primer, apply a layer of opaque.
- Firing according to Firing Chart 1.
- One layer of opaque is sufficient to cover the framework and to produce the shade required before the ceramic is applied.

Table 1 – Firing Chart:

	Temperature at start (°C / °F)	Drying time (min.)	Heating rate (°C / °F per min.)	Vacuum start (°C / °F)	Vacuum end (°C / °F)	Final temperature (°C / °F)	Holding time (min.)
ceraMotion® Me primer	500 / 932	8	75 / 167	500 / 932	980 / 1796	980 / 1796	1 min. (under vacuum)
ceraMotion® Me Paste Opaque	500 / 932	8	75 / 167	500 / 932	950 / 1742	950 / 1742	1 min. (under vacuum)

Table 2 – Technical data:

Material	
Chemical composition	Glass ceramic in paste, mainly SiO_2 , Al_2O_3 , K_2O , Na_2O , B_2O_3 , ZrO_2 , Polyvalent alcohol (glycol), silica sol and additives

Classification in accordance with DIN EN ISO 6872: Type II Class 3		
Bonding strength (Schwickerath test)	ISO 9693-1:2012 (> 25 MPa)	> 47 MPa
Coefficient of thermal expansion (25 - 500 °C / 77 – 932 °F) [10^{-6} K^{-1}]	ISO 6872:2015	2 firings ≈ 12.1 4 firings ≈ 12.0
Mechanical strength (3-point bending test)	ISO 6872:2015 (> 50 MPa)	$\approx 120 \text{ MPa}$
Chemical solubility	ISO 6872:2015 (< 100 $\mu\text{g}/\text{cm}^2$)	$\approx 20 \mu\text{g} / \text{cm}^2$
Glass transition temperature	ISO 6872:2015	$\approx 610 \text{ °C} / 1130 \text{ °F}$

Warning notices

Only for dental use by a qualified expert. The user must ensure that the intended use corresponds to the manufacturer's specifications.

Measures for protection of equipment and products

Ensure the product is handled correctly and that furnaces are maintained and calibrated regularly. Wear safety glasses and breathing protection and use a suction unit when blasting and grinding the objects. There is a risk of suffering burns when removing objects from the furnace.

General information

To protect the product against contamination, choose a suitable work area and use clean instruments and equipment. In the processing of the medical products, products in question should be stored between 10°C and 30°C (50°F and 95°F).

These non-toxic products can be disposed of without an issue if the statutory provisions are upheld.

Symbols used:

Refer to Instructions for use

Batch number (LOT)

Storable until (Year-Month)

CE mark

RxOnly

Warning: Sale and use only to/by dental technicians.

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