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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use.

3. Area of application

For optimal gingiva management and for transferring the oral situation to the working cast.

4. Medical indications/contraindications

Temporary abutment tiologic® ST

Indications

- Temporary restorations for up to 180 days in situ
- Contraindications.
- In situ longer than 180 days
- Primary bracing of abutments
- Single restoration with cantilever bridge unit
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25

5. Information on use, processing, fitting and insertion

5.1 Closure screws tiologic® ST

After implant insertion and monitoring of the operation site, the wound can be closed by suturing. Ensure that the wound closure is saliva-proof and that there is good blood circulation.

Prior to wound closure, retrieve the sterile closure screw included in the implant packaging by pressing the hex key 1.3 into the head of the screw and pulling the screw from the implant holder. The closure screw is then inserted into the implant and tightened. Closure screws should fit flush on the implant to ensure that bone tissue growth cannot penetrate into the implant. The closure screws are inscribed with S, M or L corresponding to the series of abutments.

Closure screws are intended for single use only.

Tightening torque

- Closure screw: manually or 15 Ncm
- Should a closure screw lose its sterility during the operation, it is possible to purchase closure screws individually as a non-sterile product.

IMPORTANT: Closure screws purchased separately should be sterilized before they are inserted into the implant.

5.2 Gingiva former tiologic® ST

Gingiva formers, conical or cylindrical, or bar, bridge, AngleFix abutments – particularly gentle on the tissue – are available for the implantologist to ensure optimal management of the gingiva. Gingiva formers are available for the S, M and L series of abutments.

The gingiva formers are designed to form a wide gingival contour. Depending on the type of prosthetic restoration, this can make it easier for the implantologist to fit the restoration.

The gingiva formers are selected according to the series of abutments, gingival height and insertion depth of the implant. They are available for the series of abutments S, M or L and in different gingival heights (laser-printed on gingiva former).

Tightening torque

Gingiva formers can also be used with open healing of the implant for specific indications and for preserving the soft tissue.

IMPORTANT: The gingiva formers and bridge, bar and AngleFix abutments should be cleaned and sterilized before insertion in the implant.

If a temporary restoration is fitted, the denture should be relieved during gingiva forming. The impression should not be taken until the tissue is completely free of inflammation.

Tightening torque

- Gingiva former: manually or 15 Ncm
- Bar, bridge, AngleFix abutment: 35 Ncm
- Bar, bridge, AngleFix closure screw: 15 Ncm

5.3 Impression post tiologic® ST

The impression can be taken using either the open or closed technique. Components for the series of abutments S, M and L are available for both techniques.

IMPORTANT: When positioning the impression post on the inserted implant, ensure it is seated without a gap. The impression component should NOT come into contact with the individual impression tray.

In the case of removable restorations (restorations with bar, bridge or Anglefix, ball abutments, tioloc) the impression can also be taken with other special impression components over the respective primary abutments.

Silicone or polyether impression materials are recommended for impression-taking due to their high precision and elastic recovery.

Open impression technique:

Impression posts are available in lengths of 10.0 mm and 14.0 mm with the corresponding screws to cater for different occlusal spaces. The impression posts are laser-printed with S, M or L on the retention surface and at the interface.

In order to make the abutment series more recognizable, the interface is additionally marked with dots:

- 1 dot is equivalent to abutment S
- 2 dots are equivalent to abutment M
- 3 dots are equivalent to abutment L

The screw is pushed down before fitting the impression post. This provides additional guidance when fitting the post. The inner connection is shorter with an open impression post to ensure a compression-free impression even with divergent axes. The enclosed red screw for open impression-taking has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection.

Tightening torque

- Screw impression post intra-orally: manually or 15 Ncm
- Screw impression post in the laboratory implant: manually or 15 Ncm

Closed impression technique:

Components for the closed impression technique include impression posts, retaining screws and impression caps. They are laser-printed or marked with the series of abutments S, M or L. The enclosed red retaining screw has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection. The corresponding impression aid S, M or L is placed taking the vertical retention grooves into consideration until it perceptibly and audibly clicks into place.

IMPORTANT: The impression caps are items for single use only and are NOT suitable for sterilization or multiple use.

The orientation of the retention grooves ensures that they can be positioned without coming into contact with the adjacent teeth. The impression is taken according to standard criteria.

Tightening torque

- Retaining screw, intra-orally: manually, or 15 Ncm
- Retaining screw on laboratory implant: manually, or 15 Ncm

5.4 Temporary abutment tiologic® ST

In aesthetically relevant areas, the peri-implant structures are retained with a temporary abutment.

It is possible to fit a long-term, non-functional immediate temporary restoration on implants if there is absolute primary stability and no recession of the implant site. The temporary abutment is available for the S, M and L series of abutments.

It is supplied non-sterile and is made from high-strength plastic (PEEK), which can be quickly and easily customized. The temporary abutment can be faced directly with composite or fitted with a temporary crown or bridge. In both cases the abutment is secured intra-orally with the screw for the temporary abutment; the contours are marked and adjusted extra-orally. The implant should **NOT** be subjected to mechanical stress.

The restoration should only be shortened as far as the upper edge of the screw for the temporary abutment.

Tightening torque

- Screw for temporary abutment on the model: manually
- Screw for temporary abutment in the mouth: 15 Ncm

6. Composition

- Closure screws: titanium alloy (Ti6Al4V)
- Gingiva former: titanium alloy (Ti6Al4V)
- Impression posts:
 - open – impression post and screw: titanium alloy (Ti6Al4V)
 - closed – impression post and screw: titanium alloy (Ti6Al4V)/impression cap: POM
- Temporary abutments: PEEK (polyether ether ketone)/ Screw: titanium alloy (Ti6Al4V)
- Prosthetic screws:
 - AnoTite screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw: titanium alloy (Ti6Al4V)
- Hex key:
 - Hex key - ISO shank 1.3: (1.4197)
 - Hex key - ratchet 1.3: (1.4197)

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation.All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

8. Warnings

Temporary restorations made of PEEK can be left in situ for a maximum of 180 days.

9. Storage and shelf life

Store in a dry place.

10. Information for single use products

Closure screws, gingiva formers, temporary abutments, impression posts and impression caps (closed impression-taking) are designed for single use only. Reconditioning of a closure screw, gingiva former, temporary abutment, impression post or impression cap (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

11. Information for multiple use products

The hex key 1.3 is intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Information on disinfection/cleaning/sterilizing

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

13. Product overview /scope of delivery, incl. REF

Please refer to the tiologic® ST product catalog for the tiologic® ST product range.

14. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

15. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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