DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurem.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurem.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurem GmbH & Co. KG
I Turnstr. 31
I 75228 Ispringen
I Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment. Selection aids for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

The titanium abutments tiologic® ST are prosthetic components for single tooth and bridge restorations on tiologic® ST implants.

4. Medical indications/contraindications

Titanium abutments tiologic® ST

Indications

- Individual crown and bridgework (cemented, individually screw-retained)
- Abutment on implant ø 3.3 mm only for use in the area of the upper lateral incisors and lower incisors

Contraindications

- Primary bracing of abutments
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants ø 3.3 mm
- Cast-on technique

5. Information on use, processing, fitting and insertion

Titanium abutment tiologic® ST

All titanium abutments tiologic® ST are available in conical and platform connector geometry for the abutment series S, M and L.

The straight and angled (S 15°, M 20°, L 20°) titanium abutments can be slightly altered to adapt to the gingival line (straight 1.0 mm, 2.5 mm, 4.0 mm/ angled 1.5 mm, 3.0 mm). These abutments have an exactly defined crown margin and an integrated rotational stop.

The angled titanium abutments also have different gingival heights (labial/ palatal). The universal titanium abutments are cylindrical or anatomical and can be prepared as required.

Tightening torque

- Prosthetic screw, on model: manually
- AnoTite screw intra-orally: 30 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- Titanium abutments: titanium alloy (Ti6Al4V)

- Prosthetic screws:
 - AnoTite screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw: titanium alloy (Ti6Al4V)

- Hex key:
 - Hex key - ISO shank 1.3: (1.4197)
 - Hex key - ratchet 1.3: (1.4197)

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurem.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

- Different types of alloy in the oral cavity can lead to galvanic reactions.

- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Titanium abutments are designed for single use only. Reconditioning of a titanium abutment (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex key 1.3 is intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurem.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurem generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurem.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview /scope of delivery, incl. REF

Please refer to the tiologic® ST product catalog for the tiologic® ST product range.

13. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

14. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurem.com (Explanation of symbols REF 989-313-00).

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