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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG
I Turnstr. 31
I 75228 Ispringen
I Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

AngleFix abutments tiologic® ST are prosthetic components for bridge and bar restorations with tiologic® ST implants. They can compensate implant divergences of up to 50°, thereby enabling bridge and bar restorations to have a passive fit.

4. Medical indications / contraindications

AngleFix abutments tiologic® ST

Indications

- Anchorage of implant-supported, splinted bridges and bars in the maxilla in combination with at least 6 tiologic® ST implants and, in the mandible, in combination with at least 4 or more tiologic® ST implants.

Contraindications.

- For biomechanical reasons we recommend that the following angulations are not exceeded:

AngleFix abutment with 0°:	0° - 12°
AngleFix abutment with 18°:	6° - 30°
AngleFix abutment with 32°:	20° - 44°
- Restorations with mixed retention (tooth/AngleFix)
- Combination with other anchorage elements
- Grinding and shortening AngleFix abutments

5. Information on use, processing, fitting and insertion

5.1 AngleFix abutments tiologic® ST

AngleFix abutments are available for the S, M and L series of abutments. AngleFix abutments are provided in 3 angulations: 0°, 18° and 32°. The seating surface of the AngleFix caps (ø 5.3 mm) and the cone (12°) of AngleFix abutments are always identical for the S, M and L series of abutments, so that only one size of impression posts, closure caps etc. is required. The gingival height of the straight AngleFix abutments is 1.0 mm, and 2.5 mm for the angulated AngleFix abutments.

In case of steep inclination of the implants, splinting the denture is absolutely essential.

When fitting the angulated AngleFix abutments, keep in mind that the abutments could have a larger diameter than the implants. As part of the abutments may be below the bone line, bone may have to be removed in this region, if required, so that the abutment sits flush on the implant. Prerequisite for the successful use of the AngleFix system is the best possible accurate angular position of the implant. The more accurate this angle can be maintained, the easier the prosthetic treatment, as the crowns of the abutments are then positioned parallel to one another.

The AngleFix abutments should **NOT** be ground or shortened.

When handling the AngleFix abutments, ensure that the occlusal screw is not damaged by tweezers etc.

Every construction must be checked for passive fit prior to placing it intra-orally (Sheffield test). This involves placing the construction on the AngleFix abutment and only one prosthetic screw (not included in delivery) is fixed in the AngleFix abutment. If the construction is raised and a gap is created between the construction and AngleFix abutment, stresses are present and must be eliminated.

5.2 AngleFix impression post, open and closed tiologic® ST

To achieve stable apposition of the gingiva with AngleFix abutments, it is possible to take an open or closed impression on a AngleFix abutment that is fitted in the mouth. Since the AngleFix impression posts have no rotational security, they may only be used for bridge and bar constructions.

5.3 AngleFix closure screw tiologic® ST

The AngleFix closure screw is screwed directly onto the AngleFix abutment after impression-taking and protects the surface for the caps on the AngleFix abutments.

5.4 AngleFix plastic cap tiologic® ST

The plastic cap, that burns residue-free, can be individually adapted to the space available. The processing instructions given by the investment material manufacturer must be adhered to for investing and preheating. After casting, the casting ring must be cooled to room temperature to allow devesting and to achieve stress-free cast objects. The cast object should never be blasted with abrasive agents since these can damage the inner configuration. Devesting is carried out either by steam cleaning or in an ultrasonic bath.

5.5 AngleFix titanium cap tiologic® ST

The AngleFix titanium cap is available for the passive fit technique (adhesive technique). Once a cast bar has been produced, it is bonded intra-orally to the titanium caps. This helps compensate possible inaccuracies of fit that could create tensions.

5.6 AngleFix scan cap titanium tiologic® ST

The surface of the AngleFix scan cap titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. For this purpose, the AngleFix scan cap is screwed onto the AngleFix abutment and scanned. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer.

5.7 AngleFix laboratory implant tiologic® ST for printed and cast models

The laboratory implant can be used in both analog and digital workflow. The cavities of the laboratory analogs are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory analogs with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary. In order for the counter screw supplied to fix the laboratory analog at the apical end in the printed model, the model base should be printed with a sufficient height. If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory analog in the cervical area of the printed model. The 3 retention tabs should still be visible in the printed model for monitoring purposes.

When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended series of abutments S, M or L. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The closure screw supplied with the laboratory implant enables the laboratory implant to be fixed in the end position from apical. The straight AngleFix abutments (0°) are screwed into the tiologic® ST implant with an insertion key for bar/bridge/AngleFix abutments. The angulated AngleFix abutments (18° and 32°) have a rotationally secure fit in the tiologic® ST implants and are fixed with the AnoTite screw 9.0 mm. A new AnoTite screw should be used for the final fitting in the mouth.

To fix bridge and bar constructions on AngleFix abutments, special prosthetic screws and AnoTite screws are available.

The special AnoTite screws with a reduced screw head height can be used for limited occlusal spaces. Such screws are not included in delivery.

Screws with hexagon screw head geometry are tightened manually with the hex key ratchet 1.3 or mechanically with the hex key ISO shank 1.3.

Tightening torque

- AngleFix abutments, straight on model: by hand
- AngleFix abutments, straight intra-orally: 35 Ncm
- Prosthetic screw on model: manually
- AnoTite screw intra-orally: 30 Ncm
- Prosthetic screw for AngleFix caps on model: by hand
- AnoTite screw for AngleFix caps intra-orally: 25 Ncm
- Sure-grip screw impression post on AngleFix abutment intra-orally: 15 Ncm
- Closure screw on AngleFix abutment intra-orally: 15 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw.

The prosthetic screw serves only to help the dental technician fix the abutment on the working cast. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- AngleFix abutments: titanium alloy (Ti6Al4V)
- AngleFix caps:
 - AngleFix plastic cap: polycarbonate
 - AngleFix titanium cap: titanium alloy (Ti6Al4V)
 - Closure screw: titanium alloy (Ti6Al4V)
 - AngleFix impression post, open: titanium alloy (Ti6Al4V) / screw: 1.4305
 - AngleFix impression post, closed:
 - impression post and screw: titanium alloy (Ti6Al4V) / impression cap: POM
 - AngleFix scan cap titanium: titanium alloy (Ti6Al4V)
- AngleFix laboratory implant: titanium alloy (Ti6Al4V)
- Prosthetic screws:
 - AnoTite screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw: titanium alloy (Ti6Al4V)
 - Prosthetic screw for AngleFix caps: titanium alloy (Ti6Al4V)
 - AnoTite screw for AngleFix caps: titanium alloy (Ti6Al4V)
- Insertion key bar/bridge/AngleFix ratchet: 1.4197
- Hex keys: - Hexagon key - ISO shank 1.3: 1.4197
 - Hexagon key - ratchet 1.3: 1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.
- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

AngleFix abutments and caps are designed for single use only. Reconditioning of a AngleFix abutment or cap (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex key 1.3, the insertion key bar/bridge/AngleFix and the AngleFix scan caps titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview /scope of delivery, incl. REF

Please refer to the tiologic® ST product catalog for the tiologic® ST product range.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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