---

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

## 2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function and thus improve the quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. Prosthetic abutments and components serve to produce prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components						
Series of abutments*	S		M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5	

\*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments.

The ball / tioLOC abutments are only available with a platform connector geometry.

tioLogic® ST ball and tioLOC abutments are prefabricated components for fabricating implant-retained tissue-borne prosthetic restorations. The ball and tioLOC abutments are available in different gingival heights to suit the varying gingival conditions of patients. The ball and tioLOC abutments are marked according to the S, M or L series of abutments.

Ball abutment			
Series of abutments	S	M	L
GH in mm	1.5 / 3.0 / 4.5		
Ball diameter in mm	2.25		

tioLOC abutment			
Series of abutments	S	M	L
GH in mm	1.0 / 2.0 / 3.0 / 4.0 / 5.0		

Ball abutment laboratory implants are used to create working casts. They accurately reflect the position of the tioLogic® ST ball abutment. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

The tioLogic® ST implant system contains specifically coordinated abutments and components, including instruments, for the insertion and prosthetic restoration of the implants.

Temporary abutment			
Series of abutments	S	M	L
Length in mm	12.0		

Screws for tioLogic® ST temporary abutments are used to secure temporary abutments for time-limited long-term temporary, non-functional immediate restoration on (laboratory) implants. They are supplied with the respective temporary abutment.

## 3. Intended purpose

Abutments are intended to connect the implants with the prosthetic restoration.

Ball abutment laboratory implants are designed to represent the position of the ball abutment (axial inclination and height) in a working cast.

## 4. Intended users

The tioLogic® ST product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

## 5. Medical indications / contraindications

### 5.1 Indications

#### 5.1.1 tioLogic® ST ball abutments

- Resilient anchoring of implant-retained, tissue-borne coverdentures on two implants to ensure the tangential axis of rotation
- Anchoring of implant-retained full dentures on at least two implants
- Incorporation in an existing coverdenture as an interim denture or with partial denture framework possible
- Prosthetic restoration on at least two implants

#### 5.1.2 Ball abutment laboratory implants

- precise reproduction of positions for cast models (dental stone)

#### 5.1.3 tioLogic® ST tioLOC abutments

- Anchoring of implant-retained, tissue-borne coverdentures with very low profile height
- Incorporation into an existing coverdenture is possible
- Prosthetic restoration on two or four implant posts

### 5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

#### 5.2.1 tioLogic® ST ball abutments

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- Combination with other retention elements
- Divergence of more than 20° between the implants
- Heavy mechanical strain with  $\varnothing$  3.3 mm implants

#### 5.2.2 Ball abutment laboratory implants

- Use of laboratory implants with defective or damaged connector geometry
- Use on patients

#### 5.2.3 tioLogic® ST tioLOC abutments

- Uneven number of implants per jaw
- Unfavorable positioning of implants which prevent a tangential rotation axis
- Restorations with mixed retention (tooth/implant)
- Combination with other retention elements
- Implant divergences of more than 10° using tioLogic® ST implants  $\varnothing$  3.3 mm are contraindicated
- Intermaxillary space should be at least 4.0 mm

## 6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

## 7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® ST Instructions for use.

We strongly recommend that, prior to using Dentaurum products, dental professionals attend a Dentaurum training course to familiarize themselves with the Dentaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurum recommends full clinical, radiological, photographic and statistical documentation.

The ball and tioLOC abutments and components of the tioLogic® ST implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.


**Precautions:**

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

**Important:**

**The implantologist should ensure the products cannot be aspirated during intra-oral use.**

**7.1 tioLogic® ST ball abutments**


Ball abutment

Ball anchors function in such a manner that it is advisable to support the prosthetic restoration on at least two implants. The implants must not diverge by more than 20°. Care should be taken to avoid an excessive mechanical loading when using ball head abutments together with tioLogic® ST  $\varnothing$  3.3 mm implants. The tioLogic® ST ball abutments are available with platform connector geometry in

gingival heights of 1.5 mm, 3.0 mm and 4.5 mm for S, M and L abutments. The gingival height is the distance between the uppermost edge of the implant and the lowermost edge of the cylindrical area of the ball abutments. This upper edge should be approx. 1.0 mm above the gingival line. To achieve optimum retention, the ball abutments should each be positioned at the same level. The ball head is 2.25 mm in diameter. No components of the ball abutment may be modified. Only the matrix withdrawal force can be adjusted by activating the inner matrix. The ball abutment is placed with the ball abutment insertion key.

**Working procedures**

- In the direct technique, process the matrix into an existing denture directly in the patient's mouth without fabricating a model. Neither a tioLogic® ST impression post nor laboratory implant are required. The tioLogic® ST ball abutments are attached to the tioLogic® ST implants intra-orally. During the clinical procedure, protect the undercuts with the help of a rubber dam. It is put over the matrix to prevent cold-curing resin filling the undercuts in the mouth which would prevent the denture being removed after the resin has cured. Observe the instructions for use of the resin.
- In the indirect technique, take an impression and fabricate a model. The matrix is processed into the denture in the laboratory.

Ball abutment laboratory implant

With the indirect version, the impression can also be taken using the closed technique directly over the ball abutments fixed in the mouth. To do this, select the correct tioLogic® ST ball abutments for the gingival height and S, M or L series of abutment, remove the tioLogic® ST closure screws or tioLogic® ST gingiva formers and fix the tioLogic® ST ball abutment in the implant. Apply impression material

around all ball head abutments, fill the closed impression tray with impression material and take the impression. Remove the tray after the impression material is cured. Ensure that the impression has captured the ball abutments exactly. The ball abutment laboratory implant is then repositioned in the impression in the laboratory. The flat surfaces beneath the ball ensure that the implant axes are transferred precisely. Then pour dental stone into the impression. Cast the dental arch as usual and create the base. The laboratory implant must fit precisely on the model. A ball abutment laboratory implant is available for all series of abutments as all balls have a diameter of 2.25 mm.

Processing and integration of the following matrix/inner matrix into the denture should be carried out using the corresponding components and according to the manufacturer's instructions.

REF 386-630-00

Dalbo®-PLUS matrix TE basic complete

REF 386-631-00

Dalbo®-PLUS Lamella retention insert E

REF 387-663-00

Dalbo®-PLUS screwdriver / activator

**Recall:**

Dentures and their retention units must be monitored at six-monthly intervals. Take the following points into account:

- Eliminate unfavorable movements of the denture (optimize relining of the denture, activate or replace the matrices)
- Check the fit of the ball abutments on the implant (tighten if necessary)
- Oral hygiene (remove plaque and dental calculus and, if necessary, re-instruct the patient on cleaning implants)

**Tightening torque**

Ball abutment on the model: / in the mouth

35 Ncm

**7.2 tioLogic® ST tioLOC abutments**

The tioLOC-technique involves an implant-retained, tissue-borne restoration with a very low profile height. Existing coverdentures can be modified or remade. The tioLOC abutment can be used for tissue-borne coverdentures with two to four implant abutments. The intermaxillary space should be at least 4.0 mm. The tioLOC components must not be modified. Only the tightening forces of the matrices can be adjusted as required using the exchangeable retention inserts.

1.0 mm /  
2.0 mm /  
3.0 mm /  
4.0 mm /  
5.0 mm

The tioLogic® ST tioLOC abutments are available for the S, M and L series of abutments in gingival heights of 1.0 mm, 2.0 mm, 3.0 mm, 4.0 mm and 5.0 mm. The gingival height refers to the cylindrical section of the tioLOC abutment. Its upper edge should be placed approx. 1.0 mm above the gingival line. To achieve optimum retention, the tioLOC abutments should be positioned at the same level. The head with retainer element is identical for all tioLOC

abutments. The tioLOC abutment is inserted with the insertion key tioLOC abutment. All tioLOC abutment components may not be modified.

**Working procedures**

In the direct technique, process the matrix into an existing denture directly in the patient's mouth without fabricating a model. Neither an impression post nor a laboratory implant are required. Attach the tioLogic® ST tioLOC abutments intra-orally to the corresponding tioLogic® ST implants. During the clinical procedure, protect the undercuts with the help of a rubber dam. Place this over the matrix. This prevents the undercuts in the mouth from filling with cold-curing polymer and the denture from becoming impossible to remove after curing. Observe the instructions for use of the resin.

In the indirect technique, take an impression and fabricate a model. The matrix is processed into a denture in the laboratory.

Processing and integration of the following matrices/caps into the denture should be carried out using the corresponding components and according to the manufacturer's instructions.

REF 346-512-00		Impression caps, 2 pieces	∅ 4.0 mm L 10.0 mm
REF 346-514-00		Laboratory analogs 2 pieces	∅ 4.0 mm L 10.0 mm
REF 346-516-00		Set of retention caps: 1 stainless steel case for caps, 1 protective plate, 1 laboratory cap, 3 retention caps (with internal pin)	
REF 346-518-00		Black laboratory caps (with pin), 4 pieces	GH 1.0 mm
REF 346-520-00		Yellow laboratory caps (with pin), 4 pieces	700 g 0° - 20°
REF 346-522-00		Pink laboratory caps (with pin), 4 pieces	900 g 0° - 20°
REF 346-524-00		Transparent laboratory caps (with pin), 4 pieces	1.500 g 0° - 20°
REF 346-526-00		Violet laboratory caps (with pin), 4 pieces	2.700 g 0° - 20°
REF 346-530-00		Protective plate 10 pieces	
REF 347-745-00		Insertion/removal tool for elastic caps	

**Recall:**

Dentures and their retention units must be monitored at six-monthly intervals. Take the following points into account:

- Eliminate unfavorable movements of the denture (optimize relining of the denture, activate or replace the matrices)
- Check the fit of the ball abutments on the implant (tighten if necessary)
- Oral hygiene (remove plaque and dental calculus and, if necessary, re-instruct the patient on cleaning implants)

**Tightening torque**

Ball abutment on the model: / in the mouth	35 Ncm
--	--------

**7.3 Additional components**

The prosthetic screw (not included in scope of delivery) serves only to help the dental technician fix the abutment on the working cast temporarily.

Always use the prescribed torque for the final fitting and fixation of the abutment in the mouth.

Overview of additional components:

Torque ratchet
tioLogic® ST ball abutment insertion key, ratchet
tioLOC Insertion key, ratchet

**8. Technical specification  
Composition**

Product	Composition
tioLogic® ST ball abutment	3.7165
Ball abutment laboratory implant	Al alloy
tioLogic® ST tioLOC abutment	3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-06), or see [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu). The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaurum.com](http://www.dentaurum.com).

**9. Adverse reactions and residual risks**

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

**10. Safety information and precautions**

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® ST prosthetic components in combination with tioLogic® ST implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

**Severe incidents:**

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

**11. Storage and shelf life**

Store the tioLogic® ST implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

**12. Disposal**

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

**13. Information for single use products**

Prosthetic abutments and components are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

#### 14. Information for multiple use products

The torque ratchet, ball abutment insertion key, tioLOC insertion key ratchet, hex key and the ISO shank-hex/ratchet adapter are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Scan abutments, titanium, must be replaced if they cannot be clearly identified or if their function is impaired due, for example, to poor readability of the markings or labels.

#### 15. Information on disinfection / cleaning / sterilizing

The tioLogic® ST prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaurem generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

#### 16. Product overview and scope of delivery

Please refer to the tioLogic® ST product catalog for the tioLogic® ST product range.

#### 17. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

#### 18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

#### 19. Information on MRT safety

The products within these Instructions for use are MRT safe.

#### 20. Additional information

Additional information can be found at [www.dentaurem.com/ifu](http://www.dentaurem.com/ifu)

Instructions for use	REF
tioLogic® ST Additional components	989-805-87
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurem products	989-804-43

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