DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative. As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

3. Area of application

3.1 Bridge abutment tiologic® ST

Bridge abutments are prosthetic components for bridge and bar restorations with tiologic® ST implants. They can compensate implant divergences of up to 40°, thereby enabling bridge and bar restorations to have a passive fit.

4. Medical indications/contraindications

Indications

- Anchorage of implant-supported, splinted bridges in the maxilla in combination with at least 6 tiologic® ST implants and, in the mandible, in combination with at least 4 or more tiologic® ST implants

Contraindications.

- For biomechanical reasons, we recommend that the maximum angulations for bridge abutments are not exceeded
- Reconstructions with mixed support (tooth/bridge or bar)
- Combination with other anchorage elements
- Grinding and shortening of bridge abutments

5. Information on use, processing, fitting and insertion

5.1 Bridge abutment tiologic® ST

Bridge abutments are available for the S, M and L series of abutments. The gingival heights are 1.0 mm, 2.5 mm and 4.0 mm. The seating surface of the bridge abutments should be approx. 0.5 mm above the gingiva. The seating surface and cone (20 °) for the caps on the bridge abutments is identical (ø 4.1 mm) for the S, M and L series of abutments. The same AnoTite screw (L 6.0 mm) is used for all caps. The bridge abutments should **NOT** be ground or shortened. When handling the bridge abutments, ensure that the occlusal screw is not damaged by tweezers etc. After lasering, casting or digital processing and prior to placing it intra-orally, every bridge restoration must be checked for passive fit on the model ("Sheffield Test"). This involves placing the bridge restoration on the bridge abutments and fixing it to the bridge abutment with only one prosthetic screw. If the construction is raised to create a gap between the bridge and the construction, stresses are present and must be eliminated.

5.2 Impression post bridge, open and closed tiologic® ST

To achieve stable apposition of the gingiva with bridge abutments, it is possible to take an open or closed impression on a bridge abutment that is fitted in the mouth. Since the impression posts bridge have no rotational security, they may only be used for bridge and bar constructions.

5.3 Closure screw bridge tiologic® ST

The closure screw bridge is screwed directly onto the bridge abutment after impression-taking and protects the surface for the caps on the bar abutments.

5.4 Plastic cap, bridge tiologic® ST

The plastic cap, that burns residue-free, can be individually adapted to the space available. A base structure is fabricated as a strengthener for a long-term restoration. The wax-up is fabricated taking the lingual and palatal overcasts into consideration. The waxes from Dentaurum which burn without leaving a residue (StarWax range) are suitable. This procedure guarantees that there is still sufficient space for subsequent working stages and the pre-prosthetic planning can be maintained.

5.5 Titanium cap, bridge tiologic® ST

The titanium cap bridge can also be individually adapted to the space available. Titanium must not become overheated during preparation as this can lead to differences in surface hardening (alpha-case layer). This could complicate or handicap further processing. A wax set-up is then fabricated, which can be checked using the lingual and palatal overcasts. Before waxing up the temporary restoration, ensure that there is an adequate, uniform cement gap between the bridge titanium caps and the temporary restoration by blocking out using preparation and casting wax. This guarantees stress-free fixation. Finishing, trimming and polishing should be completed according to the instructions of the acrylic manufacturer. PTFE cylinder pins are available for restorations which are bonded in the laboratory. The pins do not bond with the adhesive and so prevent the adhesive entering the screw aperture.

5.6 Scan cap titanium bridge tiologic® ST

The surface of the scan cap bridge titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. For this purpose, the scan cap bridge is screwed onto the bridge abutment and scanned. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer.

5.7 Laboratory implant bridge tiologic® ST for printed and cast models

The laboratory implant can be used in both analog and digital workflow. The cavities of the laboratory analogs are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory analogs with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary. In order for the counter screw supplied to fix the laboratory analog at the apical end in the printed model, the model base should be printed with a sufficient height. If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory analog in the cervical area of the printed model. The 3 retention tabs should still be visible in the printed model for monitoring purposes.

When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended series of abutments S, M or L. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The closure screw supplied with the laboratory implant enables the laboratory implant to be fixed in the end position from apical.

The same AnoTite screw (L 6.0 mm) is used for all caps. A new AnoTite screw should be used for the final fitting in the mouth. The bridge abutment is placed with the bar/bridge/AngleFix abutment insertion key (secure the counter screw!).

Tightening torque

- Bridge abutment on model: manually
- Bridge abutment intra-orally: 35 Nm
- Cap on bridge abutment on the model: manually
- Cap on bridge abutment intra-orally: 25 Nm
- Sure-grip screw impression post on bridge abutment intra-orally: 15 Nm
- Closure screw in bridge abutment intra-orally: 15 Nm

Attention: Do not confuse the prosthetic screw with the AnoTite screw.

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- Bridge abutment: Titanium alloy (Ti6Al4V)
- Caps for bridge abutment:
 - Plastic cap bridge: Polycarbonate (PC)
 - Titanium cap bridge: Titanium alloy (Ti6Al4V)
 - Closure screw bridge: Titanium alloy (Ti6Al4V)
 - Impression post bridge, open: Titanium alloy (Ti6Al4V) / Screw: 1.4305
 - Impression post bridge, closed:
 - Impression post and screw: Titanium alloy (Ti6Al4V) / Impression cap and cap for bite registration: POM
- Scan cap titanium bridge: Titanium alloy (Ti6Al4V)

• Laboratory implant bridge: Titanium alloy (Ti6Al4V)

• Prosthetic screw for bridge and bar abutment:

– AnoTite screw: Titanium alloy (Ti6Al4V)

– Prosthetic screw: Titanium alloy (Ti6Al4V)

- Insertion key bar/bridge/AngleFix: 1.4197

- Hex key:

– Hex key - ISO shank 1.3: 1.4197

– Hex key - ratchet 1.3: 1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.
- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

Bridge abutments and their caps are for single use only. Reconditioning of a bridge abutment or its caps (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The insertion key bar/bridge/AngleFix, the hex key 1.3 and the scan cap bridge titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview/scope of delivery, incl. REF

Please refer to the tiologic® ST product catalog for the tiologic® ST product range (REF 989-965-20).

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

 Please refer to the label. Additional information can be found at www.dentaurum.com/fu (Explanation of symbols REF 989-313-00)

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