
ceraMotion®
One Touch

MY FIRST
CERAMIC
IN PASTE.

CERAMOTION® WORKBOOK

MY FIRST
CERAMIC
IN PASTE...

ceraMotion®
One Touch

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MY FIRST
CERAMIC
IN PASTE . . .

ceraMotion®
One Touch

*"ceraMotion® One Touch
makes the impossible possible –
and it's simple and easy."*

Nancy CG

ceraMotion®
One Touch

Veneers 14 to 24 made of lithium disilicate ceraMotion® LiSi LT A1
and finished with ceraMotion® One Touch 2D pastes.

Pasta 2D
Body A

Pasta 2D
White

Pasta 2D
Grey

SKETCH

ceraMotion®
One Touch

© Nancy Ceballos G.

BEFORE AND AFTER

MY FIRST
CERAMIC
IN PASTE...

ceramotion®
One Touch

*"When the fine tip of our brush
comes in contact with the pastes,
it creates a symphony of forms and fissures.
It also produces the subtle effects."*

CHRISTIAN FERRARI

Christian Ferrari





ceraMotion®
One Touch

© Christian Ferrari

ZIRCONIA FRAMEWORK ON IMPLANTS

ceraMOTION®
One Touch

© Christian Ferrari

VENEERING WITH CERAMOTION® ZR

ceraMotion®
One Touch

© Christian Ferrari

RESULT AFTER FIRST BAKE

ceraMotion®
One Touch

© Christian Ferrari

MICRO LAYERING WITH CERAMOTION® ONE TOUCH

2D VIOLET 50%

3D OPAL HONEY 30%

2D BODY L 30%

3D IO BRIGHT 60%

3D HONEY 40%

3D BLUE 50%

3D TRANSPA 70%

3D TRANSPA 70%

3D TRANSPA 40%

3D TRANSPA 60%

MICRO LAYERING WITH CERAMOTION® ONE TOUCH

ceramotion®
One Touch

© Christian Ferrari

FINAL RESULT

MY FIRST
CERAMIC
IN PASTE...


ceraMotion®
One Touch

Waldemar Fritzer

A stylized, handwritten signature of the man in the polaroid photo, written in white ink on a black background.

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






MY FIRST
CERAMIC
IN PASTE . . .

*"ceraMotion®,
yes I love this product line
and you will like it too"*

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

MY FIRST
CERAMIC
IN PASTE . . .

"Passion for detail" is our slogan, since it embraces everything that represents us. We believe that in our hands lies an opportunity to change lives, restore self-esteem and security to a person. Therefore, it is so important for us to make unique pieces, where details make the difference. I found in ceraMotion® One Touch the perfect partner to achieve this goal.

JM
J O S É M E L O

ceraMotion®
One Touch

ceraMotion®
One Touch

ceraMotion®
One Touch

P 2D / INCISAL I2 / PASTE 2D

P 2D / BODY L

P 3D / INCISAL OPAL BRIGHT 50%
+ P 3D DENTIN BLEACH 50%

© José Melo

STEP 1

MATERIALS USED: IVOCLAR DISILICATE HTA1

ceraMotion®
One Touch

ANTES

/

DESPUÉS

STEP 2

ceraMotion®
One Touch

Our case consists of crowns in lithium disilicate, finished and characterized with ceraMotion® OneTouch No Limits offering its great potential in mass diversity, fluorescence, and enhancement of mamelons in the incisal area, opalescent effects, and halos. Also, in occlusal faces to provide dynamism, naturalness, and beauty to the pieces.

ceraMotion®
One Touch

P 2D / INCISAL I2

ST1 WHITE

P 3D DENTIN LIGHT

STEP 3

MY FIRST
CERAMIC
IN PASTE . . .

ceraMotion®
One Touch

*"ceraMotion® One Touch is an easy way
to make my restorations natural and vivid"*

Pavel Moysievich

FRAMEWORK BEFORE LAYERING

ceraMotion®
One Touch

© Pavel Moysievich


In the case you can see the characterization and final layering of 12-unit zirconium framework on 6 implants. For the framework, I use Nacera MS AL zirconium and for the individualization, I use body pasts and effects from ceraMotion® One Touch, One Touch No Limits and One Touch Pink kits.

Prior to the application of ceraMotion® One Touch products, the framework was sand blasted with 100-micron sand under 3-bar pressure.

STEP 1

LAYERING SCHEME 1

ceraMotion®
One Touch

© Pavel Moysievich

The whole framework was covered with the thin layer of ceraMotion® Paste Glaze.

At the front I use body A for the teeth bulk with the thin layer of ceraMotion® 3D opal grey and 2D blue effects at the incisal area as marked on the application scheme.


In order to intensify the translucency effect I also added the dots of ST 13 black stain to the incisal area.

At the distal areas apart from the shades mentioned above, I also apply cM body R at the gingival margin.

STEP 2

LAYERING SCHEME 2

ceramotion®
One Touch

© Pavel Moysievich

At the gum area, the dark spots were covered with Modifier past and Gingiva 3 was used for the lighter areas.
The first bake was done.

After the first bake, I apply the thin layer of Glaze to the whole prosthesis and add a touch of Gingiva 3 at the gum.
The second bake was performed.

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

MY FIRST
CERAMIC
IN PASTE . . .

ceraMotion®
One Touch

"We evaluated ceraMotion® One Touch and Pink as the best material to achieve the necessary results, due to its great aesthetic quality of completion, possibilities of anatomical and gingival customization, depth features, and its incomparable agility."

ceraMotion®
One Touch

THE SELECTION OF ZIRCONIA IS FUNDAMENTAL IN THIS TYPE OF WORK BEING NECESSARY TO COMBINE RESISTANCE AND TRANSLUCENCY.

© Clayton Schäfer

STEP 1

ceraMotion®
One Touch

© Clayton Schäfer

AFTER SINTERING

AFTER FINISHING, SANDBLASTING AND STEAM
CLEANING WE START BY APPLYING A THICK
LAYER OF GLAZE INTERNALLY TO THE PART
ENSURING AN EXTREMELY POLISHED SURFACE.

ceramotion®
One Touch

© Clayton Schäfer

STEP 2

THEN WE APPLY A THIN LAYER OF GLAZE ON THE
TEETH AND START THE MAKEUP WITH 2D PASTES.

ceramotion®
One Touch

© Clayton Schäfer

THE 2D PASTES USED WERE GLAZE, BODY A, BLUE AND GREY MIXED,
AND HONEY AND RED BROWN ALSO MIXED TO CREATE SOME VESTIBULAR
EFFECTS. WE USE DILUTING LIQUID FOR 2D PASTES.

STEP 3

APPLICATION OF 3D GINGIVAL PASTE
CHARACTERIZED WITH RED PIGMENT.

ceramotion®
One Touch

STEP 4

ceraMotion®
One Touch

ceraMot

© Clayton Schäfer

AFTER FIRST BAKE

APPLICATION OF GINGIVA ACCORDING TO
THE SELECTED COLOR AND WE COMPLEMENT
THE INCISAL WITH 3D OPAL BLUE,
OPAL GREY AND NEUTRAL PASTES MIXED
IN EQUAL PROPORTIONS TO GENERATE DEPTH
AND GREY OPALESCENCE IN THE INCISAL
REGION.

ceramotion®
One Touch

© Clayton Schäfer

STEP 5

ceraMotion®
One Touch

© Clayton Schäfer

FINISHED RESULT

MY FIRST
CERAMIC
IN PASTE . . .

ceramotion®
One Touch


Juvenal Souza Neto *oral design*

"In a few touches we can achieve an incredible result"

Master Ceramist 🇧🇷 🇺🇸
Dentogingival Concept Developer
jsnoraldesign@gmail.com

FIRST BAKE

ceraMotion®
One Touch

ENAMEL MIX

MIXING COLORS TO ACHIEVE ENAMEL QUICKLY:

P 3D OPAL BLUE + P 2D INCISAL I3 +
P 3D NEUTRAL

DESIGN STEP 1

FIRST BAKE

ceramotion®
One Touch

Camila Maia

USING THE

P 3D DENTIN DARK


TO ADD MORE VOLUME
ON THE ROOT AND
CERVICAL AREA.

ENAMEL MIX TO CREATE
A TRANSLUCENCY AREA
AND ADJUST THE SHAPE

STEP 1

SECOND BAKE

ceramotion®
One Touch

MAMELON MIX

MIXING COLORS TO ACHIEVE MAMELON'S QUICKLY:

P 2D HONEY + P 3D WHITE I3 +
P 3D DENTIN LIGHT

DESIGN STEP 2

SECOND BAKE

ceramotion®
One Touch

FIRST BAKE RESULT:
TRANSLUCENCY BASE AND ROOT SHAPE.

SECOND BAKE:
MAMELONS AND GINGIVA MODIFIER

P 3D GINGIVA 4 AND GINGIVA 2

STEP 2

GINGIVA TEXTURE


RECREATE GINGIVA TEXTURE USING A TISSUE PAPER.

ceramotion®
One Touch

STEP 3

ceramotion®
One Touch

FINAL RESULT!

MY FIRST
CERAMIC
IN PASTE...

*“You can use ceraMotion®
in every type of cases,
simple to difficult.”*

T. Wakatake

Patients who care about teeth alignment and its color tone

Since the patient requested to have as much beautiful teeth as possible without characterization, I painted the basic shade with 2D body A Incisal I3 and fired it for the first time.

Then I glazed with Paste Glaze bright and completed two firings.

In the first firing, the color tone was adjusted while controlling the thickness of the stain with the Refreshing Liquid.

In the second firing, the Paste Glaze bright was thinly applied so as not to lose the quality of its surface by using the Diluting Liquid.

SIMPLE CASE

BEFORE
AND AFTER

ceramotion®
One Touch

© Tetsuya Wakatake

S I M P L E C A S E

FULL-JAW PROSTHESIS AFTER ADJUSTMENT
OF DEFECTIVE PROSTHESIS AND RECONSTRUCTION
OF OCCLUSION

ceramotion®
One Touch

© Tetsuya Wakatake

Primary firing

The purpose of the primary firing was to provide the basic shade and the transparent edge using 2D body A 2D Incisal I3 2D violet + ST 13 black.

Secondary firing

The purpose was to give the natural view which natural teeth have. 2D white ST white 2D honey ST orange ST red brown was used for its characterization.

Third firing

I coated Dlluting Liquid thinly not to be erased the properties of its surface and baked.

Conclusion

The number of firings can be reduced by using THIS stain. It is wonderful that the stain can be thickened or thinned by using Refreshing Liquid and Dlluting Liquid appropriately.

D I F F I C U L T C A S E

DIFFICULT CASE

BEFORE
AND AFTER

ceraMotion®
One Touch

© Tetsuya Wakatake

DIFFICULT CASE

WWW.DENTAURUM.COM

D
DENTAURUM

Turnstr. 31 | 75228 Ispringen | Germany | Telefon +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaurum.com | info@dentaurum.com