## DE

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### Abbreviated instructions for use

**Dear customer,**

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

#### 1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

#### 2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient's needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

#### 3. Area of application

##### 3.1 Bar abutment tiologic® ST

Bar abutments are prosthetic components for bar restorations with tiologic® ST implants. They can compensate implant divergences of up to 80°, thereby enabling bridge and bar restorations to have a passive fit.

##### 4. Medical indications/contraindications

Indications

- Anchorage of implant-supported, splinted bars in the maxilla in combination with at least 6 tiologic® ST implants and, in the mandible, in combination with at least 4 or more tiologic® ST implants

Contraindications.

- For biomechanical reasons we recommend that the maximum angulations for bar abutments are not exceeded.

- Reconstructions with mixed support (tooth/bridge or tooth/bar)
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Combination with other anchorage elements
- Grinding and shortening bar abutments

#### 5. Information on use, processing, fitting and insertion

##### 5.1 Bar abutments tiologic® ST

Bar abutments are available for the S, M and L series of abutments. The gingival heights are 1.0 mm, 2.5 mm, 4.0 mm and 5.5 mm. The seating surface of the bar abutments should be approx. 0.5 mm above the gingiva. The seating surface for the caps on the bar abutments is identical (ø 4.5 mm) for the S, M and L series of abutments. The same AnoTite screw (L 6.0 mm) is used for all caps. The bar abutments should **NOT** be ground or shortened. When handling the bar abutments, ensure that the occlusal screw is not damaged by tweezers etc. After lasering, casting or digital processing, and prior to placing it intra-orally, every bar restoration must be checked for passive fit on the model ("Sheffield Test"). This involves placing the bar restoration on the bar abutment and fixing it to the bar abutment with only one prosthetic screw. If this raises the bar to create a gap between the bar and bar abutment, stresses are present and must be eliminated.

##### 5.2 Impression post bar, open tiologic® ST

To achieve stable apposition of the gingiva with bar abutments, it is possible to take an open impression on a bar abutment that is fitted in the mouth. Since the impression posts bar have no rotational security, they may only be used for bridge and bar constructions.

##### 5.3 Closure screw bar tiologic® ST

The closure screw bar is screwed directly onto the bar abutment after impression-taking and protects the surface for the caps on the bar abutments.

##### 5.4 Plastic cap bar tiologic® ST

The plastic cap, that burns residue-free, can be individually adapted to the space available. A prefabricated plastic bar is waxed to the plastic caps for the casting technique. Casting, trimming and polishing should be carried out according to the alloy manufacturer's instructions. A countersink or cutter are available for the post-processing of the screw aperture or the screw seating. The plastic caps are bonded to the acrylic for the milling technique. Following this, the pattern for the individual bar is fabricated and paralleled with a milling machine. For enhanced retention, horizontal and vertical attachments can be integrated. Casting, trimming and polishing should be carried out according to the alloy manufacturer's instructions. After casting, the reamer is used for finishing the screw aperture, the countersink is used for the screw seat and the cutter for the seating surface of the plastic caps.

##### 5.5 Titanium cap bar, laser weldable tiologic® ST

The prefabricated titanium bars are shortened and fixed in the correct horizontal axis for laser-welding with a parallelometer / fixator (e.g. Paralas, Dentaurum). The bar components and titanium caps are then lasered together. The only filler material

which can be used is pure titanium wire (e.g. Dentaurum). Should further questions arise on laser-welding and processing titanium, please call your local representative. After laser-welding, the titanium bar should be trimmed and polished (e.g. with a rematitan® finishing kit from Dentaurum).

##### 5.6 Titanium cap bar, adhesive technique tiologic® ST

Adhesive titanium caps are available for adhering cast bars intra-orally. Once a cast bar has been produced, it is bonded intra-orally to the titanium bonding base. Adhering the titanium caps can compensate for imprecise casting and possible differences during impression-taking or model casting.

##### 5.7 Titanium cap bar tiologic® ST

The titanium cap bar can also be individually adapted to the space available.

Titanium must not become overheated during preparation as this can lead to differences in surface hardening (alpha-case layer). This could complicate or handicap further processing. A wax set-up is then fabricated, which can be checked using the lingual and palatal overcasts. Before waxing up the temporary restoration, ensure that there is an adequate, uniform cement gap between the titanium cap bar and the temporary restoration by blocking out using preparation and casting wax. This guarantees stress-free fixation. Finishing, trimming and polishing should be completed according to the instructions of the acrylic manufacturer. PTFE cylinder pins are available for restorations which are bonded in the laboratory. The pins do not bond with the adhesive and so prevent the adhesive entering the screw aperture.

##### 5.8 Scan cap titanium bar tiologic® ST

The surface of the scan cap bar titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. For this purpose, the scan cap bar is screwed onto the bar abutment and scanned. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer.

##### 5.9 Laboratory implant bar tiologic® ST for printed and cast models

The laboratory implant can be used in both analog and digital workflow. The cavities of the laboratory analogs are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory analogs with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary. In order for the counter screw supplied to fix the laboratory analog at the apical end in the printed model, the model base should be printed with a sufficient height. If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and holding function of the laboratory analog in the neck area of the printed model. The 3 retention tabs should still be visible in the printed model for monitoring purposes. When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended series of abutments S, M or L. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The closure screw supplied with the laboratory implant enables the laboratory implant to be fixed in the end position from apical.

The same AnoTite screw (L 6.0 mm) is used for all caps. A new AnoTite screw should be used for the final fitting in the mouth. The bar abutment is placed with the bar, bridge, AngleFix abutment insertion key (secure the counter screw).

#### Tightening torque

- Bar abutment on model: manually
- Bar abutment intra-orally: 35 Ncm
- Cap on bar abutment on the model: manually
- Cap on bar abutment intra-orally: 25 Ncm
- Sure-grip screw impression post on bar abutment intra-orally: 15 Ncm
- Closure screw in bar abutment intraorally: 15 Ncm

**Attention: Do not confuse the prosthetic screw with the AnoTite screw.**
The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

#### 6. Composition

■ Bar abutment: Titanium alloy (Ti6Al4V)

■ Cap for bar abutment:

– Plastic cap bar: Polycarbonate

– Titanium cap bar (all): Titanium alloy (Ti6Al4V)

– Closure screw bar: Titanium alloy (Ti6Al4V)

– Impression post bar, open: Titanium alloy (Ti6Al4V) / Screw: 1.4305

– Scan cap titanium bar: Titanium alloy (Ti6Al4V)

■ Laboratory implant bar: Titanium alloy (Ti6Al4V)

■ Prosthetic screw for bridge and bar abutment:

– AnoTite screw: Titanium alloy (Ti6Al4V)

– Prosthetic screw: Titanium alloy (Ti6Al4V)

■ Insertion key bar/bridge/AngleFix: 1.4197

■ Hex key:

- Hex key - ISO shank 1.3: 1.4197

- Hex key - ratchet 1.3: 1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.
- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

Store in a dry place.

Bar abutments and their caps are for single use only. Reconditioning of a bar abutment or its caps (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

The insertion key bar/bridge/AngleFix, the hex key 1.3 and the scan cap bar titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again.

**11. Information on disinfection/cleaning/sterilizing**

tioLogic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentsaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

## 12. Product overview / scope of delivery, incl. REF

### 13. Quality information

for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the processing on site.

 Please refer to the label. Additional information can be found at [www.dentaurum.com](http://www.dentaurum.com) (Explanation of symbols REF 989-313-00)

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