ceraMotion® Z  
The symbol for smarter zirconia

---

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Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

## General description

ceraMotion® Z are dental zirconium oxide blanks. They have been developed specifically for the fabrication of prosthetic restorations using subtractive manufacturing. They fulfill the requirements for use in dental medicine due to their excellent chemical, mechanical and physical properties. These Instructions for use contain the main steps for processing ceraMotion® Z blanks.

### These Instructions for use are valid for the following products:

- ceraMotion® Z Hybrid (3Y/4Y-ZrO<sub>2</sub>)
- ceraMotion® Z Cubic Multishade (5Y-ZrO<sub>2</sub>)
- ceraMotion® Z HT Multishade (3Y-ZrO<sub>2</sub>)
- ceraMotion® Z HT Shade (3Y-ZrO<sub>2</sub>)
- ceraMotion® Z HT White (3Y-ZrO<sub>2</sub>)
- ceraMotion® Z White (3Y-ZrO<sub>2</sub>)

## Limitation of liability

All the information in these instructions for use takes into consideration the applicable norms and regulations, state-of-the-art technology and our many years of knowledge and experience. The manufacturer provides no warranty and will accept no liability in the following cases:

- Non-compliance with this document.
- Non-compliance with or unauthorized changes to the procedure as stated for the sintering process.

## 1. Composition

The concentrations listed refer to the respective mass fraction of the oxides and are given in mass percent (m/m). The sum of all elements equals 100%.

| Chemical composition  |                               |                                |              |
|---|-------------------------------|--------------------------------|--------------|
| ZrO <sub>2</sub> + HfO <sub>2</sub> + Y <sub>2</sub> O <sub>3</sub>                                     | Y <sub>2</sub> O <sub>3</sub> | Al <sub>2</sub> O <sub>3</sub> | other oxides |
| <b>ceraMotion® Z White, ceraMotion® Z HT White, ceraMotion® Z HT Shade, ceraMotion® Z HT Multishade</b> |                               |                                |              |
| > 98%   | 4.5% – 6%                     | < 0.5%                         | < 1.5%       |
| <b>ceraMotion® Z Cubic Multishade</b>   |                               |                                |              |
| > 98%   | 8.2% – 9.7%                   | < 0.5%                         | < 1.5%       |
| <b>ceraMotion® Z Hybrid</b>   |                               |                                |              |
| > 98%   | 5.8% – 9.7%                   | < 0.5%                         | < 1.5%       |

## 2. Technical data

The coefficient of thermal expansion (CTE) for all ceraMotion® Z blanks amounts to 10.5 ± 0.5 · 10<sup>-6</sup> K<sup>-1</sup> in the temperature range 25 °C – 500 °C / 77 °F – 932 °F and was determined according to EN ISO 6872.

| Technical data   |                              |   |
|--|------------------------------|---|
| Flexural strength (MPa)                                    | Density (g/cm <sup>3</sup> ) | Chemical solubility (µg/cm <sup>2</sup> ) |
| <b>ceraMotion® Z White, ceraMotion® Z HT White</b>         |                              |   |
| ≥ 800  | > 6.06                       | < 10                                      |
| <b>ceraMotion® Z HT Shade, ceraMotion® Z HT Multishade</b> |                              |   |
| ≥ 800  | > 6.08                       | < 10                                      |
| <b>ceraMotion® Z Cubic Multishade</b>                      |                              |   |
| ≥ 600  | > 6.04                       | < 10                                      |
| <b>ceraMotion® Z Hybrid</b>                                |                              |   |
| ≥ 800  | > 6.04                       | < 10                                      |

## 3. Intended purpose

Dental ceramic semi-finished products are intended for the production of dental restorations to restore both the chewing function and aesthetics.

## 4. Indications

ceraMotion® Z are zirconium oxide blanks for the production of fixed restorations, such as inlays, onlays, veneers, single crowns and bridges.

| Type of restoration                     |                      |                    |         |        |        |                                      |
|---|----------------------|--------------------|---------|--------|--------|--------------------------------------|
| Single crowns                           | Bridge up to 3 units | Multi-unit bridges | Veneers | Inlays | Onlays | Cantilever bridge unit premolar size |
| <b>ceraMotion® Z Hybrid **</b>          |                      |                    |         |        |        |                                      |
| •                                       | •                    | • <sup>1)</sup>    | •       | •      | •      | •                                    |
| <b>ceraMotion® Z Cubic Multishade *</b> |                      |                    |         |        |        |                                      |
| •                                       | •                    |                    | •       | •      | •      |                                      |
| <b>ceraMotion® Z HT Multishade **</b>   |                      |                    |         |        |        |                                      |
| •                                       | •                    | • <sup>2)</sup>    | •       | •      | •      | •                                    |
| <b>ceraMotion® Z HT Shade **</b>        |                      |                    |         |        |        |                                      |
| •                                       | •                    | • <sup>2)</sup>    | •       | •      | •      | •                                    |
| <b>ceraMotion® Z HT White **</b>        |                      |                    |         |        |        |                                      |
| •                                       | •                    | • <sup>2)</sup>    | •       | •      | •      | •                                    |
| <b>ceraMotion® Z White **</b>           |                      |                    |         |        |        |                                      |
| •                                       | •                    | • <sup>2)</sup>    | •       | •      | •      | •                                    |

1) Up to two pontics

2) Up to 3 sequential pontics in the anterior region, up to 2 sequential pontics in the posterior region

\* Type II Class 4 dental ceramics acc. to EN ISO 6872

\*\* Type II Class 5 dental ceramics acc. to EN ISO 6872

## 5. Contraindications and adverse reactions

There is no evidence of incompatibilities with this material when manufacturing or processing in accordance with the Instructions for use.

Contraindications for the zirconium oxide ceramics Type II Class 5 (White, HT White, HT Shade, HT Multishade, Hybrid):

- Production of abutments.
- Production of bridges with multiple segments or smaller diameters than specified.
- They must not be used on patients with bruxism.
- They must not be used if the patient is allergic to any of the constituents.

Contraindications for the zirconium oxide ceramics Type II Class 4 (Cubic Multishade):

- Production of abutments.
- Production of bridges with multiple segments or smaller diameters than specified.
- Production of cantilever bridges, Maryland bridges, inlay or onlay bridges.
- They must not be used on patients with bruxism.
- They must not be used if the patient is allergic to any of the constituents.

## 6. Patient population

Patients with defects in tooth structure for whom a restoration is necessary and indicated.

## 7. Intended users

Processing must be carried out only by dental professionals with knowledge of the production of dental restorations. Non-observance of this advice may lead to failure of the restoration and other complications.

## 8. Warnings and precautions to be taken

To avoid damage, make sure that a collision with hard objects is avoided. It must also be ensured that the material is neither crushed nor exposed to vibrations. The product must not be used if it shows cracks, splits, fractures or color irregularities.

Mechanical processing of the material can generate dust; therefore, a suitable extraction system is required when processing zirconium dioxide ceramics. In addition, processing generates heat, so that processed material may be hot. We generally recommend that gloves, protective clothing, safety glasses and especially respiratory protection equipment are worn.

Professional handling is imperative. The processing of the product is the sole responsibility of trained professionals. The patient must not come into contact with this product before the sintering process has been completed.

## 9. Reuse

Finished restorations which have been produced using semi-finished products are intended for single use only.

## 10. Storage and/or handling conditions

The product properties are not affected by normal changes in environmental conditions (such as temperature, pressure, humidity or light).

## 11. Disposal of residual material

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

## 12. Using the products

### 12.1 Construction

| ceraMotion® Z White, HT White, HT Shade, HT Multishade |                  |           |                              |
|--|------------------|-----------|------------------------------|
| Type of restoration                                    |                  | Framework | Fully anatomical restoration |
| Minimal wall thickness in mm                           | anterior region  | 0.5       | 0.8                          |
| Maximum number of connecting links                     |                  | 3         | 3                            |
| Connector cross-section in mm <sup>2</sup>             |                  | ≥ 7       | ≥ 9                          |
| Minimal wall thickness in mm                           | posterior region | 0.6       | 1.0                          |
| Maximum number of connecting links                     |                  | 2         | 2                            |
| Connector cross-section in mm <sup>2</sup>             |                  | ≥ 9       | ≥ 9                          |

| ceraMotion® Z Hybrid, Cubic Multishade     |           |                              |                  |        |
|--|-----------|------------------------------|------------------|--------|
| Type of restoration                        | Hybrid    |                              | Cubic Multishade |        |
|  | Framework | Fully anatomical restoration | Crown            | Bridge |
| <b>anterior region</b>                     |           |                              |                  |        |
| Minimal wall thickness in mm               | 0.5       | 0.8                          | 0.8              | 1.0    |
| Maximum number of connecting links         | 2         | 2                            | –                | 1      |
| Connector cross-section in mm <sup>2</sup> | ≥ 9       | ≥ 9                          | –                | ≥ 12   |
| <b>posterior region</b>                    |           |                              |                  |        |
| Minimal wall thickness in mm               | 0.6       | 1.0                          | 1.0              | 1.0    |
| Maximum number of connecting links         | 2         | 2                            | –                | 1      |
| Connector cross-section in mm <sup>2</sup> | ≥ 9       | ≥ 9                          | –                | ≥ 16   |

**Note:** Restorations must always be designed to be ceramic-compatible, i.e., without sharp edges.

### 12.2 Milling

ceraMotion® Z blanks are a semi-finished product for the production of all-ceramic dental restorations by subtractive manufacturing. Machining can be done with most commercially available CAD/CAM milling machines or with a manual milling machine. Suitable tools should be selected in accordance with the instructions in the instructions for use provided by the milling machine manufacturer or the manufacturer of the tools used. Before starting the milling process, ensure that the tools are fitted correctly and securely. Perform machining according to the parameters recommended by the respective manufacturer of the milling machine. After completion of the milling process, the framework must be removed from the blank. A diamond disc is recommended for this task. Any particles or dust remaining on the surface must be removed with the aid of compressed air.

### Shading the frameworks (optional)

The basic color of the framework can be matched to the desired tooth shade using the commercially available coloring solutions approved for zirconium oxide. An infrared lamp should be used for drying. In this context, the instructions for use of the manufacturer of the coloring liquids must be observed.

### 12.3 Sintering

Only furnaces designed for sintering these dental ceramics may be used for zirconium oxide restorations. To avoid impairing the sintering behavior, it must be ensured that no other ceramic materials are in the furnace. If other ceramic materials have previously been treated in the furnace, an empty firing at maximum furnace temperature is recommended. Only then should you begin sintering the zirconia ceramic. Use the appropriate sintering program for the sintering process (see following pages). After a temperature of 800 °C / 1472 °F has been reached during the cooling process, continue cooling to room temperature with the furnace closed. The framework can then be removed.

| Sintering program  |                  |                       |                        |                    |
|--|------------------|-----------------------|------------------------|--------------------|
| Level  | Temperature (°C) | Heating rate (°C/min) | Final temperature (°C) | Holding time (min) |
| <b>ceraMotion® Z Hybrid</b>                                |                  |                       |                        |                    |
| 1  | 30               | 8                     | 1150                   | 30                 |
| 2  | 1150             | 2                     | 1300                   | –                  |
| 3  | 1300             | 4                     | 1480                   | 120                |
| 4  | 1480             | -8                    | 800                    | –                  |
| <b>ceraMotion® Z HT Multishade, ceraMotion® Z HT Shade</b> |                  |                       |                        |                    |
| 1  | 30               | 8                     | 1200                   | 30                 |
| 2  | 1200             | 2                     | 1300                   | –                  |
| 3  | 1300             | 4                     | 1530                   | 120                |
| 4  | 1530             | -8                    | 800                    | –                  |
| <b>ceraMotion® Z HT Cubic Multishade</b>                   |                  |                       |                        |                    |
| 1  | 30               | 8                     | 1150                   | 30                 |
| 2  | 1150             | 2                     | 1300                   | –                  |
| 3  | 1300             | 4                     | 1450                   | 120                |
| 4  | 1450             | -8                    | 800                    | –                  |
| <b>ceraMotion® Z HT White, ceraMotion® Z White</b>         |                  |                       |                        |                    |
| 1  | 30               | 2.5                   | 500                    | –                  |
| 2  | 500              | 8                     | 1200                   | 30                 |
| 3  | 1200             | 2                     | 1300                   | –                  |
| 4  | 1300             | 4                     | 1530                   | 120                |
| 5  | 1530             | -8                    | 800                    | –                  |

### 12.4 Control and post processing

Check the seat and margin fit of the framework on the working model. If reworking is required, this should be performed by grinding and polishing. The restoration should be polished prior to clinical application to reduce antagonist abrasion.

Before fixing the restoration, it is recommended to clean the inner surfaces and blast them with Al<sub>2</sub>O<sub>3</sub> with a grain size of 50 µm at 2.0 bar.

### **12.5 Veneering and glazing**

The framework can be veneered with commercially available veneering ceramics specifically designed for zirconium oxide. The framework can be finished with a standard commercial glaze material suitable for zirconium oxide taking the manufacturer's instructions for use into consideration.

### **13. Shade in accordance with the manufacturer's shade guide**

ceraMotion® Z is based on the shade systems VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® (VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co KG, 79704 Bad Säckingen, Germany).

### **14. Information on combining with other products**

Only veneering ceramics with a coefficient of thermal expansion within the range specified for zirconium oxide may be used.

For best results, we recommend all materials suitable for zirconium oxide from the ceraMotion® product range: ceraMotion® Zr (classic veneering ceramic for zirconium oxide) and/or ceraMotion® One Touch (paste ceramic specifically developed for monolithic and cut-back solutions).

### **15. General information**

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident are informed accordingly. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaurum.com](http://www.dentaurum.com). For more information on our products, please visit the internet at [www.dentaurum.com](http://www.dentaurum.com). You will find these instructions for use and information on the current status at [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu). In case of questions or ideas, please contact your local representative.

### **16. Explanation of symbols used on the label**

Refer to the label. Additional information can be found on the internet at [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu) (Explanation of symbols used on the label REF 989-313-00).

### **17. Quality information**

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on processing on site by the user.

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