# rematitan® Rollendraht rematitan® Stangendraht

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# rematitan® wire on coils rematitan® straight wire

**Filler material for welding of titanium materials**

**DIN EN ISO 28319 / DIN EN ISO 5832-2**

## Composition (% by mass)

Ti grade 1 – 99,5% Ti (see table)

Ti grade 4 – 99,0% Ti (see table)

Additional elements < 1%: Fe, O

## Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaurum.com](http://www.dentaurum.com), even though you may frequently use the same product.

## 1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

## 2. General product description

rematitan® wire on coils and rematitan® straight wires are used as filler material for the laser welding technique.

## 3. Area of application

rematitan® wire on coils and rematitan® straight wire are filler materials for welding titanium and titanium alloys (Ti grade 5). They are particularly suitable for welding titanium materials which are intended for the production of crowns and bridges as well as for metal ceramics and also for the cast partial denture technique.

In addition, rematitan® straight wire in grade 4 is suitable for swinglocks for anchoring implant-supported superstructures.

## 4. Processing

The point of contact between the objects to be welded must always be clean. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine.

Welding may only be carried out in a shielding gas atmosphere (preferably argon) and is a prerequisite for good welding results.

## rematitan® wire on coils

Cross-section	Length	Material	REF
ø 0.40 mm	2 m	Ti grade 1	528-039-50
ø 0.70 mm	2 m	Ti grade 1	528-040-50

## rematitan® straight wire

Cross-section	Length	Material	REF
ø 1.00 mm	100 mm	Ti grade 1	528-041-00
ø 1.20 mm	100 mm	Ti grade 1	528-042-00
0.50 mm x 1.5 mm	100 mm	Ti grade 1	528-043-00
ø 1.50 mm	50 mm	Ti grade 4	528-050-00

## 5. Contraindications and adverse reactions

Signs of intolerance to titanium and titanium alloys are extremely rare if the manufacturer's instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Please adhere to warnings and safety data sheets.

## 6. Warnings and precautions to be taken

Metal chips and metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly also burs. Processed material may therefore be hot and/or have sharp edges. Due to the oxygen affinity of titanium, there is a risk of combustion and deflagration, particularly in the hot state. Any type of ignition source shall be kept at a distance. Metal fire powder extinguishers should always be positioned in close proximity to extinguish a fire if necessary. They should be in working order. The workpiece must be cooled during the mechanical machining in the presence of oxygen. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes.

**Please observe the relevant safety data sheet (see [www.dentaurum.com](http://www.dentaurum.com)).**

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The material should not be ground or polished in the patient's mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

## 7. Reuse

Wire that has already been processed must not be re-used.

## 8. Storage conditions

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of titanium and titanium alloys.

## 9. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

## 10. Quality


Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

## 11. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaurum.com](http://www.dentaurum.com).

These Instructions for use contain the main steps for processing titanium and titanium alloys. For more information on our products and services, visit [www.dentaurum.com](http://www.dentaurum.com). You will find these Instructions for use and information on the current status at [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu).

## 12. Explanation of symbols used on the label

 Please refer to the label. Additional information can be found at [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu) (Explanation of symbols REF 989-313-00).

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Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germania