NiCr-Schweißdraht

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NiCr welding wire

Filler material for welding of nickel-based alloys, DIN EN ISO 28319

Composition (% by mass)

Ni	Cr	Mo	Nb
63.0	22.0	9.0	3.5

Other elements < 1%: Al, Co, Fe, Ti

This product contains nickel and cobalt.

Beryllium, lead and cadmium free.

Dear customer,

Thank you for choosing a quality product from Dentaaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

2. General product description

NiCr welding wire is used as a filler material in conjunction for the laser welding technique.

3. Area of application

NiCr welding wire is suitable as a filler material for welding nickel-based alloys, for the production of crowns and bridges which are intended for the production of crowns and bridges and for metal ceramics..

4. Procedure

Suitable methods, such as laser welding, produce a mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Instructions for the preparation of the welding process and the process itself can be found in the brochure „Prosthetics. Laser welding examples“ (989-818-12) on the Dentaaurum website.

5. Contraindications and adverse reactions

Signs of intolerance to NiCr welding wire are extremely rare if the manufacturer's Instructions for use and processing instructions are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Individual alloy components may have a carcinogenic effect. Warnings and safety data sheets must be observed. In the State of California elements of this alloy are known to be carcinogenic.

6. Warnings and precautions to be taken

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burrs. Processed material may therefore be hot and/or have sharp edges. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

Please observe the relevant safety data sheet (see www.dentaaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

7. Reuse

Welding wire that has already been processed must not be re-used.

8. Storage conditions

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of cobalt-based alloys.

9. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

10. Quality

Dentaaurum ensures faultless quality of its products. These recommendations are based upon Dentaaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaaurum, have no influence on the processing on site.

11. General information

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaaurum.com. These Instructions for use contain the main steps for processing NiCr welding wire. Further details can be found in the processing instructions: Prosthetics, Laser welding examples (REF 989-818-12) and Orthodontics, Laser welding examples (REF 989-629-20). For more information on our products and services, please visit www.dentaaurum.com. You will find these Instructions for use and information on the current status at www.dentaaurum.com/ifu.

12. Explanation of symbols used on the label

⚠ Please refer to the label. Additional information can be found at www.dentaaurum.com/ifu (Explanation of symbols REF 989-313-00).

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Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germania