DE

PL CoCr rod

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EN

PL CoCr rod

CoCrMo - Dental alloy based on cobalt for metal/ceramic bonding, type 5 according to DIN EN ISO 22674 / DIN EN ISO 9693 for milling

Co	Cr	Mo
64.0	28.0	6.0

Other elements < 1%: Mn, Si, Fe, N
This product contains cobalt.

Nickel, beryllium, lead, cadmium-free.

Technical data

Yield strength*	R _{p0,2}	1.150 MPa
Tensile strength*	R ^m	1.350 MPa
Hardness*	H	420
Elongation at rupture*	A ₅	20%
Modulus of elasticity*	E	230 GPa
Density	ρ	8.3 g/cm ³
Solidus temperature	T _s	1.380 °C / 2.516 °F
Liquidus temperature	T _l	1.420 °C / 2.588 °F
CTE (25 °C - 500 °C / 77 - 932 °F)	α	14.3 x 10 ⁻⁶ K ⁻¹

* * Variations of ± 10 % are possible depending on the batch.

Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG I Turnstraße 31 I 75228 Ispringen I Germany

2. General product description

Non-precious alloy for the manufacture of medical devices by means of the subtractive manufacturing process.

3. Area of application

PL CoCr rods for subtractive manufacturing are cobalt-based dental alloys. They are particularly suitable for the production of constructions to be produced by milling in the form of crowns, abutments and bars, as well as for metal ceramics.

4. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction must be in line with dental regulations in order to fulfil mechanical requirements.

5. Milling

The machining can be carried out with most commercially available CAD/CAM milling machines or hand-guided copy milling machines. Suitable tools should be selected according to the machine manufacturer’s operating instructions. Before starting the milling process, ensure that the tools are clamped correctly and securely. Machining should be carried out with the milling parameters recommended by the machine manufacturer.

6. Post treatment

After completion of the milling operation, the framework must be removed from the blank. Support residues can be removed by means of finely-toothed hard metal burs. Any particles or dust remaining on the surface must be removed with compressed air.

7. Finishing

The finishing of the framework and the surfaces should be carried out by means of finely-toothed hard metal burs or diamond grinding points.

8. Welding

Suitable welding processes such as laser welding and TIG pulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler material for this alloy, e.g. CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

9. Soldering

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting temperature of the alloy used., e.g. CoCrMo-Sold 1 (REF 102-306-00) in conjunction with rema® Flux 1 (REF 102-304-00).

10. Ceramic veneering

Non-precious cobalt-based bonding alloys can be veneered with all suitable ceramic masses, e.g. ceraMotion® Me (Dentaurum), as long as the CTE is suitable for the alloy. The highest recommended firing temperature is 980 °C / 1796 °F. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

11. Contraindications and adverse reactions

Signs of intolerance to non-precious cobalt-based alloys are extremely rare if the manufacturer’s Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California elements of this alloy are known to be carcinogenic.

12. Warnings and precautions

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly also burs. Processed material may therefore be hot and/or have sharp edges. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient’s mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

13. Reuse

Fully machined rods must not be further processed for further use as medical devices. In principle, the use of new material for the production of high-quality dental prosthetics is recommended, so that clear batch traceability is also possible.

14. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of cobalt-based alloys.

15. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

16. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

17. Additional tips

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

18. Explanation of symbols used on the label

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).










IN	Gebrauchsanweisung I Instructions for use I Mode d’emploi I Modo de empleo I Modalità d’uso
----	---

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