DE

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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG
| Turnstr. 31
| 75228 Ispringen
| Germany

2. General device description, intended purpose and characteristics, technical description

The tiologic® ST implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® ST implants. Only original components of the tiologic® ST implant system should be combined in accordance with the Instructions for use. The abutment is selected by the dentist/dental technician to suit the individual patient’s needs.

The diameter of the implant and the gingival thickness determine the choice of abutment. Selection aids for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

The tiologic® ST CAD/CAM components are used for the transfer and manufacture of digital prosthetic solutions such as single-tooth and bridge restorations on tiologic® ST implants.

4. Medical indications/contraindications

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- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants ø 3.3 mm

- Cast-on technique

5. Information on use, processing, fitting and insertion

The tiologic® ST CAD/CAM components are used for the transfer and manufacture of digital prosthetic solutions such as single-tooth and bridge restorations on tiologic® ST implants.

tiologic® ST CAD/CAM data sets for 3shape™, dental wings and exocad can be downloaded free of charge from www.dentaurum.com/digital. They include all the information needed for all types of restorations as a complete package, i.e. for printed models, hybrid abutments, one-piece abutments, bridges and bars.

5.1 Scan abutment titanium tiologic® ST

The surface of the scan abutment titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. The long cylinder at the connection point of the scan abutment titanium serves for easy integration into the (laboratory) implant and for exact axial detection of the position. The enclosed red retaining screw has a shortened thread, which only grips in the (laboratory) implant when the scan abutment has been inserted into the connection point in the correct position. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer

5.2. Laboratory implant tiologic® ST for printed and cast models

The laboratory implant can be used in both analog and digital workflow. The cavities of the laboratory analogs are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory analogs with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary. In order for the counter screw supplied to fix the laboratory analog at the apical end in the printed model, the model base should be printed with a sufficient height. If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory analog in the cervical area of the printed model. The 3 retention tabs should still be visible in the printed model for monitoring purposes.

When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended series of abutments S, M or L. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The closure screw supplied with the laboratory implant enables the laboratory implant to be fixed in the end position from apical.

5.3 CAD/CAM titanium block, tiologic® ST, PreForm

By using original Dentaurum CAD/CAM PreForm titanium blocks, it is possible to create highly precise abutments that are specific to each patient, have a precise fit and can be directly screw-retained. The original PreForm titanium blocks were designed for the Medentika® PreFace® abutment holders. These ensure a perfect connection to the tiologic® ST implant interface, both for in-house production and for production in a manufacturing center. CAD/CAM PreForm titanium blocks are available for all three series of abutments S, M and L.

Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

5.4 CAD/CAM titanium base tiologic® ST

The geometry of the titanium bases is designed to ensure a reliable, aesthetic bond with the ceramic mesostructure. tiologic® ST CAD/CAM titanium bases are available in all three series of abutments S, M and L.

Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition the surfaces to be bonded. Before bonding, ensure the surfaces are dry and free from grease.Adhesives must be used according to the manufacturer’s instructions (e.g.: "PANAVIA™ F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurum GmbH & Co. KG gives no warranty for the stability of the meso/ superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the rotational security). The all-ceramic crown is completed according to the manufacturer’s instructions.

Important information

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.

Additional components

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. A new AnoTite screw should be used for the final fitting in the mouth.

Tightening torque

- Prosthetic screw, on model: manually

- AnoTite screw intra-orally: 30 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw.
The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

Contraindications.

- Primary bracing of abutments

6. Composition

- CAD/CAM titanium base: Titanium alloy (Ti6Al4V)

- CAD/CAM titanium block: Titanium alloy (Ti6Al4V)

- Scan abutment titanium: Titanium alloy (Ti6Al4V)

6 Screws:

- AnoTite screw: Titanium alloy (Ti6Al4V)
- Prosthetic screw: Titanium alloy (Ti6Al4V)
- Retaining screw: Titanium alloy (Ti6Al4V)

- Hex key:**
 - Hex key - ISO shank 1.3: (1.4197)
 - Hex key - ratchet 1.3: (1.4197)

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.

- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.

- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

- Different types of alloy in the oral cavity can lead to galvanic reactions.

- Only use tiologic® ST prosthetic components in combination with tiologic® ST implants!

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

CAD/CAM titanium blocks and titanium bases are intended for single use only. Reconditioning of a CAD/CAM titanium block or titanium base (recycling) or reuse on patients are not permitted since it cannot be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex key 1.3 and the scan abutment titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again.

Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tiologic® ST prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs.

Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview /scope of delivery, incl. REF

Please refer to the tiologic® ST product catalog for the tiologic® ST product range.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

i Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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