tioLogic®ST

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Prosthetic abutments and components replace missing teeth, restore chewing function and thus improve the quality of life. The implants of the tioLogic® implant family are used to fixate corresponding abutment components, depending on the indication. The abutment components serve to fabricate prosthetic meso-constructions, which are then restored with a prosthetic superstructure.

Abutment components					
Series of abutments*	S	M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5

*all of the abutment components listed in these Instructions for use are available in the S, M and L series of abutments.

tioLogic® ST scan abutments, titanium, are intended for digital recording of the intraoral implant position of the tioLogic® ST implants (axis inclination and position of the rotational safety lock). This can be performed either on the implant in the mouth or on the laboratory implant on the working cast. The scan abutments, titanium, are supplied with a retaining screw. The tioLogic® ST CAD/CAM scan abutments, titanium, are reusable. The retaining screw is intended for single use and is marked on the blister packaging with .

Scan abutment			
Series of abutments	S	M	L
Length in mm	12.0		

tioLogic® ST laboratory implants for printed and cast models are used to create working casts. They precisely reproduce the position of the tioLogic® ST implants (axis inclination and position of the rotational safety lock). The laboratory implants are provided together with a counter screw for model fixation. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

Laboratory implant			
Series of abutments	S	M	L
Length in mm	12.0		

tioLogic® ST CAD/CAM titanium blocks, PreForm are prefabricated blanks/components for the machined production of individual, single-piece abutments (and gingiva formers). The titanium blocks are designed for fixation in the Medentika PreFace abutment holder and subsequent computer-controlled machined production (CAM). The CAD/CAM titanium blocks are each supplied with an AnoTite screw.

Titanium block			
Series of abutments	S	M	L
Length in mm	19.9		

tioLogic® ST CAD/CAM titanium bases serve as an adhesive base for customized, implant-supported dental restorations or for the production of two-part abutments. The CAD/CAM titanium bases are each supplied with an AnoTite screw. In addition, modeling sleeves are available specifically for CAD/CAM titanium bases.

Titanium base			
Series of abutments	S	M	L
Length in mm	3.5		

The modeling sleeves for CAD/CAM titanium bases are intended for the production of mesostructures and crown frameworks that are waxed up with a modeling wax or acrylic, or are ground with a milling cutter for acrylic to achieve the final shape. The wax-up is then transferred to a framework using the casting technique and bonded to the CAD/CAM titanium base. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

The tioLogic® ST implant system contains specifically coordinated abutments and components, including instruments, for the insertion and prosthetic restoration of the implants.

3. Intended purpose

Impression components are intended for three-dimensional recording and transfer of the intraoral position of the implants.

tioLogic® ST laboratory implants for printed and cast models are designed to represent the position of the implant (axial inclination and position of the rotational safety lock) in a working cast.

Abutments are intended to connect the implants with the prosthetic restoration.

The modeling sleeve for CAD/CAM titanium bases is designed to be processed as wax-ups for the production of cast mesostructures and crown frameworks on CAD/CAM titanium bases.

4. Intended users

The tioLogic® ST product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

5.1.1 tioLogic® ST scan abutments

- Localization of the exact implant position in the intra-oral scan
- Localization of the exact implant position in the model scan

5.1.2 tioLogic® ST laboratory implants for printed and cast models

- precise reproduction of positions when printing models
- precise reproduction of positions for cast models (dental stone)

5.1.3 tioLogic® ST CAD/CAM titanium blocks

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- For the production of customized gingiva formers
- Adhesive bond with ceramic mesostructures
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

5.1.4 tioLogic® ST CAD/CAM titanium bases

- For the production of individual crown and bridge techniques as a connecting piece between the implant and the prosthetic abutment (cemented, individually screw-retained)
- Adhesive bond with ceramic mesostructures
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper lateral incisors and lower incisors

5.1.5 Modeling sleeves for CAD/CAM titanium bases

- Only use in conjunction with a CAD/CAM titanium base
- Wax-up for mesostructures and crown frameworks

5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

5.2.1 tioLogic® ST scan abutments

- Use as restoration

5.2.2 tioLogic® ST laboratory implants for printed and cast models

- Use of laboratory implants with defective or damaged connector geometry
- Use on patients

5.2.3 tioLogic® ST CAD/CAM titanium blocks

- Primary bracing of titanium blocks
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper and lower posterior teeth

5.2.4 tioLogic® ST CAD/CAM titanium bases

- Primary bracing of titanium bases
- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants \varnothing 3.3 mm
- Cast-on technique
- Direct laser welding
- Restorations with a length ratio to the length of the implant which exceeds 1:1.25
- Single crowns with a cantilever bridge unit
- Abutment on implant \varnothing 3.3 mm only for replacement in the area of the upper and lower posterior teeth

5.2.5 Modeling sleeves for CAD/CAM titanium bases

- Use on patients
- Use on other abutments than on CAD/CAM titanium bases
- Use as direct restoration (without cast)

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® ST Instructions for use.

We strongly recommend that, prior to using Dentaureum products, dental professionals attend a Dentaureum training course to familiarize themselves with the Dentaureum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaureum recommends full clinical, radiological, photographic and statistical documentation.

The CAD/CAM abutments and components of the tioLogic® ST implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

Important:

The implantologist should ensure the products cannot be aspirated during intra-oral use.

The CAD/CAM data sets for 3shape™, dental wings and exocad can be downloaded free of charge from www.dentaureum.de/cadcam. The data sets include all the information needed for all types of restorations as a complete package, i.e. for printed implant models, hybrid abutments, one-piece abutments, as well as bridges and bars.

7.1 tioLogic® ST CAD/CAM scan abutments, titanium, platform

The tioLogic® ST scan abutments, titanium, are available in the S, M or L series of abutments. Place the scan abutment intra-orally on the tioLogic® ST implant or extra-

orally on the laboratory implant in the master model (pay attention to the rotational security) and secure it with the hex key AF 1.3 and the enclosed red fixing screw. The long cylinder at the connection point of the scan abutment, titanium, serves for an easy integration into the (laboratory) implant and for exact axial detection of the position. The enclosed red retaining screw features a shortened thread, which only grips in the (laboratory) implant when the scan abutment has been inserted into the connection point in the correct position. The surface of the scan abutments titanium is optimized for digital capture, both intra-orally as well as on the model without scanning spray. The scanning and matching process of the scan data within the respective CAD software is performed according to the specifications of the individual software manufacturer. The scan abutments, titanium, are intended for multiple use while maintaining the integrity of the scan body and the connection to the (laboratory) implant. The retaining screw is intended for single use and is marked on the blister packaging with .

Tightening torque

Retaining screw with shortened thread on the model / in the mouth	15 Ncm
---	--------


7.2 tioLogic® ST laboratory implants for printed and cast models

The laboratory implant can be used in both analog and digital workflows. The cavities of the laboratory implant are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory implant with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary.

To ensure that the supplied counter screw can fixate the laboratory implant at the apical end in the printed model, the model base should be printed with a sufficient height **1**.

If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory implant in the neck area of the printed model **2**. The 3 retention tabs should still be visible in the printed model for monitoring purposes.

When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended S, M or L series of abutments or the bridge, bar orAngleFix. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The counter screw **3** supplied with the laboratory implant can be used to fixate the laboratory implant in the end position from apical.

Laboratory implants for printed and cast models are intended for multiple uses in the laboratory, provided that the outer contour and the connection to the impression and abutment components remain intact.


Tightening torque

Counter screw on the model	15 Ncm
----------------------------	--------

7.3 tioLogic® ST CAD/CAM titanium blocks, PreForm

By using original Dentaureum CAD/CAM PreForm titanium blocks, it is possible to create highly precise abutments for the prosthetic restoration that are specific to each patient, have a precise fit and can be directly screw-retained. The titanium blocks were designed for the Medentika® PreFace® abutment holders. These ensure a perfect connection to the tioLogic® ST implant interface, both for in-house production and for production in a manufacturing center.

CAD/CAM titanium block M. PreForm

CAD/CAM PreForm titanium blocks are available for all three series of abutments S, M and L. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.


Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.4 tioLogic® ST CAD/CAM titanium bases

The geometry of the tioLogic® ST CAD/CAM titanium base is designed to ensure a reliable, aesthetic bond with the ceramic mesostructure. tioLogic® ST CAD/CAM titanium bases are available in all three series of abutments S, M and L. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

The angulated screw aperture offers the optimum solution for the complete segment that is occlusally screw-retained. The screw aperture can be inclined at an angle of 20° to the implant axis, allowing discrete emergence in the palatal area, particularly in the anterior region.

Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm / 2 bar). It is also advisable to use the silanization method to condition the surfaces to be bonded. Before

bonding, ensure the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer's instructions (e.g.: "PANAVIA™ F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurem GmbH & Co. KG gives no warranty for the stability of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the rotational safety lock). The all-ceramic crown is completed according to the manufacturer's instructions.

Important information:

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- The height and diameter of the area to be bonded must not be reduced

Tightening torque	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm
AnoTite screw for angled screw aperture on the model	25 Ncm
AnoTite screw for angled screw aperture in the mouth	25 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.5 Modeling sleeve for CAD/CAM titanium bases

The modeling sleeve can be ground with a milling cutter for acrylics and shaped into the desired form of the mesostructure or the crown framework with standard burn-out modeling waxes or acrylics.

The wax-up is then transferred to a framework using the casting technique and bonded to the CAD/CAM titanium base. For additional information on bonding, see 7.4 CAD/CAM titanium base.

7.6 Additional components

The prosthetic screw (not included in scope of delivery) serves only to help the dental technician fix the abutment on the working cast temporarily.

A new AnoTite screw should always be used for the final fitting and fixation of the abutment in the mouth at the prescribed torque.

Overview of additional components:

Torque ratchet
Hex key
Hex ball point key for angulated screw apertures
Insertion key – ratchet
PentaGrip insertion key – ISO shank
Adapter, ISO shank, hexagon/ratchet
Prosthetic screw on the model

Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

8. Technical specification

Composition

Product	Composition
tioLogic® ST CAD/CAM scan abutment, titanium	3.7165
tioLogic® ST laboratory implant, for printed and cast models	3.7165
tioLogic® ST CAD/CAM titanium blocks, PreForm	3.7165
tioLogic® ST CAD/CAM titanium bases	3.7165
Modeling sleeve for CAD/CAM titanium bases	PC
AnoTite screw	3.7165
Retaining screw	3.7165
Counter screw	3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-06), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® ST prosthetic components in combination with tioLogic® ST implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurem GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® ST implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

Prosthetic abutments and components are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties.

14. Information for multiple use products

The scan abutments, titanium, are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Scan abutments, titanium, must be replaced if they cannot be clearly identified or if their function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The tioLogic® ST prosthetic abutments and components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaurem generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Please refer to the tioLogic® ST product catalog for the tioLogic® ST product range.

17. Quality information

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaurem.com/ifu

Instructions for use	REF
Screws	989-804-98
tioLogic® ST Additional components	989-805-87
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurem products	989-804-43

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