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Haltbarkeit

24 Monate bei Raumtemperatur.

Instructions for use RK 3 and RK-DUR fine investment

Dear customer

Thank you for choosing a quality product from Dentaureum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaureum.com, even though you frequently use the same product.

1. Manufacturer

Dentaureum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. With reference to our quality standards

Dentaureum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

3. Instructions for use

Before use, shake the bottle well until the content has a uniform and creamy consistency. Apply a **thin** layer of fine investment material to the cast partial denture structure. The layer must not be too thick. Do not paint over areas of fine investment which have already dried.

Allow the surface to **dry briefly** before covering with investment material. Do not allow the fine investment layer to dry fully. Mixing without vacuum improves its flow properties. For further processing steps, follow the instructions for the investment material.

Important: To prevent the investment material from becoming viscous, the bottle must be tightly sealed after use.

The effect of RK-DUR

Crystalization of a molten metal starts from nuclei on which atoms, ions or molecules of the surrounding melt accumulate. During the course of solidification, in terms of space and time, it is characterized by two processes, nuclei formation and crystal growth.

Crystal growth is **promoted by the nuclei** – and the number of nuclei thus determines the grain size. The grain size in turn largely influences the mechanical properties of a casting. Partial castings, with a fine-grained structure, have considerably better mechanical properties, such as 0.2% yield strength, tensile strength and modulus of elasticity, than those with a coarse grain.

RK-DUR fine investment contains metal oxides as nuclei-forming agents which additionally promote the formation of crystals and cause the melt to solidify to a fine-grained casting structure. More reliability with RK-DUR due to improved 0.2% yield strength, tensile strength and modulus of elasticity.

Example

Alloy	Vickers hardness H (HV10)	Elongation at fracture A ₅	0.2% Yield strength R _{p0.2} (MPa)	Tensile strength R _m (MPa)	Modulus of elasticity E (GPa)
remanium® GM 380+ without RK-DUR	360	6.5%	640	900	220
remanium® GM 380+ with RK-DUR	338	6.5%	654	963	220

Greatly improved mechanical properties result when RK-DUR is used.

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