remanium[®]
star 

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remanium® star CL 10 - 30 µm

Powder

CoCrW dental alloy based on cobalt for the firing and partial denture technique, Type 5 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 for processing with laser melting machines

1. Composition (% by mass)

| Co | Cr | W | Si |
|------|------|-----|-----|
| 60.5 | 28.0 | 9.0 | 1.5 |

Further elements < 1%: Mn, N, Nb. This product contains cobalt. Free from nickel, beryllium, lead, cadmium and iron.

2. Technical data

| Technical data | | Heat treatment | | |
|--|------------|---|---|---|
| | | classic/speed** | eco** | sprint*** |
| Yield strength* | $R_{p0.2}$ | 800 MPa | 840 MPa | 800 MPa |
| Tensile strength* | R_m | 1.170 MPa | 1.170 MPa | 1.170 MPa |
| Hardness* | H | 395 HV 10 | 395 HV 10 | 395 HV 10 |
| Elongation at rupture* | A_5 | 11% | 10.5% | 13% |
| Modulus of elasticity* | E | 230 GPa | 230 GPa | 230 GPa |
| Density | ρ | 8.6 g/ cm ³ | 8.6 g/ cm ³ | 8.6 g/ cm ³ |
| Solidus temperature | T_s | 1.320 °C / 2921 °F | 1.320 °C / 2921 °F | 1.320 °C / 2921 °F |
| Liquidus temperature | T_L | 1.420 °C / 3002 °F | 1.420 °C / 3002 °F | 1.420 °C / 3002 °F |
| CTE (25 °C - 500 °C / 77 °F - 932 °F) α | | 14.4 x 10 ⁻⁶ K ⁻¹ | 14.4 x 10 ⁻⁶ K ⁻¹ | 14.4 x 10 ⁻⁶ K ⁻¹ |

The technical data represent typical values that are achieved when using remanium® star CL with the heat treatment recommended in these Instructions for use.

* Variations of $\pm 10\%$ are possible depending on the batch.

** Typical values when using remanium® star CL via Concept Laser Mlab cusing R 100 and a layer thickness of 25 μm .

*** Typical values when using remanium® star CL via TruPrint 2000 and a layer thickness of 20 μm .

3. Intended purpose

Non-precious alloy for the manufacture of medical devices by means of an additive powder bed fusion process, e.g. selective laser melting.

4. Indication

remanium® powders for additive manufacturing are dental alloys based on cobalt. remanium® star CL is particularly suitable for the manufacture of fixed and removable prosthetic restorations, appliances and metal-ceramic frameworks.

5. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction is to be in line with dental regulations in order to fulfill the necessary mechanical requirements. The finished product must, however, have a material thickness of at least 0.3 mm.

6. Processing in the powder bed system

Before processing in the powder bed system, ensure that the system is suitable for processing the material. Product-optimized process parameters are provided for the following machine manufacturers: ZoneLab GmbH, ATLIX S.r.l, Concept Laser GmbH, Guangdong Hanbang 3D Tech Co., Ltd. (HBD), imes-icore GmbH, Nanjing Chamliion Laser Technology

Co., Ltd, Nanjing Profeta Intelligent Technology Co., Ltd, One Click Metal GmbH, Trumpf Laser- und Systemtechnik GmbH.

Further information can be found in our remanium® brochure "Processing instructions. Non-precious alloy for laser melting or milling" (REF 989-814-00). If you would like to work with a different supplier, please contact your machine manufacturer in order to clarify this or contact us. We will be pleased to help you check the system requirements.

7. Post treatment

Powder residues on the components manufactured by the additive manufacturing process must be removed with due care, taking safety precautions into consideration.

8. Stress relief heat treatment

The process outlined below describes stress relief heat treatment for products manufactured by additive manufacturing technology. The process ensures that the material properties are as indicated and guarantee that the components are of high quality and lasting durability. If the heat treatment program is altered, results may differ. This, in turn, may void the guarantee. Stress relief heat treatment optimizes the material structure of the CoCrW alloy. Inner stress is relieved and the accuracy of fit of the dental component is improved. We recommend that heat treatment of the components is carried out under a protective atmosphere on a build plate in a suitable furnace. When using a gas box, observe the gas flow rates recommended by the manufacturer. Argon or nitrogen are suitable shielding gases. Heat treatment under ambient air causes a thicker reaction layer on the surface of the components which may negatively impact the optics and possibly restrict the function of the products. The layer can, however, be removed by subsequent ablation.

Heat treatment classic:

1. Preheat the furnace to 600 °C / 1112 °F and then place the build plate with the supported components into the furnace.
2. Heat the furnace at maximum power to 1150 °C / 2102 °F.
3. Hold the temperature at 1150 °C / 2102 °F for 60 minutes, then switch off the heating.
4. Remove the build plate when the temperature has cooled down to 600 °C / 1112 °F.
5. Processing of the components can continue once they have cooled down.

Heat treatment speed:

1. Preheat the furnace to 1050 °C / 1922 °F and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 1050 °C / 1922 °F for 60 minutes.
3. Then, remove the build plate and let it cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

Heat treatment eco:

1. Preheat the furnace to 850 °C / 1562 °F and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 850 °C / 1562 °F for 40 minutes.
3. Then, remove the build plate and let it cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

Heat treatment sprint (TruPrint 1000 and 2000):

1. Preheat the furnace to 1100 °C / 1922 °F and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 1100 °C / 2012 °F for 45 minutes.
3. Then, remove the build plate and let it cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

9. Post treatment

The components can be removed from the build plate using a separating disc or band saw. Reaction layers and support residues can be removed and surfaces smoothed using a finely-toothed hard metal bur. The surface can be conditioned either by polishing manually or mechanically. When using a polishing unit with a dry suspension, the surfaces are processed in preparation using blasting abrasive (REF 128-016-00) and, if required, using finely-toothed hard metal burs.

10. Powder processing

Excess powder can, for example, be passed through a vibratory sieve shaker from the AS 200 series (Retsch) with a mesh of approx. 40 µm. For compact screening machines with ultrasound stimulation and outlet tube, a screening mesh width of 63 µm is recommended, e.g. CSM 500 (Assonic). The separated powder can undergo visual inspection and, if considered in good order, be further processed in the powder bed system. Alternatively it can be stored, dry and under seal, in its container.

11. Welding

Suitable welding processes such as laser welding and TIG pulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Suitable welding wires are available as filler materials for all remanium® alloys. Examples: CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

12. Soldering

Soldering should be avoided if possible to prevent a mixture of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting range of the alloy used. Example: CoCrMo-Sold 1 (REF 102-306-00) in connection with rema® Flux 1 (REF 102-304-00).

13. Ceramic veneering

remanium® bonding alloys can be veneered with all suitable ceramic masses, e.g. ceraMotion® Me (Dentaurum), as long as the CTE is suitable for the alloy. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

14. Contraindications and adverse reactions

Signs of intolerance to non-precious remanium® alloys are extremely rare if the manufacturer's Instructions for use are adhered to. In the case of a proven allergy to any component of the alloy, it must not be used for safety reasons. There have been individual reports of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California elements of this alloy are considered to be carcinogenic.

15. Warnings and precautions

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burrs. Processed material may therefore be hot and/or have sharp edges. It is possible that dust may be generated when the package is opened or the remanium® powder is transferred. Metal dust is combustible. Keep away from ignition

sources. Always keep metal fire extinguishing powder within easy reach and ready for use. Processing at temperatures above the solidus temperature may generate vapors. Equip the workplace with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and especially respiratory protection equipment are worn. Should there be contact with the skin, wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Consult a doctor if you are exposed to individual elements.

Please observe the safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

16. Reuse

The alloy may change its composition and properties if reused. In principle, we recommend using new material to manufacture high-quality dental restorations.

17. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the remanium® alloys. remanium® powder should be stored in a dry place in its original packaging with the lid closed tightly.

18. Disposal

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

19. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, they should inform the manufacturer and the competent authority in the country in which the user and/or the patient are resident. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and at www.dentaurum.com.

These Instructions for use contain the main steps for processing the non-precious remanium® alloys. Further information can be found in our remanium® star brochure "Processing instructions. Non-precious alloy for laser melting or milling", REF 989-814-00.

For more information on our products, please visit www.dentaurum.com. You will find these Instructions for use and information on the current status at www.dentaurum.com/ifu.

Our dental technical team is available to answer any questions that may arise. Dental Technical Hotline Tel No: +49-7231-803-410.

20. Explanation of symbols used on the label

Refer to the label. Additional information can be found on the internet at www.dentaurum.com/ifu (Explanation of symbols used on the label REF 989-313-00).

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