Gebrauchsanweisung | Instructions for use | Mode d'emploi |
Modo de empleo | Modalità d'uso | Instrukcja używania

cera**motion**[®]
LiSi

Lithium-Disilikat Presskeramik

Lithium disilicate press ceramic

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Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germany

Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

1. Composition (% by mass)

Shaded lithium disilicate glass-ceramic, made of

59%-70%	SiO ₂
10%-19%	Li ₂ O
1%-3%	Al ₂ O ₃
0-3%	K ₂ O
0-6%	P ₂ O ₅
0-3%	ZnO
0-8%	ZrO ₂
0-2%	Pigments from inorganic oxide-ceramics

2. Technical data

The technical data was determined on the basis of DIN EN ISO 6872.

	Type	Class	Coefficient of thermal expansion α (25 °C to 500 °C / 77 - 932 °F) in $1 \cdot 10^{-6} \text{ K}^{-1}$	Chemical solubility <i>L</i> in $\mu\text{g}/\text{cm}^2$	Flexural strength σ in MPa
ceraMotion® LiSi	II	3	10.2	20	400

3. Intended purpose

Dental ceramic semi-finished products are intended for the production of dental restorations to restore both the chewing function and aesthetics.

4. Indication

Press ceramic for the production of veneers, inlays, onlays, partial crowns, anterior and posterior tooth crowns, three-unit anterior bridges, three-unit bridges in the premolar area whereby the second premolar acts as the terminal abutment.

5. Contraindications and adverse reactions

There is no evidence of incompatibility with these materials when manufacturing or processing in accordance with the Instructions for use. The product should not be used if there is a known allergic reaction to one or more of the material components.

Ceramic restorations are not recommended or are only recommended to a limited extent in the case of bruxism or other parafunctions. Consultation with the attending dentist is strongly recommended.

The following are generally contraindicated: inadequate oral hygiene, extreme dystrophy of the teeth, non-suitable or self-adhesive cementation, use or combination of the product with materials which are not included in the indications. Production of inlay bridges, cantilever bridges and adhesive bridges is not permitted.

6. Patient population

Patients with defects in tooth structure for whom a restoration is necessary and indicated.

7. Intended users

Processing must be carried out only by dental professionals with knowledge of the production of dental restorations. These Instructions for use alone are insufficient for personnel without adequate training to professionally and correctly process prosthetic restorations and related components. Non-observance of this advice may lead to failure of the restoration and other complications.

8. Warnings and precautions to be taken

Dust may also be generated during mechanical processing. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. Vapors may also form if liquid substances are used during processing. Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract.

Equip the workplace with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Consult a doctor if you are exposed to critical substances.

Adhere to warnings and safety data sheets.

If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary. The prosthetic restoration should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers, or children.

9. Reuse

Finished restorations which have been produced using semi-finished products are intended for single use only.

10. Storage and/or handling conditions

The product properties are not affected by normal changes in environmental conditions (such as temperature, pressure, humidity or light). Close the original packaging after use to prevent contamination. In general, the product should be kept dry and at room temperature in the tightly closed original packaging.

11. Disposal

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

12. Using the products

12.1 Preparation

In preparation for using ceramics, the tooth structure should be reduced to a reproduction of the tooth shape that is smaller and less detailed. Sharp edges and interior angles should be avoided. The design of the restoration must be adapted to the static and dynamic occlusion. Observe the minimum wall thicknesses. These depend on the application or the restoration in question and the ensuing processing technique (Figs. 1 to 4 and Table 1).


Fig. 1: Preparation instructions, crown, with minimum wall thicknesses


Fig. 2: Preparation instructions, crown, with minimum wall thicknesses


Fig. 3: Preparation instructions, inlay/onlay with minimum wall thicknesses


Fig. 4: Preparation instructions, veneer (thin) with minimum wall thicknesses

Table 1 Minimum wall thicknesses (in mm)

Design	Veneer		Inlay	Onlay	Crown		Three-unit bridge	
					anterior	posterior	anterior area	posterior region*
Staining technique	cervical	0.3 to 0.6	1.0	1.0	1.0	1.0	1.2	1.5
	incisal / occlusal	0.4 to 0.7	1.0	1.0	1.0	1.0	1.5	1.5
Cut-back technique	cervical	0.6	—	1.5	1.0	1.0	1.2	1.5
	incisal / occlusal	0.4	—	0.8	0.5	0.5	0.8	0.8
Layering technique*	cervical	—	—	—	0.5	0.5	0.6	0.8
	incisal / occlusal	—	—	—	0.5	0.5	0.8	0.8
Cross-section connector (mm ²)		—	—	—	—	—	16**	16**

*The overall wall thickness of the restoration must be made of at least 50% ceraMotion® LiSi. A bridge in the posterior region is limited to the second premolar.

** Height of the connector must always be the same or more than its width.

The master models are produced analogous to the models for the manufacture of metal ceramic work and inlays. The preparation margin is established. A layer of sealer is recommended to harden the surface.

In the case of crowns and veneers, the application of a suitable die spacer up to 1 mm apical to the preparation margin is recommended for a cement gap of approx. 20 µm. For inlays and onlays, apply the die spacer up to shortly before the preparation margin for a cement gap of approx. 30 µm.

12.2 Preparation

Modelling of the wax objects is carried out to suit the finishing technique in use (staining technique, micro-layering, cut back and layering technique).

12.3 Attaching sprues


Fig. 5: Sprue attachment single crown and 3-unit bridge

One ingot is sufficient for a mass of sprued models $\leq 0,75$ g. Use two ingots for a mass of sprued models up to 1.6 g.

The press sprue diameter of the investment ring base must correspond to the diameter of the press ingots and the press plunger (≥ 13.0 mm). The sprues should be attached in line with the boundary conditions outlined in Fig. 5 and in Table 2.

Table 2 Boundary conditions

Investment ring system	100 g and 200 g
Diameter pressing channel of sprue	$\varnothing$ 3.0 mm to 3.5 mm
Length pressing channel of sprue	4.0 mm – 6.0 mm, attach sprues to objects at the same height
Attaching sprues on the waxed-up object	at the most voluminous point / veneers -> incisal Inlay and partial crowns -> interproximal
Angle of sprue to the waxed-up object	in flow direction of the pressing
Angle of sprue to the investment ring base	45 ° to 60 °
Distance between press objects	min. 5.0 mm
Minimum thickness between object and investment ring	10.0 mm

12.4 Investing

Apply a thin layer of Vaseline to the ring base and the ring leveller. If using a wax surface reducing agent, remove this from the press object with oil-free compressed air. We recommend the press investment material ceraMotion® press invest. Refer to the Instructions for use (REF 989-301-03) for more detailed information. Follow the steps in Fig. 6 for investing.

Fig. 6a) Pour investment material into the casting ring with slight vibration.

Fig. 6b) Position ring leveler after sufficient filling.

Fig. 6c) Once the investment material has set, remove the casting ring.

Fig. 6d) Using a plaster knife, smooth out any unevenness.

12.5 Procedure for investing

Select a suitable press program (Table 3) and preheat the press furnace to the start temperature. When moving the casting ring, ensure that it does not cool down too much and that the ingot is quickly placed into the sprue with the marking facing upwards. Then use a cold press plunger, place the casting ring placed in the pressing furnace and start the press program. After the program has finished, remove the casting ring and allow it to cool down to room temperature.

Table 3: Pressing recommendations

Casting ring	Start-temp. °C (°F)	Speed of temperature increase (°C(°F)/min)	Final temperature (°C/°F)	Holding time (mins)	Max. pressing time (mins)	Pressure (bar)	Vacuum
100 g	800 (1472)	60 (140)	900 (1652)	20	2	2 to 3	Yes
200 g	800 (1472)	60 (140)	905 (1661)	20	2	2 to 3	Yes

Parameters are guide values and must be adjusted if necessary.

12.6 Devesting

Follow the steps in Figs. 7) – e) for correct devesting.

Fig. 7a) After cooling, mark the length of the press plunger on the casting ring.

Fig. 7b) Using a separating disc, create a predefined breaking point all around.

Fig. 7c) Break at the incision using a plaster knife.

Fig. 7d) Separated casting ring at transition point from plunger to pressing cone.

Fig. 7e) Rough devesting with Al_2O_3 at max. 4 bar pressure. Careful, precise blasting with glass beads (medium grain size) at max. 2 bar pressure.

12.7 Separating and finishing

To prevent the ceramic from overheating, it is recommended that work is carried out at a low pressing pressure and water is used to cool if necessary. Cut off the press sprue with a thin diamond disc without exerting pressure.

Use grinding tools that are suitable for ceramic for the finishing process. Carefully fit the object onto the stump. Before further processing (glazing/staining/layering), blast surface at 1 bar pressure with Al_2O_3 medium grain size, then steam-clean.

13. Shade in accordance with the manufacturer's shade guide

ceraMotion® LiSi is based on the shade systems VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER®

(VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® are registered trademarks of Vita Zahnfabrik H.Rauter GmbH & Co KG, 79704 Bad Säckingen, Germany).

14. Information on combining with other products

The dental ceramic may only be combined with other products if recommended by Dentaaurum. Veneering ceramics must be suitable for the material composition and the CTE of ceraMotion® LiSi (e.g. ceraMotion® One Touch and ceraMotion® Zr).

15. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaaurum.com.

For more information on our products and services, visit www.dentaaurum.com. You will find these Instructions for use and information on the current status at www.dentaaurum.com/ifu.

In case of questions or ideas, please contact your local representative.

16. Explanation of symbols used on the label

Refer to the label. Additional information can be found at www.dentaaurum.com/ifu (Explanation of symbols REF 989-313-00).

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