cera^{motion}_{zr}[®]

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Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

1. Composition (% by mass)

Colored glass ceramic, made of

35%-60%	SiO_2
9%-19%	Al_2O_3
0-9%	K_2O
0-20%	Na_2O
0-5%	Li_2O
0-10%	B_2O_3
0-5%	CaO
0-3%	MgO
0-5%	BaO
0-21%	ZrO_2
0-8%	BaTiO_3
0-15%	pigments from inorganic oxide-ceramics

Fluids used to mix pastes, to adjust viscosity and to reinforce contrasts consist mainly of alcohols.

2. Technical data

The technical data was determined on the basis of DIN EN ISO 6872.

ceraMotion® Zr	Type	Class	Coefficient of thermal expansion α (25 °C to 500 °C/77 °F to 932 °F or T_g) in $1 \cdot 10^{-6} K^{-1}$	Glass transition temperature T_g in °C (°F)	Chemical solubility L in $\mu g/cm^2$	Flexural strength σ in MPa
Liner	I	1	9.2	565 (1049)	33	140
Dentin	I	1	9.2	530 (986)	30	125
Incisal	I	1	9.2	530 (986)	30	125
Modifier	I	1	9.2	530 (986)	30	125
Touch Up	I	1	8.8	520 (968)	20	135
Glaze	I	1	9.0	480 (896)	40	140
Stains	I	1	8.4	530 (986)	45	130
2D Paste	I	1	9.0	490 (914)	40	120
3D Paste	I	1	9.1	490 (914)	40	120

3. Intended purpose

Dental ceramic semi-finished products are intended for the production of dental restorations to restore both the chewing function and aesthetics.

4. Indication

Bonding ceramic to veneer frameworks, single crowns and bridges made of zirconium dioxide with a coefficient of thermal expansion (CTE) from $9.7 \cdot 10^{-6} K^{-1}$ to $10.6 \cdot 10^{-6} K^{-1}$ (25 °C to 500 °C/77 °F to 932 °F) and of lithium disilicate.

5. Contraindications and adverse reactions

There is no evidence of incompatibility with these materials when manufacturing or processing in accordance with the Instructions for use. The product should not be used if there is a known allergic reaction to one or more of the material components.

Ceramic restorations are not recommended or are only recommended to a limited extent in the case of bruxism or other parafunctions. Consultation with the attending dentist is strongly recommended.

The following are generally contraindicated: inadequate oral hygiene, extreme dystrophy of the teeth, non-suitable or self-adhesive cementation, use or combination of the product with materials which are not included in the indications.

6. Patient population

Patients with defects in tooth structure for whom a restoration is necessary and indicated.

7. Intended users

Processing must be carried out only by dental professionals with knowledge of the production of dental restorations. These Instructions for use alone are insufficient for personnel without adequate training to professionally and correctly process prosthetic restorations and related components. Non-observance of this advice may lead to failure of the restoration and other complications.

8. Warnings and precautions to be taken

It is possible that dust may be generated when the package is opened or the powder is transferred. Dust may also be generated during mechanical processing. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. Vapors may also form if liquid substances are used during processing. Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract.

Adhere to warnings and safety data sheets!

Equip the workplace with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Consult a doctor if you are exposed to critical substances.

If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

The prosthetic restoration should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers, or children.

9. Reuse

Finished restorations which have been produced using semi-finished products are intended for single use only.

10. Storage and/or handling conditions

The product properties are not affected by normal changes in environmental conditions (such as temperature, pressure, humidity or light). Close the original packaging after use to prevent contamination. In general, the product should be kept dry and at room temperature in the tightly closed original packaging. Do not return pastes that have been mixed to the jar.

11. Disposal

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

12. Characteristics of the products

Use for	Product	Explanation
	Bonder B	Only to be used in combination with ceraMotion® Ti. The Powder Bonder creates an optimal bonding between titanium and ceramic. The powder Bonder should be mixed with the Powder BOL Liquid. The Bonder is not fluorescent.
	Paste Bonder PB	Only to be used in combination with ceraMotion® Ti. The Paste Bonder creates an optimal bonding between titanium and ceramic. The consistency of the paste can be adjusted by adding Paste Liquid. The Bonder is not fluorescent.
	Primer CAD/CAM	Primer especially for use with CAD/CAM manufactured crowns and bridges made of CoCr alloys.
	PO Base	Especially developed for precious metal alloys to completely cover the framework, ensuring there will be no discoloration.
	Paste Opaque PO	Best bonding with low layer thickness and high opacity. The consistency of the paste can be adjusted by adding Paste Liquid.
	Paste Opaque Modifier POM	The Paste Opaque Modifier materials can be used to individually characterize the Paste Opaque. The consistency of the paste can be adjusted by adding Paste Liquid.
	Opaque O	The powder Opaque. Optimal bonding in a thin layer, dense coverage of the metal framework, high opacity. The powder Opaque is mixed using the Powder BOL Liquid.
	Opaque Modifier OM	The powder Opaque Modifier materials can be used to individualize the powder Opaque. The powder Opaque Modifier materials are mixed with the Powder BOL Liquid
	Liner L	Only use in connection with ceraMotion® Zr. For use on non-shaded zirconium oxide frameworks to ensure shade harmonization. The Liner can also be used to protect the shade on shaded zirconium oxide frameworks. The Liner is mixed with the Powder BOL Liquid.
	Liner Modifier LM	The Liner Modifier materials can be used to individualize the Liner. Mix with the Powder BOL Liquid

Use for	Product	Explanation
Ti Me Zr Lf	Base Dentin BD	The Base Dentin materials correspond to the main dentin shade. They can be used alone, mixed with other materials or as an undercoat. The Base Dentin has a higher degree of opacity than Dentin and creates the appearance of depth. Base Dentin can be used to enhance the basic tooth shade where there is little space.
Ti Me Zr	Base Dentin Modifier BDM	The Base Dentin Modifier materials have the same degree of opacity as Base Dentin. They are used to create shade effects with depth.
Ti Me Zr Lf	Dentin D	The Dentin materials A-D, in combination with the Incisal materials, are used to reproduce the shades within the Vita Classic shade guide. The Dentin materials have a natural fluorescent appearance.
Ti Me Zr Lf	Dentin Modifier Chroma DM C	The Dentin Modifier Chroma materials are used to intensify the Dentin color. Applied in a thin layer they can underline the tooth shade within the tooth body. They can create or replicate shade effects.
Ti Me Zr Lf	Dentin Modifier Fluo DM F	The Dentin Modifier Fluo materials have a higher degree of fluorescence than the Dentin materials and can create or replicate shade effects.
Ti Me Zr Lf	Incisal I	The Incisal materials, combined with the Dentin materials, are used to reproduce the shades within the Vita Classic shade guide.
Ti Me Zr Lf	Incisal Opal IO	The Incisal Opal materials show natural opalescence. Bluish effects in reflecting light, orange effects in transmitting light.
Ti Me Zr	Incisal Transpa IT	The Incisal materials Transpa have a higher degree of translucency than the (standard) Incisal materials.
Ti Me Zr Lf	Transpa T	The Transpa material allows the translucency and density to be controlled.
Ti Me Zr Lf	Incisal Modifier IM	The Incisal Modifier materials are used for the individual characterization of the incisal area.
Ti Me Zr Lf	Chroma Concept CC	The Chroma Concept materials underline the Dentin shade and can be used to create or replicate shade effects. The shade spectrum of these materials allows the reproduction of reddish, yellowish and whitish shades. This means that all natural shades can be reproduced, including Vita-3D - Master® or Ivoclar-Chromascop® shades. Specific opaque materials and an incisal material round off this assortment.
Me Zr	Value Concept IV, IM, VM, TM, TV	The Value Concept masses regulate the brightness without changing the shade. Modifiers for the cervical, dentine and incisal areas.

Use for	Product	Explanation
Ti Me Zr Lf	Shoulder SM	The Shoulder materials are divided into the shade groups A-B-C-D and allow the reproduction of the Vita Classic shade guide. The Shoulder material "white" can be mixed accordingly to produce all shade nuances from A1 to D4. Adding the Shoulder material "transparent" enables the translucence and the saturation within the shoulder area to be controlled.
Ti Me Zr Lf	Gingival G	The Gingival materials enable an aesthetical and anatomical reconstruction of the gingiva. They are not fluorescent.
Ti Me Zr	Touch Up TU	The Touch Up materials have a lower firing temperature which allows the glaze effects on the surface to be controlled. Corrections and a glaze firing can be performed in one firing sequence. They are available as Base Dentin, Dentin and Incisal materials. There is a repair opaque in addition for the Me ceramics.
Ti Me Zr Lf	Correction C	The Correction material is non-shaded, has a lower firing temperature and can be used for smaller corrections. It can be used either alone or mixed with other materials for small corrections.
Ti Me Zr	Glaze GL	The powder Glaze material increases the shine during the glaze firing.
Zr	Paste Glaze PGL/ Paste Glaze bright	The paste Glaze material increases the shine during the glaze firing.
Ti Me Zr	Body Stains BST/ Stains ST	Body Stains and Stains can be used to create shade effects. They can either be inlaid or mixed. They can be mixed with the ceramic materials up to a maximum of 10 wt %. In addition, the Stains can be used to enhance the shades on the ceramic surface.
Lf	Lf Body Stain LF BST/ Lf Stains LF ST	Body Stains and Stains can be used to create shade effects. They can either be inlaid or mixed. They can be mixed with the ceramic materials up to a maximum of 10 wt %. In addition, the Stains can be used to enhance the shades on the ceramic surface. The CTE is higher, especially for the high-expanding ceraMotion® Lf ceramic materials.
Ti Me Zr Lf	Modelling Liquid	Standard modelling liquid.
Me	Me Standard Modelling Liquid	Standard liquid for mixing ceramic materials, especially for ceraMotion® Me.
Ti Me Zr Lf	Modelling Liquid +	Alternative to Modelling Liquid and Me Standard Modelling Liquid. It has a slightly higher viscosity in comparison to the Me Standard Modelling Liquid and can be diluted with distilled water.
Ti Me Zr Lf	Modelling Liquid MV universal	Modelling liquid to increase plasticity and extend the modelling properties.

Use for	Product	Explanation
	Paste Liquid	For adjusting the viscosity of the Paste Opaque, the Paste Bonder and the Me Primer CAD/CAM.
	Powder BOL Liquid	For mixing the ceraMotion® powder Bonder, powder Liner und powder Opaque.
	Shoulder Liquid	For mixing the Shoulder materials, to increase the stability.
	Stains Liquid	For mixing the glaze material Glaze and the stains material Stains/Body Stains.
	Contrast Marker	The shaded liquids can be used to show a better contrast between the mixed materials during build-up. The viscosity is equivalent to the Modelling Liquid.

13. Using the products

13.1 Wax-up

Always keep to the layer thicknesses stipulated for veneering (Fig. 1 and Fig. 2). At the same time, pay attention to the relation in layer thicknesses between the framework and the ceramic. The framework replicates the tooth in an anatomical form reduced in size; it is important to avoid corners and edges in the framework. The veneering ceramic can only be fired at a maximum layer thickness of 2 mm.

13.2 Processing steps for build-up

Step 1 Framework preparation: The framework should be finished and conditioned in accordance with the manufacturer's instructions.

Step 2 Connecting or liner firing: Connecting firing for shaded frameworks is achieved by a thin layer of Dentin or Base Dentin in the required shade. White frameworks must be given a thin, masking layer of Liner.

Step 3 Dentin firing 1: The complete anatomical shape of the tooth is built up with dentin and incisal masses. Modifiers are used as per individual requirements. The layers should be somewhat over-contoured to compensate for shrinkage during sintering. To control shrinkage on bridges, the applied ceramic should be separated interdentally down to the framework before the first dentin firing. Carry out the first dentin firing according to the Firing Chart.

Step 4 Dentin firing 2: Carry out finishing and shape corrections after the first dentin firing. Carry out the second dentin firing according to the Firing Chart.

Step 5 Correction firing: Repeat the process for the second dentin firing until the desired result is achieved.

Step 6 Glaze firing: Grind the entire surface evenly and clean thoroughly before the glaze firing. Individual shade effects can be achieved by applying Stains/Body Stains to the surface. Apply glaze material if required.


Fig. 1: Framework design, anterior crown


Fig. 2: Framework design, posterior tooth crown

13.3 Firing Chart

The firing temperatures indicated were determined in a calibrated furnace. Use only calibrated furnaces to ensure the required product quality. Firing temperatures may deviate nonetheless depending on the furnace. It may be necessary to individually adapt the specifications given. A test firing is recommended, whereby Transpa T is mixed with Modelling Liquid. The sample is placed on a platinum sheet. The correct furnace temperature is visually assessed and determined by increasing or decreasing the temperature in 10 °C/50 °F steps. The furnace temperature is optimal if the sample is translucent and has sharp edges. If the temperature is too high, the sample has no sharp edges and appears to be very shiny. If the temperature is too low, the sample is white and milky.

Zirconium oxide

	Temperature at start °C (°F)	Drying time min	Heating rate °C (°F)/min	Vacuum start °C (°F)	Vacuum end °C (°F)	Firing temperature °C (°F)	Holding time**
Connecting Liner 1 + 2	500 (932)	4	55 (131)	500 (932)	810 (1490)	810 (1490)	2 min (under vacuum)
Connecting Base Dentin/Dentin	500 (932)	6	55 (131)	500 (932)	780 (1436)	780 (1436)	2 min (under vacuum)
Shoulder 1 + 2	500 (932)	6	55 (131)	500 (932)	780 (1436)	780 (1436)	2 min (under vacuum)
Dentin 1	500 (932)	6	55 (131)	500 (932)	750 (1382)	750 (1382)	2 min (under vacuum)
Dentin 2	500 (932)	4	55 (131)	500 (932)	750 (1382)	750 (1382)	2 min (under vacuum)
Correction material, pure	500 (932)	4	55 (131)	500 (932)	715 (1319)	715 (1319)	1 min (under vacuum)
Correction material 1:1 Dentin/Incisal**	500 (932)	4	55 (131)	500 (932)	730 (1346)	730 (1346)	1 min (under vacuum)
Glaze***	500 (932)	4	55 (131)	***	***	750 (1382)	1 min
Shine with glaze	500 (932)	6	55 (131)	500 (932)	750 (1382)	750 (1382)	1 min
Touch Up Glaze and Correction	500 (932)	6	55 (131)	500 (932)	730 (1346)	730 (1346)	1 min

* Reducing the heat rate improves the firing quality in the case of large restorations.

** Extend the hold-time in the case of large restorations to compensate for the poor thermal conductivity of ZrO₂.

*** The glaze firing can be carried out with or without vacuum.

**** The Correction material can be mixed 1:1 with Base Dentin, Dentin or Incisal.

Lithium disilicate

	Temperature at start (°C/°F)	Drying time min	Heating rate °C (°F)/min	Vacuum start °C (°F)	Vacuum end °C (°F)	Firing temperature °C (°F)	Holding time
Base Dentin Connecting	500 (932)	6	55 (131)	500 (932)	760 (1400)	760 (1400)	1 min (under vacuum)
Dentin 1	500 (932)	6	55 (131)	500 (932)	760 (1400)	760 (1400)	1 min (under vacuum)
Glaze*	500 (932)	4	55 (131)	*	*	750 (1382)	1 min
Shine with glaze	500 (932)	6	55 (131)	500 (932)	750 (1382)	750 (1382)	1 min

* The glaze firing can be carried out with or without vacuum.

Please take note of the manufacturer's instructions when creating the framework.

14. Shade in accordance with the manufacturer's shade guide

ceraMotion® Zr is based on the shade systems VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® (VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co KG, 79704 Bad Säckingen, Germany).

15. Information on combining with other products

The dental ceramic may only be combined with other products if recommended by Dentaaurum. Only frameworks may be veneered where the coefficient of thermal expansion lies within the given range for veneering ceramics.

The dental ceramic may not be used for frameworks made of lithium disilicate, zirconium oxide, titanium, titanium alloys or dental alloys with a silver content exceeding 30%.

16. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaaurum.com.

For more information on our products and services, visit www.dentaaurum.com. You will find these Instructions for use and information on the current status at www.dentaaurum.com/ifu.

In case of questions or ideas, please contact your local representative.

17. Explanation of symbols used on the label

Refer to the label. Additional information can be found at www.dentaaurum.com/ifu (Explanation of symbols REF 989-313-00).

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