
Triceram®

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Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

1. Composition (% by mass)

Colored glass ceramic, made of

35%-60%	SiO_2
9%-19%	Al_2O_3
0-9%	K_2O
0-17%	Na_2O
0-3%	Li_2O
0-10%	B_2O_3
0-5%	CaO
0-3%	MgO
0-5%	BaO
0-21%	ZrO_2
0-8%	BaTiO_3
0-15%	pigments from inorganic oxide-ceramics

Fluids used to layer ceramic masses, to adjust viscosity and to reinforce contrasts consist mainly of alcohols.

2. Technical data

The technical data was determined on the basis of DIN EN ISO 6872.

Triceram®	Type	Class	Coefficient of thermal expansion α (25 °C to 500 °C/77 °F to 932 °F or T_g) in $1 \cdot 10^{-6} K^{-1}$	Glass transition temperature T_g in °C (°F)	Chemical solubility L in $\mu g/cm^2$	Flexural strength σ in MPa
Bonder	I	1	9.4	575 (1067)	35	125
Opaque	I	1	8.9	565 (1049)	33	140
Dentin	I	1	8.5	545 (1013)	30	125
Incisal	I	1	8.5	545 (1013)	30	125
Glaze	I	1	8.3	525 (977)	40	140
Stains Universal	I	1	8.4	530 (986)	45	130

3. Intended purpose

Dental ceramic semi-finished products are intended for the production of dental restorations to restore both the chewing function and aesthetics.

4. Indication

Bonding ceramic to veneer frameworks, single crowns and bridges made of titanium and titanium alloys and made of zirconium dioxide with a coefficient of thermal expansion (CTE) from $9.6 \cdot 10^{-6} K^{-1}$ to $10.5 \cdot 10^{-6} K^{-1}$ (25 °C to 500 °C/77 °F to 932 °F).

5. Contraindications and adverse reactions

There is no evidence of incompatibility with these materials when manufacturing or processing in accordance with the Instructions for use. The product should not be used if there is a known allergic reaction to one or more of the material components.

Ceramic restorations are not recommended or are only recommended to a limited extent in the case of bruxism or other parafunctions. Consultation with the attending dentist is strongly recommended.

The following are generally contraindicated: inadequate oral hygiene, extreme dystrophy of the teeth, non-suitable or self-adhesive cementation, use or combination of the product with materials which are not included in the indications.

6. Patient population

Patients with defects in tooth structure for whom a restoration is necessary and indicated.

7. Intended users

Processing must be carried out only by dental professionals with knowledge of the production of dental restorations. These Instructions for use alone are insufficient for personnel without adequate training to professionally and correctly process prosthetic restorations and related components. Non-observance of this advice may lead to failure of the restoration and other complications.

8. Warnings and precautions to be taken

It is possible that dust may be generated when the package is opened or the powder is transferred. Dust may also be generated during mechanical processing. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. Vapors may also form if liquid substances are used during processing. Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract. Equip the workplace with a suitable suction or ventilation system.

Adhere to warnings and safety data sheets!

Equip the workplace with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Consult a doctor if you are exposed to critical substances.

If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

The prosthetic restoration should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers, or children. During dental treatment, known cross-reactions or interactions of the medical device with other medical devices or materials present in the mouth must be taken into account.

9. Reuse

Finished restorations which have been produced using semi-finished products are intended for single use only.

10. Storage and/or handling conditions

The product properties are not affected by normal changes in environmental conditions (such as temperature, pressure, humidity or light). Close the original packaging after use to prevent contamination. In general, the product should be kept dry and at room temperature in the tightly closed original packaging. Do not return pastes that have been mixed to the jar.

11. Disposal

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

12. Characteristics of the products

Device	Explanation
Bonder B	Only use in connection with Triceram®. The powder Bonder creates an optimal bonding between titanium and ceramic. The powder Bonder is mixed with the Universal B.O.L. Liquid.
Paste Bonder PB	Only use in connection with Triceram®. The Paste Bonder creates an optimal bonding between titanium and ceramic. The consistency of the paste can be adjusted by adding Universal Paste Liquid.
Opaque O	Optimal bonding in a thin layer, dense coverage of the metal framework, high opacity. The powder opaque is mixed with the Universal B.O.L. Liquid.
Opaque Modifier	The powder Opaque Modifier materials can be used to individualize the powder opaque. The powder Opaque Modifier is mixed with the Universal B.O.L. Liquid.
Fluo Liner FL	Use only for zirconium oxide frameworks to achieve an optimum bond and reliable shade coordination. Only use in connection with Triceram®. The Liner is mixed with the Universal B.O.L. Liquid.
Opaque Dentin OD	The opaque dentins correspond to the dentin shade. They can be used alone, mixed with other materials or as an undercoat. The Base Dentin has a higher degree of opacity than Dentin and creates the appearance of depth. Base Dentin can be used where there is insufficient space to enhance the basic tooth shade.
Dentin D	The Dentin materials A-D, in combination with the Incisal materials, are used to reproduce the shades within the Vita Classic shade guide. The Dentin materials have a natural fluorescent appearance.
Chroma Dentin CD	The Chroma Dentin CD materials are used to intensify the Dentin shade. Applied in a thin layer they can underline the shade intensity within the tooth body. They can create or replicate shade effects.
Incisal transparent IT	Standard incisal materials. The Incisal materials are used in combination with the Dentin materials to reproduce the shades within the Vita Classic shade guide.
Incisal opalescent IO	The opalescent incisal materials show natural opalescence. Bluish effects in reflecting light, orange effects in transmitting light.
Neutral transparent NT	The Transpa material allows the translucency and intensity to be controlled.
Transpa Effect TE	Individualization and characterization in the incisal region.

Device	Explanation
Shoulder Material SM and Modifier SM	The Shoulder Materials are divided into the shade groups A-B-C-D and allow the reproduction of the Vita Classic shade guide. The Shoulder Material "white" can be mixed accordingly to produce all shade nuances from A1 to D4. Adding the Shoulder Material "transparent" enables the translucency and the saturation within the shoulder area to be controlled.
Gingival G	The Gingival materials enable an aesthetical and anatomical reconstruction of the gingiva.
Correction Material CM	The Correction Material is non-shaded, has a lower firing temperature and can be used for smaller corrections. It can be used either alone or mixed with other materials for small corrections.
Stains Universal ST	The Stains Universal products can be used to modify the shades. They can also be mixed with the ceramic materials up to a maximum of 10 wt%. In addition, the Stains can be used to enhance the shades on the ceramic surface.
Modelling Liquid MV Universal	To increase plasticity and extend the modelling properties.
Modelling Liquid LV + Universal	Shrinking of the materials is reduced, the ceramic remains moist longer and can be processed longer.
Universal Paste Liquid	For adjusting the viscosity of the Paste Bonder.
Universal B.O.L. Liquid	For mixing the Powder Bonder, Powder Liner und Powder Opaque.
Shoulder Materials Liquid Universal	For mixing the Shoulder Materials, to increase the stability.
Stains Universal Liquid	For mixing Stains Universal.

13. Using the products

13.1 Wax-up

Always keep to the layer thicknesses stipulated for veneering (Fig. 1 and Fig. 2). At the same time, pay attention to the relation in layer thicknesses between the framework and the ceramic. The framework replicates the tooth in an anatomical form reduced in size; it is important to avoid corners and edges in the framework. The veneering ceramic can only be fired at a maximum layer thickness of 2 mm.

13.2 Processing steps for build-up

Step 1 Framework preparation: The framework should be finished and conditioned in accordance with the manufacturer's instructions. Application and firing of the Liner for zirconium oxide frameworks. Continue with Step 4.

Step 2: Application and firing of Paste and/or Powder Bonder on titanium frameworks is vital.

Step 3 Opaque firing 1 + 2: Application and firing of powder Opaque material

Apply the powder Opaque evenly to the entire framework. No wash firing is required (observe the instructions of the alloy manufacturer). Repeat the process once.

Step 4 Dentin firing 1: The complete anatomical shape of the tooth is built up with dentin and incisal masses. Modifiers are used as per individual requirements. The layers should be somewhat over-contoured to compensate for shrinkage during sintering. To control shrinkage on bridges, the applied ceramic should be separated interdentally down to the framework before the first dentin firing. Carry out the first dentin firing according to the Firing Chart.

Step 5 Dentin firing 2: Carry out finishing and shape corrections after the first dentin firing. Carry out the second dentin firing according to the Firing Chart.

Step 6 Correction firing: Repeat the process for the second dentin firing until the desired result is achieved.

Step 7 Glaze firing: Grind the entire surface evenly and clean thoroughly before the glaze firing. Individual shade effects can be achieved by applying Stains/Body Stains to the surface. Apply glaze material if required.


Fig. 1: Framework design, anterior crown


Fig. 2: Framework design, posterior tooth crown

13.3 Firing chart for Triceram on zirconium oxide

The firing temperatures indicated were determined in a calibrated furnace. Use only calibrated furnaces to ensure the required product quality. Firing temperatures may deviate nonetheless depending on the furnace. It may be necessary to individually adapt the specifications given. A test firing is recommended, whereby Transpa T is mixed with Modelling Liquid. The sample is placed on a platinum sheet. The correct furnace temperature is visually assessed and determined by increasing or decreasing the temperature in 10 °C/50 °F steps. The furnace temperature is optimal if the sample is translucent and has sharp edges. If the temperature is too high, the sample has no sharp edges and appears to be very shiny. If the temperature is too low, the sample is white and milky.

	Temperature at start °C (°F)	Drying time min	Heating rate °C (°F)/min	Vacuum start °C (°F)	Vacuum end °C (°F)	Firing temperature °C (°F)	Holding time	Cooling time mins
Liner 1 + 2	500 (932)	4	65 (149)	500 (932)	800 (1472)	800 (1472)	1 min (under vacuum)	—
Shoulder	500 (932)	6	55 (131)	500 (932)	790 (1454)	790 (1454)	1 min (under vacuum)	—
Dentin 1	500 (932)	6	55 (131)	500 (932)	760 (1400)	760 (1400)	1 min (under vacuum)	—
Dentin 2	500 (932)	4	55 (131)	500 (932)	760 (1400)	760 (1400)	1 min (under vacuum)	—
Glaze	500 (932)	2	55 (131)	*	*	760 (1400)	1 min*	—
Correction material	500 (932)	4	55 (131)	500 (932)	720 (1328)	720 (1328)	1 min (under vacuum)	—
Gingival	500 (932)	6	55 (131)	500 (932)	760 (1400)	760 (1400)	1 min (under vacuum)	—

13.4 Firing chart for titanium

	Temperature at start °C (°F)	Drying time min	Heating rate °C (°F)/min	Vacuum start °C (°F)	Vacuum end °C (°F)	Firing temperature °C (°F)	Holding time	Cooling time mins
Bonder firing (Powder Bonder)	500 (932)	4	65 (149)	500 (932)	795 (1463)	795 (1463)	1 min (under vacuum)	—
Bonder firing (Paste Bonder)	500 (932)	6	65 (149)	500 (932)	795 (1463)	795 (1463)	1 min (under vacuum)	—

Opaque firing 1+2	500 (932)	4	65 (149)	500 (932)	795 (1463)	795 (1463)	1 min (under vacuum)	—
Shoulder	500 (932)	6	55 (131)	500 (932)	785 (1445)	785 (1445)	1 min (under vacuum)	—
Dentin firing	500 (932)	6	55 (131)	500 (932)	755 (1391)	755 (1391)	1 min (under vacuum)	—
Correction firing	500 (932)	4	55 (131)	500 (932)	755 (1391)	755 (1391)	1 min (under vacuum)	—
Glaze	500 (932)	2	55 (131)	*	*	755 (1391)*	1 min *	—
Correction material	500 (932)	4	55 (131)	500 (932)	715 (1319)	715 (1319)	1 min (under vacuum)	—
Gingival	500 (932)	6	55 (131)	500 (932)	715 (1319)	715 (1319)	1 min (under vacuum)	—

* The glaze firing can be carried out with or without vacuum. The glaze level is increased by extending the holding time.

Entry temperature must be right before objects are placed.

To this end, mix Triceram® Neutral transparent with Modelling Liquid. Place the sample on a platinum foil and fire on the basis of the firing parameters for Triceram® in the dentin firing 1.

14. Shade in accordance with the manufacturer's shade guide

Triceram® is based on the shade systems VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® (VITA Classical A1-D4® and VITA SYSTEM 3D-MASTER® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co KG, 79704 Bad Säckingen, Germany).

15. Information on combining with other products

The dental ceramic may only be combined with other products if recommended by Dentaaurum. Only frameworks may be veneered where the coefficient of thermal expansion lies within the given range for veneering ceramics.

The dental ceramic may not be used for frameworks made of lithium disilicate, zirconium oxide, titanium, titanium alloys or dental alloys with a silver content exceeding 30%.

16. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaaurum.com.

For more information on our products and services, visit www.dentaaurum.com. You will find these Instructions for use and information on the current status at www.dentaaurum.com/ifu.

In case of questions or ideas, please contact your local representative.

17. Explanation of symbols used on the label

Refer to the label. Additional information can be found at www.dentaaurum.com/ifu (Explanation of symbols REF 989-313-00).

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