discovery® | discovery® smart |
topic | equilibrium® mini |
equilibrium® 2 | equilibrium® ti

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo
Modalità d'uso

D
DENTAURUM

- ☞ Informationen zu Produkten finden Sie unter www.dentaurum.com
- ☞ For more information on our products, please visit www.dentaurum.com
- ☞ Vous trouverez toutes les informations sur nos produits sur www.dentaurum.com
- ☞ Descubra nuestros productos en www.dentaurum.com
- ☞ Informazioni su prodotti sono disponibili nel sito www.dentaurum.com

Stand der Information | Date of information | Mise à jour |
Fecha de la información | Data dell'informazione: 2023-02

Änderungen vorbehalten | Subject to modifications |
Sous réserve de modifications | Reservado el derecho de modificación |
Con riserva di apportare modifiche

D
DENTAURUM

Turnstr. 31 | 75228 Ispringen | Germany | Tel. +49 7231/803-0 | Fax +49 7231/803-295
www.dentaurum.com | info@dentaurum.com

989-804-44 Printed by Dentaurum Germany 2023-02

DE

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General description

2.1 discovery®/topic – Premium metal brackets

Our premium metal brackets include the discovery® bracket family and the topic bracket. They are the result of many years of experience. discovery® metal brackets are designed as a one-piece bracket. They are manufactured on state-of-the-art machines using the MIM technique (Metal Injection Molding) with a special dental alloy. A nickel-free alloy is used for topic brackets. They are optimally adapted to the tooth contour using CAD software (Computer Aided Design) and are therefore easy to handle.

2.2 equilibrium® - Standard metal brackets

Our equilibrium® brackets are part of Dentaurum's standard metal bracket range. They are formed and milled out of rod material and are processed in many tried-and-tested steps.

2.3 Bracket base

The laser-structured base of all named Dentaurum metal brackets ensures the best possible bonding and guarantees easy debonding.

3. Intended purpose

For the treatment of tooth misalignments as part of orthodontic treatment.

4. Intended users

All products are intended for dental professionals only.

5. Indications

For tooth and jaw misalignments.

6. Contraindications

Stainless steel brackets must not be used in patients with a nickel allergy. Warnings are found on the respective product labels. Also see 9. Composition.

7. Patient target group

All patients for whom a fixed appliance with orthodontic braces can be placed as part of orthodontic treatment.

8. Processing

8.1 Informing the patient

Inform the patient about the possible consequences and risks of bracket treatment.

8.2 Conditioning teeth

Check each tooth thoroughly for any existing damage to the enamel to avoid damaging the tooth during debonding.

8.3 Selection of bracket adhesive

discovery®, topic and equilibrium® brackets can be bonded using all commercially available bracket adhesives.

8.4 Removing and handling

Grip the brackets mesially-distally using crimping tweezers to remove them from the packaging. The laser-structured base of the brackets should not be touched with your fingers or contaminated in any way as this could otherwise have a detrimental effect on the bond. Detailed instructions on removing the brackets from the bracket blister can be found in „Notes on removing the brackets“ (REF 989-802-39) at www.dentaurum.com.

8.5 Using the bracket adhesive

Prepare the bracket adhesive according to the manufacturer's instructions. Apply a small amount of adhesive to the center of the bracket base, enough to wet the entire bracket base. This guarantees an optimum bond.

8.6 Placing the brackets - direct bonding

After applying the adhesive, place the bracket in the best possible position on the tooth and press it down on the tooth. Please completely remove any excess adhesive before it has cured.

8.7 Placing the brackets - indirect bonding

After applying the adhesive, place the transfer tray with the brackets as usual in the mandible and/or maxilla. Completely remove any excess adhesive before it has cured.

8.8 Ligating archwires

Once the adhesive has cured, ligate the archwires as usual with wire or plastic ligatures.

8.9 Bracket debonding

For the debonding of discovery®, topic and equilibrium® brackets, we recommend using Weingart pliers (e.g. REF 003-120-00). To this end, grasp the bracket from mesial-distal on the tie wings. Remove the base carefully from the tooth by pressing the pliers together and bending the bracket away from the tooth. If necessary, turn the bracket briefly to the right or the left, without pulling.

9. Composition

- discovery® and discovery® smart brackets are made of the following material: DIN 1.4404.
- equilibrium® brackets are made of the following materials: DIN 1.4301, 1.4303 and 1.4404.
- topic brackets are made of a CoCr alloy.
- equilibrium® ti brackets are made of the following materials: DIN 3.7035 and 3.7165.

The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

10. Safety information

The product should not be used if there is a known allergic reaction to one or more of the material components. Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly. topic brackets (CoCr alloy) must not be ground in the patient's mouth or without a suction apparatus to prevent the release of cobalt. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

11. Storage and shelf life

If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

12. Information for single use products

The product is intended for single use only. Reconditioning of brackets that have already been used (recycling) and their reuse on a patient is not permitted.

13. Product range

Please refer to our Orthodontics Catalog for information on the products available in the bracket range.

14. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

15. Explanation of symbols used on the label

Please refer to the label. These Instructions for use and additional information can be found at www.dentaurum.com (explanation of symbols REF 989-313-00).

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