CITO mini®

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

The CITO mini® O-ring matrices and the O-rings for CITO mini® matrices are prefabricated appliance retaining elements used to anchor implant-supported dentures in the mandible and/or maxilla in conjunction with CITO mini® implants. The one-piece implants are fitted with a ball Ø 1.8 mm, on top of which matrices of the system are fixated once they have been polymerized into the patient's dental prosthesis. Depending on the indication, there are matrices with different withdrawal forces and freedom of movement for divergencies between the implants.

CITO mini® implant data			
Implant diameters in mm	1.8	2.2	2.5
Implant lengths in mm	11/13	11/13	11/13
Ball diameter in mm*	1.8		

*all O-ring matrices listed in these instruction for use can be processed with all CITO mini® implants, as the ball head diameter is identical at 1.8 mm.

CITO mini® ball abutment laboratory implants are used to create working casts. They precisely reproduce the position of the CITO mini® implants, including the ball head. Processing and use take place exclusively in the (practice) laboratory or in the manufacturing center and not on the patient. This is not a medical device.

CITO mini® PentaGrip insertion key, ISO shank for handpiece insertion, tightening, or removal of CITO mini® implants. The insertion key can be used with the handpiece set to the relevant torque. It is intended for multiple use and reprocessing.

	Length in mm
CITO mini® PentaGrip insertion key, ISO shank	22.3

The adapters – ISO shank / ratchet – allow manual handling of the handpieces and components (ISO shank) with the torque ratchet or the sure-grip wheel. They are intended for multiple use and reconditioning.

	Length in mm
Adapter, ISO shank, hexagon/ratchet	15.0 / 20.0

To better access hard-to-reach areas in the mouth with handpiece instruments, they can be connected to the drill extension via the ISO shank. The drill extension is intended for multiple use and reconditioning.

	Length in mm
Drill extension – ISO shank hexagon	21.0

The CITO mini® implant system contains specially coordinated abutments and accessories, including instruments, for the insertion and prosthetic restoration of the implants.

3. Intended purpose

O-ring matrices are intended for connecting the implants in the maxilla and mandible with the prosthetic restoration and restoration of the chewing function.

O-rings for CITO mini® O-ring matrices are components regularly replaced in the CITO mini® O-ring matrices to maintain the optimum retention force of the denture.

CITO mini® ball abutment laboratory implants are designed to represent the position of the abutment (axial inclination and height) in a working cast.

Instruments for screwing implants and screws for abutments and impression components in and out.

The adapters - ISO shank, hexagon/ratchet - are intended for manually handling the motorized rotary instruments and components (ISO shank) with the ratchet.

Drill extensions are intended to extend rotary instruments with a handpiece connection.

4. Intended users

The CITO mini® product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

5. Medical indications / contraindications

5.1 Indications

The CITO mini® O-ring matrices are used to stabilize implant-supported prostheses in the edentulous or partially edentulous maxilla and mandible on CITO mini® implants. The indication for an implant should be made with consideration of possible advantages and disadvantages, as well as risks of implant treatment and alternative treatments.

To fixate full dentures, Dentaurum recommends at least four implants in the mandible and six implants in the maxilla. For restorations in the maxilla, use CITO mini® implants Ø 2.2 mm or Ø 2.5 mm.

When positioning the implants, pay attention to keeping an adequate distance between the implants and to ensuring that the static force on the bone is optimally applied. The size of the matrix should also be considered.

Immediate implant loading is only possible if a torque of 35 Ncm was reached during insertion. Torque should not exceed 40 Ncm during insertion.

5.1.1 CITO mini® O-ring matrices

- Restoration of full/partial and free-end dentures in combination
- Removable dental restorations on CITO mini® implants
- Immediate restoration of CITO mini® implants with an insertion torque of at least 35 Ncm
- Incorporation in an existing coverdenture as an interim denture or with partial denture framework possible

5.1.2 O-rings for CITO mini® O-ring matrix

- Replacement article for mounted O-rings in CITO mini® O-ring matrices
- Use of O-rings for matrix 1 in CITO mini® O-ring matrices 1
- Use of O-rings for matrix 2 + 3 in CITO mini® O-ring matrices 2 + 3

5.1.3 CITO mini® ball abutment laboratory implant

- precise reproduction of positions for cast models (dental stone)

5.1.4 CITO mini® PentaGrip insertion key, ISO shank

- contact-free removal of the CITO mini® implants from the implant packaging
- handpiece or manual insertion or removal of the CITO mini® implants into/from the prepared drill hole, taking torque specifications into account

5.1.5 Adapter ISO shank / ratchet

- Adapters for instruments with ISO shank connection enable manual insertion or removal
- For use by hand, with the sure-grip wheel or the torque ratchet

5.1.6 Drill extension

- Extension of instruments with ISO shank connection for handpiece or manual insertion or removal

5.2 Contraindications

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically poor prosthetic restoration can cause premature implant loss and fatigue fractures in implants and abutments under biomechanical loading. CAUTION: The insertion of CITO mini® implants in Ø 1.8 mm in soft bone and/or the maxilla is not indicated.

5.2.1 CITO mini® O-ring matrices

- Use with third-party products
- Use with incorrect ball diameter
- Exceeding the specified angle tolerances
- Failure to observe the matrix diameters
- Immediate restoration of CITO mini® implants with an insertion torque of at least 35 Ncm
- Heavy mechanical strain with Ø 1.8 mm implants

Matrices	Angle tolerance	Retentive force	Height	∅
CITO mini®O-ring matrix 1	0° - 25°	3 N	3.55 mm	4.7 mm
CITO mini® O-ring matrix 2	0° - 15°	4 N	3.3 mm	4.1 mm
CITO mini® O-ring matrix 3	0° - 5°	4 N	3.2 mm	4.1 mm

5.2.2 O-rings for CITO mini® O-ring matrix

- Use of O-rings for matrix 1 in CITO mini® O-ring matrices 2 + 3
- Use of O-rings for matrix 2 + 3 in CITO mini® O-ring matrix 1
- Exceeding the specified angle tolerances
- Use with incorrect ball diameter

5.2.3 CITO mini® ball abutment laboratory implant

- Use of laboratory implants with defective or damaged connector geometry
- Use on patients

5.2.4 CITO mini® PentaGrip insertion key, ISO shank

- Use as an insertion key for other abutments and components
- Use with incorrect torque
- Use without torque limit
- Use without the insertion key properly engaging with the implant

5.2.5 Adapter – ISO shank / ratchet

- Use as insertion key or hex key
- Use with incorrect torque
- Use without torque limit

5.2.6 Drill extension

- Use as insertion key or hex key
- Use with incorrect torque
- Use without torque limit

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the CITO mini® Instructions for use.

We strongly recommend that, prior to using Dentaaurum products, dental professionals attend a Dentaaurum training course to familiarize themselves with the Dentaaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaaurum recommends full clinical, radiological, photographic and statistical documentation.

Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Important:

The implantologist should ensure the products cannot be aspirated during intra-oral use.

7.1 CITO mini® ball abutment laboratory implant

7.1.1 Impression-taking:

With the indirect processing version, the impression can be taken directly via the ball heads. Apply impression material around all ball heads, fill the closed impression tray with impression material and take the impression. Remove the tray after the impression material is cured. Ensure that the impression has captured the ball heads exactly.

7.1.2 Casting the model

The ball abutment laboratory implant is then repositioned in the impression in the laboratory. The indent beneath the ball ensures that the implant axes are transferred precisely. An impression is cast in dental stone and given a base to make a model. The laboratory implant must fit precisely on the model. After the dental stone has cured, carefully remove the impression from the model.

A ball abutment laboratory implant is available for all CITO mini® implant types, as all balls have a diameter of 1.8 mm.

7.2 CITO mini® O-ring matrices

Ball-anchored dentures are tissue-borne restorations. When using ball abutments existing coverdentures can be used as temporary dentures or modified with a cast partial denture framework – alternatively a new coverdenture can be fabricated.

CITO mini® implants function in such a manner that it is advisable to support the prosthetic restoration on at least four implants in the mandible and six implants in the maxilla.

7.2.1 Indirect processing version:

The CITO mini® O-ring matrices are inserted on the CITO mini® ball abutment laboratory implants and their angle of insertion must be aligned parallel and axially in the insertion direction. Ensure that the dimensions and angle tolerances of the selected O-ring matrix are observed. In order to protect the balls against acrylic seeping in and to ensure that the matrices are parallel in the insertion direction, block out the undercuts with wax. Grind the prosthesis in the area of the ball heads such that the matrices have sufficient clearance. To check the exact fit of the matrices, drill a window from the lingual or palatal aspect. Secure the matrices with cold-curing resin through the opening. Observe the instructions for use of the resin. Once the resin has cured, remove the denture from the model. Remove excess resin on the lower margins of the matrices. Smooth and polish the denture. Check that the ball abutments function correctly on the working cast.

If a new full denture is fabricated, use exactly the same procedures, except that the matrices are incorporated into the denture base.

7.2.2 Direct processing version:

In the direct technique, incorporate the matrices into an existing denture directly in the patient's mouth without fabricating a model. No laboratory implant is needed. During the clinical procedures, protect the undercuts with the help of spacers. Position them between the gingiva and the edge of the matrix to prevent cold-curing resin filling the undercuts in the mouth which would prevent the denture being removed after the resin has cured. Follow the instructions for use from the polymer and spacer manufacturers.

Important:

The implants must not diverge by more than 25°.

The lower edge of the ball head should be approx. 1.8 mm above the gingival line. To achieve optimum retention, position all ball heads at the same level. The ball head is 1.8 mm in diameter. The components of the ball head may not be modified.


Matrices	Angle tolerance	Retentive force	Height	∅
CITO mini®O-ring matrix 1	0° - 25°	3 N	3.55 mm	4.7 mm
CITO mini® O-ring matrix 2	0° - 15°	4 N	3.3 mm	4.1 mm
CITO mini® O-ring matrix 3	0° - 5°	4 N	3.2 mm	4.1 mm

7.2.3 Monitoring:

Dentures and their retention units must be monitored at six-monthly intervals. Take the following points into account:

- Eliminate unfavorable movements of the denture (optimize relining of the denture or replace the matrices)
- Oral hygiene (remove plaque and dental calculus and, if necessary, re-instruct the patient on cleaning implants)
- Regularly changing the O-rings in the matrices guarantees optimum retention of the denture

7.3 CITO mini® PentaGrip insertion key, ISO shank


PentaGrip insertion key

Grasp the CITO mini® implant with the CITO mini® PentaGrip insertion key. The insertion key snaps into place on the implant when turned slightly. The insertion key must engage fully on the antirotational mechanism. The implant is otherwise not secure and may be damaged during insertion. The CITO mini® PentaGrip insertion key is designed to ensure contact-free insertion, regardless of the indication. The PentaGrip insertion key can be extended using a drill extension and can be used manually by attaching the adapter.

Maintain a maximum motor speed of 10 rpm during handpiece insertion. Excessive torque or motor speed can damage the implant site.

7.4 Adapter ISO shank / ratchet

Adapter, ISO shank, hexagon - ratchet

The adapter – ISO shank / ratchet allows manual handling of the handpieces and components (ISO shank) with the torque ratchet or the sure-grip wheel. Depending on the available space and the situation, two different adapter lengths can be selected. Instruments with mechanical connection are coupled via the ISO shank to the ISO shank/ratchet adapter and locked in place by the spring on the

adapter. They can then be used with the torque ratchet or sure-grip wheel, taking into account the respective torque values for the different screws (AF 1.3). They are intended for multiple use and reconditioning.

7.5 Drill extension

Mechanical ISO shank hexagon extension

To better access hard-to-reach areas in the mouth with handpiece instruments, they can be connected to the drill extension via the ISO shank. They are intended for multiple use and reconditioning.

8. Technical specification

Composition

Product	Composition
CITO mini® O-ring matrices	3.7165 / NBR
CITO mini® O-ring for O-ring matrices	NBR
CITO mini® laboratory implant	Al alloy
CITO mini® PentaGrip insertion key, ISO shank	1.4197 / 1.4310
Adapter, ISO shank, hexagon/ratchet	1.4197
Drill extension, ISO shank hexagon	1.4197

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during surgical treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during surgical treatment (permanent adverse reactions) are:

Chronic pain in connection with the implant, permanent paresthesia, dysesthesia, nerve damage, exfoliation, hyperplasia, local or systemic infections, oroantral or oronasal fistulas, alveolar bone loss in maxilla/mandible, irreversible damage to adjacent teeth, fractures to the implant, jaw, bone or restoration, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- All CITO mini® prosthetic components may only be used in combination with CITO mini® implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the CITO mini® implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

CITO mini® matrices and the O-rings for matrices are delivered non-sterile. All prosthetic abutment components are for single use on one patient only.

Reprocessing of CITO mini® abutment components (recycling) or reuse on patients are not permitted as safe reprocessing and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

The CITO mini® PentaGrip insertion keys, ISO shank/ratchet adapters and drill extensions are delivered non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. Drills, instruments and components must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The CITO mini® matrices and the O-rings for matrices are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Refer to the CITO mini® order form for the CITO mini® product range.

17. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaurum.com/ifu

Instructions for use	REF
CITO mini® Planning Scheme	989-502-45
CITO mini® implants	989-365-98
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurum products	989-804-43

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