tiologic[®]
TWINFIT
tiologic[®]ST

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EN Instructions for use · Screws

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1. Hersteller

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

The various screws fix prosthetic abutments, caps, as well as components for implant closure and components for analog and digital impressions on the tioLogic® ST and tioLogic® TWINFIT implants.

Abutment components					
Series of abutments	S	M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5

Retaining screws for tioLogic® ST and tioLogic® TWINFIT impression posts, closed and scan abutments fix the tioLogic® ST and tioLogic® TWINFIT impression post, closed during impression taking with a closed impression tray and the scan abutment on the (laboratory) implants.

tioLogic® ST and tioLogic® TWINFIT impression components, closed and scan abutments			
Series of abutments	S	M	L
Length of retaining screws in mm	14.0		

AnoTite screws 9.0 mm fixate tioLogic® ST and tioLogic® TWINFIT abutments, which are screw-retained directly onto the implant with the corresponding tioLogic® ST and tioLogic® TWINFIT implant.

For final fitting, an unused AnoTite screw is recommended!

The 9.0 mm prosthetic screw is available for processing the tioLogic® ST and tioLogic® TWINFIT abutments in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

AnoTite screw and prosthetic screw 9.0 mm			
Series of abutments	S	M	L
Length of AnoTite screw and prosthetic screws in mm	9.0		

AnoTite 6.0 mm screws fixate tioLogic® ST bar, bridge, and AngleFix caps and tioLogic® TWINFIT 4Base caps, that are screw-retained directly to the respective tioLogic® ST bar, bridge, and AngleFix caps and tioLogic® TWINFIT 4Base abutment.

For final fitting, an unused AnoTite screw is recommended!

The 6.0 mm prosthetic screw is available for processing tioLogic® ST bar, bridge, AngleFix caps and tioLogic® TWINFIT 4Base caps in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

AnoTite screw and prosthetic screw 6.0 mm				
Series of abutments	Bar abutment	Bridge abutment	AngleFix abutment	4Base Abutment
Length of AnoTite screw and prosthetic screws in mm	6.0			

AnoTite screws for angulated screw apertures 8.5 mm fixate tioLogic® TWINFIT CAD/CAM titanium bases via an angulated screw aperture directly screw-retained on the tioLogic® TWINFIT implant.

For final fitting, an unused AnoTite screw for angulated screw aperture is recommended!

The 8.5 mm prosthetic screw for angulated screw apertures is available for processing the tioLogic® TWINFIT CAD/CAM titanium bases in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

AnoTite screw and prosthetic screw, for angulated screw apertures 8.5 mm			
Series of abutments	S	M	L
Length of AnoTite screw and prosthetic screws in mm	8.5		

AnoTite screws for angulated screw apertures 5.5 mm fixate bridges and bar constructions via an angulated screw aperture screw-retained directly onto the tioLogic® TWINFIT 4Base abutments.

For final fitting, an unused AnoTite screw for angulated screw aperture is recommended!

The 5.5 mm prosthetic screw for angulated screw apertures is available for processing the tioLogic® TWINFIT 4Base caps in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

AnoTite screw and prosthetic screw, for angulated screw apertures 5.5 mm	
Series of abutments	4Base Abutment
Length of AnoTite screw and prosthetic screws in mm	5.5

3. Intended purpose

Screws are used to fix abutments on the implants.

Prosthetic screw 9.0 mm is intended to fixate the tioLogic® ST and tioLogic® TWINFIT abutments during production of the dental restoration on the model.

Prosthetic screw 6.0 is intended to fixate the tioLogic® ST bar, bridge and AngleFix caps and tioLogic® TWINFIT 4Base caps during production of the dental restoration on the model.

Prosthetic screw 8.5 mm for angulated screw apertures is intended to fixate the tioLogic® TWINFIT CAD/CAM titanium bases to the model during production of the dental restoration.

Prosthetic screw 5.5 mm for angulated screw apertures is intended to fixate the tioLogic® TWINFIT 4Base caps to the model during production of the dental restoration

4. Intended users

The tioLogic® ST and tioLogic® TWINFIT product ranges are supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by surgeons, dentists or dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restorations.

5. Medical indications / contraindications

5.1 Indications

5.1.1 Retaining screw

- For use with the tioLogic® ST and tioLogic® TWINFIT impression post, closed
- For use with the tioLogic® ST and tioLogic® TWINFIT scan abutments

5.1.2 AnoTite screw 9.0 mm

- For use with tioLogic® ST and tioLogic® TWINFIT abutments for direct fixation on tioLogic® ST and tioLogic® TWINFIT implants

5.1.3 Prosthetic screw 9.0 mm

- For use with tioLogic® ST and tioLogic® TWINFIT abutments for direct fixation on tioLogic® ST and tioLogic® TWINFIT laboratory implants on the model for the production of dental restorations
- Use exclusively in the laboratory

5.1.4 AnoTite screw 6.0 mm

- Fixation of tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base caps on the tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base abutments

5.1.5 Prosthetic screw 6.0 mm

- Fixation of tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base caps on the tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base abutments or the respective laboratory implants on the model for the production of dental restorations
- Use exclusively in the laboratory

5.1.6 AnoTite screw 8.5 mm for angulated screw apertures

- Fixation of tioLogic® TWINFIT abutments directly over the tioLogic® TWINFIT implant via an angulated screw aperture

5.1.7 Prosthetic screw 8.5 mm for angulated screw apertures

- For use with tioLogic® TWINFIT CAD/CAM titanium bases for direct fixation on tioLogic® TWINFIT laboratory implants on the model for the production of dental restorations
- Use exclusively in the laboratory

5.1.8 AnoTite screw 5.5 mm for angulated screw apertures

- Fixation of bridge and bar constructions directly over the tioLogic® TWINFIT 4Base abutment via an angulated screw aperture

5.1.9 Prosthetic screw 5.5 mm for angulated screw apertures

- Fixation of tioLogic® TWINFIT 4Base caps on the tioLogic® TWINFIT 4Base abutments or 4Base laboratory implants on the model for the production of dental restorations
- Use exclusively in the laboratory

5.2 Contraindications

5.2.1 Retaining screw with shortened thread

- For use with other abutments except on tioLogic® ST and tioLogic® TWINFIT impression posts, closed as well as tioLogic® ST and tioLogic® TWINFIT scan abutments

5.2.2 AnoTite screw 9.0 mm

- For use with other abutments except tioLogic® ST and tioLogic® TWINFIT abutments for direct fixation on tioLogic® ST and tioLogic® TWINFIT implants

5.2.3 Prosthetic screw 9.0 mm

- For use with other abutments except tioLogic® ST and tioLogic® TWINFIT abutments for direct fixation on tioLogic® ST and tioLogic® TWINFIT laboratory implants on the model for the production of dental restorations
- Not for use on or in patients

5.2.4 AnoTite screw 6.0 mm

- Fixation of other abutments except tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base caps on the tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base abutments

5.2.5 Prosthetic screw 6.0 mm

- Fixation of other abutments except tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base caps on the tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT 4Base abutments or laboratory implants on the model for the production of dental restorations
- Not for use on or in patients

5.2.6 AnoTite screw 8.5 mm for angulated screw apertures

- Fixation of other abutments except tioLogic® TWINFIT CAD/CAM titanium bases directly over the tioLogic® TWINFIT implant via an angulated screw aperture

5.2.7 Prosthetic screw 8.5 mm for angulated screw apertures

- For use with other abutments except tioLogic® TWINFIT CAD/CAM titanium bases for direct fixation on tioLogic® TWINFIT laboratory implants on the model for the production of dental restorations
- Not for use on or in patients

5.2.8 AnoTite screw 5.5 mm for angulated screw apertures

- Fixation of other abutments except bridge and bar constructions directly over the tioLogic® TWINFIT 4Base abutment via an angulated screw aperture

5.2.9 Prosthetic screw 5.5 mm for angulated screw apertures

- Fixation of other abutments except tioLogic® TWINFIT 4Base caps on the tioLogic® TWINFIT 4Base abutments or 4Base laboratory implants on the model for the production of dental restorations
- Not for use on or in patients

6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

7. Information on use, processing, fitting and insertion

Prior to use, dental professionals should ensure that they have carefully read and understood the full tioLogic® ST and tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaureum products, dental professionals attend a Dentaureum training course to familiarize themselves with the Dentaureum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaureum recommends full clinical, radiological, photographic and statistical documentation.

The implantologist should ensure the products cannot be aspirated during intra-oral use.

Precautions:



Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.

Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

7.1 Retaining screw with shortened thread

7.1.1 tioLogic® ST and tioLogic® TWINFIT impression posts, closed


Impression post M, closed Impression cap M

Components for the closed impression technique include impression posts, retaining screws and impression caps. Impression posts and caps are laser-marked or marked with the S, M or L series of abutments.

After the impression has been taken, an individual tray is fabricated. It is designed to be large enough in the region of the implants so it does not come into contact with the impression posts, including the attached impression caps.

Remove the cover screws, the gingiva formers and the temporary restoration prior to taking the impression.

For easier fitting, the impression posts for closed-tray impressions can be mounted on positioning aids in the corresponding series of abutments (S, M, or L). Fit the impression post for closed impression-taking with the long tube into the implant. Then align the antirotational mechanism so the PentaStop interlocks. The screw connection can be made occlusally through the positioning aid. The enclosed red retaining screw has a shortened thread which will only grip in the (laboratory) implant if the impression post has been inserted in the correct position into the connection. Align the impression post again if necessary and check to ensure that it fits

correctly (X-ray check of the proper fit of the impression post). Tighten the retaining screw with the hex key AF 1.3.

Fit the corresponding impression cap S, M or L according to the vertical retention grooves on the impression post until it perceptibly and audibly clicks into place.

IMPORTANT:

The impression caps are items for single use only and are NOT suitable for sterilization or multiple use. Multiple use results in transfer inaccuracies.

The orientation of the retention grooves ensures that they can be positioned without coming into contact with the adjacent teeth. The impression is to be taken with a material based on silicone or polyether. The impression caps are secured in the impression material by the retention. Ensure that the peri-implant region is accurately reproduced in the impression. Remove the impression tray after curing the impression compound. Send the impression posts with fixing screws together with the impression to the laboratory. Use the practice record card (REF 989-966-22) to provide the dental technician with all the relevant information. The respective gingiva formers are refitted after impression-taking.

Tightening torque	
Retaining screw with shortened thread on the model / in the mouth	15 Ncm

7.1.2 tioLogic® ST and tioLogic® TWINFIT scan abutments

The tioLogic® TWINFIT scan abutments, titanium, are available in the S, M or L series of abutments. Place the scan abutment intra-orally on the tioLogic® TWINFIT implant or extra-orally on the laboratory implant in the master model (pay attention to the rotational security) and secure it with the hex key AF 1.3 and the enclosed red fixing screw. The long cylinder at the connection point of the scan abutment, titanium, serves for an easy integration into the (laboratory) implant and for exact axial detection of the position. The enclosed red retaining screw features a shortened thread, which only grips in the (laboratory) implant when the scan abutment has been inserted into the connection point in the correct position. The surface of the scan abutments titanium is optimized for digital capture, both intra-orally as well as on the model without scanning spray. The scanning and matching process of the scan data within the respective CAD software is performed according to the specifications of the individual software manufacturer. The scan abutments, titanium, are intended for multiple use while maintaining the integrity of the scan body and the connection to the (laboratory) implant. The retaining screw is intended for single use and is marked on the blister packaging with .

Tightening torque	
Retaining screw with shortened thread on the model / in the mouth	15 Ncm

7.2 AnoTite screw 9.0 mm

AnoTite screws 9.0 mm are long-term connector elements which fixate tioLogic® ST and tioLogic® TWINFIT abutments which are screw-retained directly onto the implant with the corresponding tioLogic® ST and tioLogic® TWINFIT implant.

For final fitting, an unused AnoTite screw is recommended!
The 9.0 mm prosthetic screw is available for processing the tioLogic® ST and tioLogic® TWINFIT® abutments in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

The AnoTite 9.0 mm screw and the 9.0 mm prosthetic screw are fixed to the respective (laboratory) implant using the AF 1.3 hex key, taking into account the specified torque.

Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

Tightening torque	
tioLogic® ST and tioLogic® TWINFIT titanium abutments, CAD/CAM abutments and bases	
Prosthetic screw on the model	30 Ncm
AnoTite screw intra-orally	30 Ncm
tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT® 4Base abutments	
Prosthetic screw on the model	35 Ncm
AnoTite screw intra-orally	35 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.3 AnoTite screw 6.0 mm

AnoTite 6.0 mm screws are long-term connector elements to fixate tioLogic® ST bar, bridge and AngleFix caps and tioLogic® TWINFIT® 4Base caps that are screw-retained directly to the respective tioLogic® ST bar, bridge, AngleFix caps and tioLogic® TWINFIT 4Base abutment.

For final fitting, an unused AnoTite screw is recommended!

The 6.0 mm prosthetic screw is available for processing tioLogic® ST bar, bridge, AngleFix caps and tioLogic® TWINFIT 4Base caps in the laboratory.

The prosthetic screw is intended for laboratory use only and is not a medical device!

The AnoTite 6.0 mm screw and the 6.0 mm prosthetic screw are fixed to the respective (laboratory) implant using the AF 1.3 hex key, taking into account the specified torque.

Caution:

Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

Tightening torque	
tioLogic® ST bar, bridge and AngleFix and tioLogic® TWINFIT® 4Base caps	
Prosthetic screw on the model	25 Ncm
AnoTite screw intra-orally	25 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.4 AnoTite screws for angulated screw apertures

The angulated screw aperture offers the optimum solution for the complete segment that is occlusally screw-retained. The screw aperture can be inclined at an angle of 20° to the implant axis, allowing discrete emergence in the palatal area, particularly in the anterior region.

7.4.1 AnoTite screws and prosthetic screws 8.5 for angulated screw apertures

AnoTite screws 8.5 mm for angulated screw apertures are long-term connector elements and fixate tioLogic® TWINFIT CAD/CAM titanium bases which are screw-retained directly onto the implant with the tioLogic® TWINFIT implant.


For final fitting, an unused AnoTite screw is recommended!

The 8.5 mm prosthetic screw with an angulated screw aperture is available for processing the tioLogic® TWINFIT® CAD/CAM titanium bases in the laboratory.

Note:


The prosthetic screw is intended for laboratory use only. This is not a medical device!

7.4.2 AnoTite screws and prosthetic screws 5.5 for angulated screw apertures

AnoTite screws 5.5 mm for angulated screw apertures are long-term connector elements and fixate tioLogic® TWINFIT 4Base caps which are directly screw-retained on the tioLogic® TWINFIT 4Base abutment.

For final fitting, an unused AnoTite screw is recommended!

The 5.5 mm prosthetic screw with an angulated screw aperture is available for processing the tioLogic® TWINFIT 4Base caps with an angulated screw aperture in the laboratory.

Note:

The prosthetic screw is intended for laboratory use only. This is not a medical device!

The AnoTite screw and prosthetic screw for angulated screw apertures 8.5 mm and 5.5 mm are fixed to the respective (laboratory) implant or 4Base abutment using the hex ball point key, ISO shank AF 1.3, taking into account the specified torque.

Caution:

Do not confuse the prosthetic screw with the AnoTite screw for angulated screw apertures. The prosthetic screws are colored silver and the AnoTite screws are anodized blue.

Tightening torque	
tioLogic® TWINFIT CAD/CAM bases	
AnoTite screw 8.5 mm for angulated screw apertures intra-orally	25 Ncm
Prosthetic screw 8.5 mm for angulated screw apertures on the model	25 Ncm
tioLogic® TWINFIT® 4Base caps	
AnoTite screw 5.5 mm for angulated screw apertures intra-orally	25 Ncm
Prosthetic screw 5.5 mm for angulated screw apertures on the model	25 Ncm

To achieve optimum preloading of the AnoTite screw, the AnoTite screw must be retightened with the same torque approx. 5 minutes after insertion.

7.5 Additional components

Overview of additional components:

Torque ratchet
Hex key
Adapter, ISO shank, hexagon/ratchet

8. Technical specification

Composition

Product	Composition
Fixing screws with shortened thread for impression posts, closed and scan abutment	3.7165
AnoTite screw 9.0 and 6.0 mm	3.7165
Prosthetic screw 9.0 and 6.0 mm	3.7165
AnoTite screws for angulated screw apertures 8.5 and 5.5 mm	3.7165
Prosthetic screws for angulated screw apertures 8.5 and 5.5 mm	3.7165

Information about the composition is included in the material list for our implantology products (REF 989-801-06), or see www.dentaurum.com/ifu. The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use all tioLogic® ST prosthetic components in combination with tioLogic® ST implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Only use all tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

11. Storage and shelf life

Store the tioLogic® TWINFIT implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

13. Information for single use products

Prosthetic abutments and screws are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

14. Information for multiple use products

The hex key AF 1.3 and hex ball point key, ISO shank AF 1.3, are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used initially and before being used again. All reprocessed items must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

15. Information on disinfection / cleaning / sterilizing

The screws for securing prosthetic abutments, caps, as well as the components for implant closure and the components for analog and digital impressions are supplied non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaureum generally recommends that the product be sterilized in addition, as non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

16. Product overview and scope of delivery

Please refer to the tioLogic® ST product catalog for the tioLogic® ST product range. Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

17. Quality information

Dentaureum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaureum has no influence on the processing on site.

18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

19. Information on MRT safety

The products within these Instructions for use are MRT safe.

20. Additional information

Additional information can be found at www.dentaureum.com/ifu

Instructions for use	REF
tioLogic® ST closure screws, gingiva formers, impression posts, temporary abutments	989-801-86
tioLogic® TWINFIT closure screws, gingiva formers, impression posts, temporary abutments	989-365-90
tioLogic® ST CAD/CAM abutments and components	989-802-88
tioLogic® ST titanium abutments and components	989-801-87
tioLogic® ST bar abutments and components	989-802-38
tioLogic® ST bridge abutments and components	989-801-91
tioLogic® ST AngleFix abutments and components	989-801-89
tioLogic® TWINFIT CAD/CAM abutments and components	989-365-89
tioLogic® TWINFIT titanium abutments and components	989-365-91
tioLogic® TWINFIT 4Base abutments and components	989-365-93
tioLogic® ST Additional components	989-805-87
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaureum products	989-804-42

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