**D**  
DENTAURUM  
1886

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---

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

### 2. General device description

The pOPosition for tioLogic® drills and components are intended for the template-guided preparation of the implant site for tioLogic® ST and tioLogic® TWINFIT implants.

The pOPosition for tioLogic® drills are supplied non-sterile and packaged individually. The pOPosition for tioLogic® drills feature an integrated depth stop. The ADVANCED for tioLogic® thread tap is used for both the pOPosition for tioLogic® drill system as well as the ADVANCED for tioLogic® drill system. In addition, the instruments are color-coded according to their specific diameter.

Implant diameters in mm	3.3	3.7	4.2	4.8	5.5
Color coding	brown	yellow	red	blue	green
Implant lengths in mm	9/11/13/15	7/9/11/13/15			9/11/13/15

The sleeves for the pOPosition for tioLogic® system are manufactured from titanium grade 5 and guarantee exact guidance of the corresponding pOPosition for tioLogic® drills. For easier identification, the sleeves are color-coded according to their diameter. They have a grip section which, due to its special bending zone, can be adapted three-dimensionally according to the amount of space available orally. This ensures optimal use of the inner sleeves, particularly in the distal region.

The overview shows the allocation of the implant diameters to the corresponding sleeves (basic sleeves / inner sleeves):

	S	M	L		
Implant diameter	∅ 3.3 mm	∅ 3.7 mm	∅ 4.2 mm	∅ 4.8 mm	∅ 5.5 mm
Basic sleeve					
Inner diameter	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm
Outer diameter	5.1 mm	5.1 mm	5.6 mm	6.2 mm	6.9 mm
Inner sleeve Depth drill					
Inner diameter	2.05 mm	2.05 mm	2.05 mm	2.05 mm	2.05 mm
Outer diameter	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm
Inner sleeve Stepped countersink					
Inner diameter	2.85 mm	3.15 mm	3.65 mm	4.25 mm	4.95 mm
Outer diameter	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm

### 3. Intended purpose

Drills are intended for atraumatic preparation of the implant site for the implants in the maxilla and mandible.

The sleeves for the pOPosition for tioLogic® system are intended for precision guiding of the pOPosition for tioLogic® drills to the previously planned position using a surgical stent.

### 4. Intended users

The pOPosition for tioLogic® product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

### 5. Medical indications / contraindications

When planning implants, observe and take into account the indications and contraindications for the tioLogic® ST implants (REF 989-701-27) and tioLogic® TWINFIT implants (989-365-95) given in the respective IFUs.

#### 5.1 Indications

##### 5.1.1 pOPosition gingiva cutter

- diameter-specific and minimally invasive removal of the gingiva
- safe separation of the periosteum
- guidance and depth stop via the basic sleeve
- can be reused 15 to 20 times in hard bone

##### 5.1.2 pOPosition depth drill

- determining the planned implant length
- guidance and depth stop via the inner sleeve of the depth drill
- can be reused 15 to 20 times in hard bone

##### 5.1.3 pOPosition stepped drill

- preparing the implant site according to diameter and length
- guidance and depth stop via the inner sleeve of the stepped countersink
- collection of bone chips as autologous augmentation material
- can be reused 15 to 20 times in hard bone

##### 5.1.4 pOPosition expander

- preparing the implant site according to diameter and length
- guidance and depth stop via the basic sleeve
- can be reused 15 to 20 times in hard bone

##### 5.1.5 ADVANCED thread tap

- pre-cutting of the threads to reduce the implant insertion torque
- can be reused 15 to 20 times in hard bone

##### 5.1.7 pOPosition basic sleeves

- guidance and depth stop for gingiva cutter and expander
- insertion and fixation of the inner sleeves depth drill and stepped countersink

##### 5.1.8 pOPosition inner sleeves

- guidance and depth stop for depth drill and stepped countersink

#### 5.2 Contraindications

##### 5.2.1 pOPosition gingiva cutter

- use for implant site preparation
- use of an incorrect diameter
- use without basic sleeve
- use without internal cooling
- reuse exceeding 15 to 20 times in hard bone

##### 5.2.2 pOPosition depth drill

- use without inner sleeve depth drill
- use with inner sleeve stepped countersink
- use of an incorrect length
- use without internal cooling
- reuse exceeding 15 to 20 times in hard bone

##### 5.2.3 pOPosition stepped countersink

- use without inner sleeve stepped countersink
- use with inner sleeve depth drill
- use of an incorrect diameter
- use of an incorrect length
- use without internal cooling
- exceeding the maximum speed
- reuse exceeding 15 to 20 times in hard bone

##### 5.2.4 pOPosition expander

- use without basic sleeve
- use without internal cooling
- exceeding the maximum speed

- use of an incorrect diameter
- reuse exceeding 15 to 20 times in hard bone

#### 5.2.5 ADVANCED thread tap

- exceeding the maximum speed
- use of an incorrect diameter
- reuse exceeding 15 to 20 times in hard bone

#### 5.2.6 pOPosition basic sleeves

- positioning at the incorrect height and the wrong orientation
- failure to observe the sleeve height and the distance between the sleeve and the planned implant position
- not for multiple use

#### 5.2.7 pOPosition inner sleeves

- use of the inner sleeve depth drill with the stepped countersink
- use of the inner sleeve stepped countersink with the depth drill
- use without or with defective silicone ring
- use of an incorrect diameter
- not for multiple use

## 6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

## 7. Information on use, processing, fitting and insertion

Prior to use, dental professionals should ensure that they have carefully read and understood the full pOPosition for tioLogic®, tioLogic® ST and tioLogic® TWINFIT Instructions for use.

We strongly recommend that, prior to using Dentaurem products, dental professionals attend a Dentaurem training course to familiarize themselves with the Dentaurem products. Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaurem recommends full clinical, radiological, photographic and statistical documentation.

### Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.


Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

### Important:

**The implantologist should ensure the products cannot be aspirated during intra-oral use.**

### 7.1 Preparation for surgery

Ensure that hygiene measures are carefully followed, documented and validated throughout the surgical procedure. Prepare the treatment room, instruments and patient under these aspects. All surgical instruments required for the operation should be checked to ensure that they are complete, functional and sterile. We recommend having several implants and surgical instruments available as a precaution.

All pOPosition for tioLogic® drills, ADVANCED for tioLogic® thread taps, and the components for guided surgery, are supplied non-sterile. They should be thoroughly cleaned, disinfected, conditioned and sterilized before being used for the first time (factory new) and immediately after each use. Check the sharpness, condition, and markings of the rotary instruments after each use, as they have a limited service life. Instruments can become blunt as a result of use and cleaning. Only instruments that are sharp and free from defects should be used. The rotary instruments can be reused 15 to 20 times in hard bone, provided they are undamaged, not contaminated, and handled with appropriate care. Avoid any further reuse or the use of damaged and / or contaminated instruments. The dental professional bears responsibility. No liability is accepted if these instructions are disregarded.

## 7.2 Notes on operation planning

The diagnostic documentation and the surgical stent prepared in advance must be available for the surgical procedure with the pOPosition for tioLogic® drills.

Implants need a distance of at least 1.5 mm to the neighboring teeth. The distance between the implants must be at least 3.0 mm. Select the implant diameter and length so that there is adequate bone around the implant (min. 1.0 mm).

For technical reasons, the pOPosition for tioLogic® depth drills and the pOPosition for tioLogic® stepped countersinks are 1.0 mm longer apically than the given preparation length. This should be taken into account during diagnostics, planning and preparation.

The distance to the mandibular nerve and inferior alveolar nerve should be at least 1.5 mm!

The sleeve height (5.0 mm) and the sleeve distance to the implant shoulder (4.0 mm) must be observed and strictly adhered to.

## 7.3 General information on implant site preparation

### 7.3.1 pOPosition for tioLogic® rotary instruments

Careful preparation of the hard and soft tissue in the bone is important for optimum healing of the implant site. Avoid thermal or mechanical trauma at all costs. To this end, keep the temperature as low as possible when preparing the site and do not exceed the maximum implant insertion torque (max. 40 Ncm).

For this reason, a low-speed preparation is recommended for certain drills.

As a general rule, the speed should be reduced as the drill diameter increases, as the speed is directly related to the development of heat.

Article		Maximum torque
pOPosition gingiva cutter		max. 800 rpm
pOPosition depth drill		max. 800 rpm
pOPosition stepped countersink		max. 500 rpm
pOPosition expander		max. 500 rpm
ADVANCED thread tap		max. 10 rpm or manually

### Important:

**To prevent chip formation on the basic and inner sleeves caused by the pOPosition for tioLogic® drills, the pOPosition for tioLogic® drills should only be started, i.e. set in rotation, after they have been completely inserted into the respective guide sleeve.**

## 7.4 Implant site preparation

### 7.4.1 pOPosition for tioLogic® rotary instruments

The handpiece instruments can also be used manually as required using an adapter (max. permitted torque 40 Ncm). The instruments should be inserted rotationally secure and the fit checked. The manual hex keys and adapter can be used with the torque ratchet set at the relevant torque. The respective components should be secured with a sterile safety cord to prevent aspiration during use.

Select the required rotary instruments and their drilling sequence and depth for the preparation of the implant site depending on the bone quality. A drilling protocol for different bone qualities (soft, medium, hard) is available for the implantologist. Determination of the bone quality rests on the implantologist. Drilling should be performed using an intermittent drilling technique, in accordance with the specified drilling sequence and with sufficient cooling. The pOPosition for tioLogic® drills feature internal cooling. If additional external cooling is possible, this should be done with pre-cooled (5°C/41°F), sterile, physiological saline solution.

Before starting surgery, check that the diagnostic information is complete, including the surgical stent that was fabricated and disinfected beforehand. When inserting the disinfected surgical stent check that it is correctly positioned, e.g. on the residual dentition of the partially edentulous jaw. Positioning of the surgical stent in the patient's mouth must correspond exactly with planned positioning of the surgical stent in the planning software. The gingiva or other structures should not interfere with the correct fit of the stent. If this is the case, the gingiva must be mobilized or the structure adapted so that it is possible to fit the stent without any interference. Raising a gingival flap also makes it easier to view the operation site.

Preparation protocol – pOPosition for tioLogic®.

Adapt preparation protocol according to the indication and situation of each patient if required.  
 The initial situation is the preparation protocol  
 Adapt if required

	Illustration	Soft bone quality				
		∅ 3.3	∅ 3.7	∅ 4.2	∅ 4.8	∅ 5.5
Gingiva cutter pOPosition <sup>1</sup>		(X)	(X)	(X)	(X)	(X)
Depth drill pOPosition <sup>2</sup>		X	X	X	X	X
Stepped countersink pOPosition <sup>1</sup>		X	X	X	X	X
Expander pOPosition <sup>1</sup>		X	X	X	X	X
Thread tap ADVANCED <sup>13</sup>						

X Preparation depth in accordance with implant length


( ) Optional application  
 (taking into account the respective bone quality)

	Illustration	Medium bone quality				
		∅ 3.3	∅ 3.7	∅ 4.2	∅ 4.8	∅ 5.5
Gingiva cutter pOPosition <sup>1</sup>		(X)	(X)	(X)	(X)	(X)
Depth drill pOPosition <sup>2</sup>		X	X	X	X	X
Stepped countersink pOPosition <sup>1</sup>		X	X	X	X	X
Expander pOPosition <sup>1</sup>		X	X	X	X	X
Thread tap ADVANCED <sup>13</sup>						

	Illustration	Hard bone quality				
		∅ 3.3	∅ 3.7	∅ 4.2	∅ 4.8	∅ 5.5
Gingiva cutter pOPosition <sup>1</sup>		(X)	(X)	(X)	(X)	(X)
Depth drill pOPosition <sup>2</sup>		X	X	X	X	X
Stepped countersink pOPosition <sup>1</sup>		X	X	X	X	X
Expander pOPosition <sup>1</sup>		X	X	X	X	X
Thread tap ADVANCED <sup>13</sup>		(X)	(X)	(X)	(X)	(X)

<sup>1</sup> Exemplary illustration of rotary instruments ∅ 4.2 mm (red).  
<sup>2</sup> The insertion depth/length of the depth drill depends on the implant length.  
<sup>3</sup> Only use the ADVANCED thread tap for tioLogic® ST and tioLogic® TWINFIT implants. Thread taps must be used if implant insertion torque exceeds 40 N cm. The depth scales must be observed.

7.4.1.1 pOPosition gingiva cutters


Gingival incision pOPosition

The internally cooled gingiva cutter removes the mucosa depending on the diameter and is used for all implant lengths of a given diameter. The integrated depth stop ensures that it does not exceed the predefined insertion depth. The gingiva cutter is color-coded to indicate the diameter as well as laser-marked with the implant diameter.

As an alternative to the method for opening the mucosa described above, the gingiva can be removed minimally invasively using the pOPosition gingiva cutter. It features a color coding, which indicates the planned final implant diameter.

The internally cooled gingiva cutter is inserted in the basic sleeve and cuts through the mucosa to the bone. The stent must be removed to allow excision of the separated mucosa. The green contra angle (800 rpm) is used for drilling with internal and external cooling using sterile, cooled physiological saline solution (5°C/41°F). To allow the tip of the gingiva cutter to cool down, drill intermittently without applying pressure.

7.4.1.2 pOPosition depth drills

Inserted inner sleeve for depth drilling

Depth drilling pOPosition

The correct color-coded inner sleeves should be used for application of the depth drills. The internally cooled depth drill prepares the defined depth and direction of the implant site independent of the implant diameter. The respective inner sleeve for depth drilling should be used for application of the pOPosition for tioLogic® depth drill. This is inserted into the corresponding basic sleeve so that the 2.0 mm diameter depth drill is guided accurately. Select the appropriate depth drill for the planned length and insert it through the corresponding inner sleeve as far as the integrated depth stop. This ensures that the planned insertion length is not exceeded, even when the site is difficult to see. All pOPosition for tioLogic® depth drills feature internal cooling and a diameter of 2.0 mm. They are laser-marked with the respective implant length. For technical reasons, the depth drill is 1.0 mm longer than the given preparation length. This should be taken into account during planning and preparation. The green contra angle (800 rpm) is used for drilling with internal and external cooling using sterile, cooled physiological saline solution (5°C/41°F). To allow the tip of the depth drill to cool down, drill intermittently without applying pressure.

7.4.1.3 pOPosition stepped countersinks


Inserted inner sleeve for stepped countersinking

Stepped countersinking pOPosition

Use the correct color-coded inner sleeves for application of the stepped countersinks. The internally cooled stepped countersink widens the implant site to correspond to the implant contour. The respective inner sleeve for stepped countersinking should be used for application of the pOPosition stepped countersink. This is inserted into the corresponding basic sleeve so that the stepped countersink is guided accurately. Select the appropriate stepped countersink for the planned implant diameter and length and insert it through the corresponding inner sleeve as far as the integrated depth stop. This ensures that the planned insertion depth is not exceeded, even when the site is difficult to see. All stepped countersinks have internal cooling and a special cavity for collecting bone chips. They are laser-marked with the respective diameter and length and color-coded to indicate the diameter. For technical reasons the stepped countersink is 1.0 mm longer than the given preparation length. This should be taken into account during planning and preparation. The green contra angle (500 rpm) is used for drilling with internal and external cooling using sterile, cooled physiological saline solution (5°C/41°F). To allow the tip of the stepped countersink to cool down, drill intermittently without applying pressure.

7.4.1.4 pOPosition expanders


Expanding pOPosition

After use of the stepped countersink, the fine thread section of the implant site is prepared to the diameter of the implant using the internally cooled expander. Insert the internally cooled expander through the basic sleeve as far as the integrated depth stop. The expander is laser-marked with the respective diameter and color-coded to indicate the diameter.

The green contra angle (500 rpm) is used for drilling with internal and external cooling using sterile, cooled physiological saline solution (5°C/41°F). To allow the tip of the expander to cool down, drill intermittently without applying pressure.

7.4.1.5 ADVANCED thread taps

Depending on the bone quality, finishing the preparation with a thread tap is recommended as an option. It is available in the same diameters as the implants, and features a color-coded groove and an additional marking on the shank. When using the ADVANCED thread tap, the depth reached can be read at bone level. Remove the surgical stent for optimum visibility. Corresponding depth markings (7.0 mm, 9.0 mm, 11.0 mm, 13.0 mm


and 15.0 mm) on the thread tap indicate when the planned implant length has been attained. In addition, the ADVANCED thread tap features a hexagon chucking system for the transmission of high torques.

Manual preparation using a torque ratchet adapter or sure-grip wheel may be helpful in very narrow bone conditions. Alternatively, the thread tap can also be used with a handpiece (max 10 rpm) with the relevant attachments in the same way as for the manual technique.

After thread tapping is completed, the alveolus should be thoroughly rinsed with sterile, cooled physiological saline solution (5 °C/41°F).

#### 7.4.1.6 Implant insertion

##### 7.4.1.6.1 tioLogic® ST

Implant insertion

manual insertion key  
S-M-L tioLogic®


The implant holder and the placement aid attached to the implants are designed to ensure contact-free insertion with all indications.

Manual insertion: insert the tioLogic® ST hex key, ratchet AF 2.5 (in 3 different lengths) into the insertion post on the implant, loosen the implant from the implant holder by turning the insertion key by a ¼ turn, and insert it manually through the basic sleeve into the prepared implant site. The planned implant position is attained when the groove on the placement aid is level with the upper edge of the basic sleeve.

Handpiece insertion: insert the tioLogic® ST hex key, ISOs shank AF 2.5 (in 2 different lengths) into the insertion post on the implant, loosen the implant from the implant holder by turning the insertion key by a ¼ turn, and insert it using the handpiece through the basic sleeve into the prepared implant site. The planned implant position is attained when the groove on the placement aid is level with the upper edge of the basic sleeve.

The tioLogic® ST implant insertion key, ratchet can be used for final positioning (insertion height) and alignment (PentaStop) of the implants. To do this, insert the implant insertion key, ratchet into the implant and secure it in the implant using the supplied matching counter screw.

##### 7.4.1.6.2 tioLogic® TWINFIT

Implant insertion

manual insertion key  
S-M-L tioLogic® TWINFIT

Remove the silicone seal with the closure screw. Note: Do not tilt as the implant lies loose in the inner container. Then grasp the implant with the handpiece insertion key. The insertion key is designed to ensure contact-free insertion, regardless of the indication.

Manual insertion: the ISO shank/ratchet adapter can be used to convert the ISO shank implant insertion key into a manual instrument. To do this, fit the ISO shank/ratchet adapter onto the ISO shank implant insertion key. After making a slight turn, the insertion key snaps into place in the implant and can be inserted manually through the basic sleeve into the prepared implant site. The planned implant position is attained when the groove on the tioLogic® TWINFIT insertion key is level with the upper edge of the basic sleeve.

Handpiece insertion: the tioLogic® TWINFIT implant insertion key, ratchet can be used for final positioning (insertion height) and alignment (PentaStop) of the implants. To do this, insert the implant insertion key, ratchet into the implant and secure it in the implant using the supplied matching counter screw. The planned implant position is attained when the groove on the tioLogic® TWINFIT insertion key is level with the upper edge of the basic sleeve.

A torque of 40 Ncm should not be exceeded with any insertion procedure. Maintain a maximum motor speed of 10 rpm during handpiece insertion. Excessive torque or motor speed can damage the implant site.

The additional 5 marks on the insertion key correspond to the 5 antirotational mechanisms in the implant and allow alignment of these with regard to the subsequent prosthetic restoration. The implant position can be finely readjusted if a point in the final position is not aligned exactly to vestibular or in the direction of the alignment mark placed on the basic sleeve by the dental technician.

## 8. Technical specification

### Composition

Product	Composition
pOPosition gingiva cutter	1.4197
pOPosition depth drill	1.4197
pOPosition stepped countersink	1.4197
pOPosition expander	1.4197
ADVANCED thread tap	1.4197
pOPosition basic sleeves	3.7165
pOPosition inner sleeves depth drill	3.7165 / 1.4301 / Silicone
pOPosition inner sleeves stepped countersink	3.7165 / 1.4301 / Silicone
tioLogic® ST hex key ISO shank hex AF 2.5	1.4197 / VMQ
tioLogic® ST hex key ratchet AF 2.5	1.4197 / VMQ
tioLogic® TWINFIT implant insertion key ISO shank hex	1.4197 / 1.4310
tioLogic® ST / tioLogic® TWINFIT implant, insertion key, ratchet	1.4197 / 1.4305

Information about the composition is included in the material list for our implantology products (REF 989-801-06), or see [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu). The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaurum.com](http://www.dentaurum.com).

## 9. Adverse reactions and residual risks

Possible complications and risks during surgical treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

Chronic pain in connection with the dental implant, permanent paresthesia, dysesthesia, nerve damage, exfoliation, hyperplasia, local or systemic infections, oroantral or oronasal fistulas, alveolar bone loss in maxilla/mandible, irreversible damage to adjacent teeth, fractures to the implant, jaw, bone or restoration, aesthetical problems.

## 10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- All pOPosition for tioLogic® drills and the ADVANCED for tioLogic® thread tap may only be used in conjunction with tioLogic® ST and tioLogic® TWINFIT implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Rotary instruments – used with proper care and provided that they are not damaged or contaminated – can be reused in dense bone 15 to 20 times. Any further reuse or the use of damaged and/or contaminated instruments must be avoided. The dental professional bears responsibility. No liability is accepted if these instructions are disregarded.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

- Use of the pOPosition for tioLogic® system in combination, requires special knowledge and skills in 3D planning and the implementation of planning data in the production of a surgical stent.

#### Severe incidents:

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

### 11. Storage and shelf life

Store the pOPosition for tioLogic® drills, the ADVANCED thread taps and components for guided surgery strictly in their original packaging, in a dry, dark place and at room temperature. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

### 12. Disposal

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

### 13. Information for single use products

The pOPosition for tioLogic® basic sleeves and the pOPosition for tioLogic® inner sleeves are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reconditioning of abutments and components (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

### 14. Information for multiple use products

The pOPosition for tioLogic® drills and the ADVANCED thread taps are supplied non-sterile and are intended for multiple use. The product must be disinfected, cleaned, and sterilized before initial use and each subsequent use. pOPosition for tioLogic® drills and the ADVANCED thread tap must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

### 15. Information on disinfection / cleaning / sterilizing

The pOPosition for tioLogic® drills and the ADVANCED thread tap are supplied non-sterile. They are intended for multiple use. They must be cleaned, disinfected and sterilized before and after being used on a patient (e.g. transfer to dental laboratory). For sterilization purposes, the pOPosition for tioLogic® drills and the ADVANCED thread tap are packaged in disposable sterilization packaging intended for sterilization in accordance with ISO 11607-1/2.

### 16. Product overview and scope of delivery

For the pOPosition for tioLogic® product range, please refer to the tioLogic® ST product catalog.

For the pOPosition for tioLogic® product range, please refer to the tioLogic® TWINFIT product catalog.

### 17. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

### 18. Explanation of symbols

Please refer to the label. See chapter "Additional information".

### 19. Additional information

Additional information can be found at [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu)

Instructions for use	REF
tioLogic® ST implants	989-701-27
tioLogic® TWINFIT implants	989-365-95
tioLogic® ST Additional components	989-805-87
tioLogic® TWINFIT Additional components	989-805-88
Torque ratchet	989-719-00
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00

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The diagram illustrates the correspondence between implant diameters and guide sleeves. It is organized into three main sections: Base Sleeve (Boccola di base), Internal Pilot Drill Sleeve (Boccola interna fresa pilota), and Internal First Drill Sleeve (Boccola interna prima fresa). Each section shows five columns corresponding to implant diameters: 3.3 mm (Brown), 3.7 mm (Yellow), 4.2 mm (Red), 4.8 mm (Blue), and 5.5 mm (Green). The letters S, M, and L are placed above the columns to indicate size categories.

	S		M		L	
diametro dell'impianto	∅ 3.3 mm	∅ 3.7 mm	∅ 4.2 mm	∅ 4.8 mm	∅ 5.5 mm	
Boccola di base						
Diametro interno	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm	
Diametro esterno	5.1 mm	5.1 mm	5.6 mm	6.2 mm	6.9 mm	
Boccola interna fresa pilota						
Diametro interno	2.05 mm	2.05 mm	2.05 mm	2.05 mm	2.05 mm	
Diametro esterno	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm	
Boccola interna prima fresa						
Diametro interno	2.85 mm	3.15 mm	3.65 mm	4.25 mm	4.95 mm	
Diametro esterno	4.3 mm	4.3 mm	4.8 mm	5.4 mm	6.1 mm	

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