

Drilling protocol tioLogic® TWINFIT and tioLogic® ST.

Adapt preparation protocol to suit indication and patient if necessary.





























































- ⊘ = max. implant length
- ⑦ = min. 7.0 mm preparation depth
- ⊕ = optional application
-  = follow the instructions for use

Illustration Surgical Tray



tion and		 = max. implant length						
		 = min. 7.0 mm preparation depth						
		 = optional application						
		 = follow the instructions for use						
								
			ø 3.3	ø 3.7	ø 4.2	ø 4.8	ø 5.5	
		Marking drill						
		Depth drill						
		Surface cutter						
		Preceding stepped countersink						
		Stepped countersink						
Preparation according to bone quality			Preparation according to bone quality	Preparation according to bone quality	Preparation according to bone quality	Preparation according to bone quality	Preparation according to bone quality	Preparation according to bone quality
			D4-D3 D3-D2 D2-D1 soft medium hard	D4-D3 D3-D2 D2-D1 soft medium hard	D4-D3 D3-D2 D2-D1 soft medium hard	D4-D3 D3-D2 D2-D1 soft medium hard	D4-D3 D3-D2 D2-D1 soft medium hard	
		Expander						
		Thread tap						
			ø 3.3	ø 3.7	ø 4.2	ø 4.8	ø 5.5	

Guide sleeves are only available for tioLogic® TWINFIT.

Drilling speed.







Prepare the hard and soft bone tissue carefully to support optimal healing of the implant.
Avoid thermal or mechanical trauma at all costs. Keep the temperature during preparation of the implant site as low as possible and adhere to the maximum insertion torque (max. 40 Ncm) during insertion.

Low-speed preparation is therefore recommended for certain drills.

This can be carried out without cooling to prevent the bone chips from being washed out.

In general, the speed should be reduced as the drill diameter increases.



Product		Low-speed preparation procedure	Maximum speed
Marking drill		max. 800 rpm	max. 800 rpm
Depth drill		max. 800 rpm	max. 800 rpm
Surface cutter		50-80 rpm	max. 500 rpm
Stepped countersink		50-80 rpm	max. 500 rpm
Expander		50-80 rpm	max. 500 rpm
Thread tap		max. 10 rpm or manually	max. 10 rpm or manually