*tiologic®ST*

---

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Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

## 2. General device description

tioLogic® ST additional components are instruments and components that support and enable the planning and implementation of implants and their prosthetic restoration.

Abutment components					
Series of abutments	S	M		L	
Implant diameters in mm	3.3	3.7	4.2	4.8	5.5

X-ray reference spheres with a diameter of 5.0 mm are polymerized into the plastic template for the X-ray reference foil and used as X-ray reference points to determine the implant position. If they are positioned directly on the mucosa, the thickness of the mucosa can be calculated. The X-ray reference spheres are intended for multiple use and reconditioning.

The guide sleeve is available with a length of 6.0 mm and is incorporated into the plastic template for the radiographic or surgical stent in the prosthetically and anatomically ideal implant position and direction. This enables the precise positioning of implants. The guide sleeves are intended for single use.

The paralleling pin can be inserted following depth drilling. This is used as orientation for subsequent depth drilling and as a depth gauge of the depth drilling performed. It is available in two diameters: 1.4 mm (for pilot drills) and 2.0 mm (for depth drills). In addition, the paralleling pin with a short thread can be fixed in an implant that has already been placed, where it also serves as a guide for further implants. The paralleling pin should be secured with a sterile cord. The paralleling pin is intended for multiple use and reprocessing.

The implant hex keys, ratchet, AF 2.5, and implant insertion keys, ratchet are instruments intended for manual implant insertion, the hex keys, ISO shank, AF 2.5 are used for handpiece insertion. Using the AF 2.5 hex keys, the implant can be removed directly from the packaging via the placement aid without being touched and then inserted. Depending on the implant diameter (S, M, or L), the implants can be connected directly to the inner connection using the appropriate insertion key and inserted into the planned final position. The implant hex keys and insertion keys are designed for repeated use and reprocessing.

Series of abutments			
Series of abutments	S	M	L
tioLogic® ST Implant hex key, ISO shank, hex AF 2.5, length in mm		19.0 / 25.0	
tioLogic® ST Implant hex key, ratchet, AF 2.5, length in mm		8.0 / 13.0 / 23.0	
Implant insertion key, ratchet, length in mm		26.6	

Hex keys AF 1.3 are designed for handpiece or manual tightening or unscrewing of the corresponding screws (AF 1.3). The manual hex keys can be used with the torque ratchet or the sure-grip wheel taking into account the respective torque values. The hex key ISO shank can be used with the handpiece taking into account the respective torque. They are intended for multiple use and reprocessing.

	Length in mm
Hex key, ISO shank, AF 1.3	20.0 / 26.0
Hex key, ratchet, AF 1.3	16.0 / 26.0
Hex ball point key ISO shank for angulated screw aperture, AF 1.3	26.0

Insertion key – ratchet for manually inserting, tightening, or unscrewing various abutment variants. The insertion keys can be used with the torque ratchet or the sure-grip wheel taking into account the respective torque values. They are intended for multiple use and reprocessing.

	Length in mm
Bar/bridge/AngleFix insertion key, ratchet	15.0 / 20.0
Ball abutment insertion key, ratchet	15.0
tioLOC Insertion key, ratchet	15.5

The adapters – ISO shank / ratchet – are used for manual handling of handpieces and components (ISO shank) with the torque ratchet or the sure-grip wheel. They are intended for multiple use and reprocessing.

	Length in mm
Adapter, ISO shank, hexagon / ratchet	15.0 / 20.0

To better access hard-to-reach areas in the mouth with handpiece instruments, they can be connected to the drill extension via the ISO shank. The drill extension is intended for multiple use and reprocessing.

	Length in mm
Drill extension – ISO shank hex	21.0

If there is sufficient intraoral space, the sure-grip wheel can be used instead of the torque ratchet for manually handling manual instruments and components and it is intended for repeated use and reprocessing.

The holding key placement aid is used to fixate the placement aid when loosening the retention screw in an implant in case of unfavorable bone conditions.

A silicone O-ring is pre-mounted on the screws for impression posts, open and on the implant hex keys, AF 2.5, and prevents the screw from dropping out when positioning the impression post in the implant or fixates the placement aid with the implant during insertion on the hex key AF 2.5. The silicone O-rings are intended for single use.

The ProstheticSet is equipped with the torque ratchet and two hex keys with lengths 16.0 mm and 26.0 mm and provides safe storage and sterilization of these items.

The tioLogic® ST implant system contains specially coordinated abutments and accessories, including instruments, for the insertion and prosthetic restoration of the implants.

## 3. Intended purpose

X-ray reference spheres are intended for use in planning and diagnostics to provide a two-dimensional reference point in the orthopantomogram and to determine the mucosal thickness.

Guide sleeves are intended for two-dimensional model analysis in planning and diagnostics (distance to adjacent teeth, opposing jaw and implant alignment).

Paralleling pins are used in planning and diagnostics for orientation for additional depth drilling and as a depth gauge for depth drilling.

tioLogic® ST Implant hex keys with ratchet AF 2.5 are instruments intended for removing the tioLogic® ST implants without touching from the packaging and for manual insertion or removal from the jaw.

Instruments for screwing implants and screws for abutments and impression components in and out. tioLogic® ST Implant insertion keys are instruments intended for removing the tioLogic® ST implants without touching from the packaging and for manual insertion or removal from the jaw.

Hex keys AF 1.3 are instruments intended for the insertion or removal of screws AF 1.3 with a handpiece or by hand.

Insertion keys – ratchets are instruments intended for manually screwing in or unscrewing various types of abutments.

The adapter - ISO shank hexagon/ratchet - is intended for manual handling of the motorized instruments and components (ISO shank) with the ratchet.

Drill extensions are intended to extend rotary instruments with a handpiece connection.

The sure-grip wheel is intended for manual screwing in and out with instruments and components.

The tioLogic® ST placement aid holding key is a manual instrument for securing the placement aid when unscrewing the retaining screw in the implant.

The silicone O-rings are retaining elements for instruments for screwing in or out.

ProstheticSet is used for secure storage and sterilization of the torque ratchet and the hex key.

#### 4. Intended users

The tioLogic® ST product range is supplied exclusively to qualified personnel such as surgeons, dentists and dental technicians. It should only be used by qualified personnel who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

### 5. Medical indications / contraindications

#### 5.1 Indications

##### 5.1.1 X-ray reference spheres

- Fixation in X-ray or surgical stent as intraoral reference point for the orthopantomogram
- Determination of the mucosal thickness

##### 5.1.2 Guide sleeves

- Planning the implant position and alignment
- Guidance for depth drilling with  $\varnothing$  2.0 mm

##### 5.1.3 Paralleling pin

- Monitors drilling depth with the depth drill
- Orientation aid for further deep drilling and implants
- Determination of the mucosal thickness

##### 5.1.4 tioLogic® ST implant hex key AF 2.5 and implant insertion key

- Contact-free removal of the tioLogic® ST implants from the implant packaging
- Handpiece or manual insertion or removal of the tioLogic® ST implants into/from the prepared drill hole, taking torque specifications into account

##### 5.1.5 Hex key ISO shank hex AF 1.3

- Handpiece and manual insertion or removal of the screws with AF 1.3, taking torque specifications into account

##### 5.1.6 Insertion key, ratchet

- Manual insertion or removal of various abutment variants, taking torque specifications into account

##### 5.1.7 Adapter – ISO shank / ratchet

- Adapters for instruments with ISO shank connection for manual insertion or removal
- For use by hand, with the sure-grip wheel or the torque ratchet

##### 5.1.8 Drill extension

- Extension of instruments with ISO shank connection for handpiece or manual insertion or removal

##### 5.1.9 Sure-grip wheel

- Engaging instruments with manual ratchet connection

##### 5.1.10 Holding key for placement aid

- Fixation of the placement aid when unscrewing the retaining screw in the implant

##### 5.1.11 Silicone O-rings

- Retaining element for screws for impression posts, open, and implant hex key AF 2.5

##### 5.1.12 ProstheticSet

- Container for the secure storage of the included torque ratchet and hex keys
- Container for the sterilization of the included torque ratchet and hex keys

#### 5.2 Contraindications

- In patients with a known intolerance to individual components.
- The contact surface between the implant and abutment must not be blasted or otherwise mechanically processed
- Insufficient oral hygiene
- Bruxism

##### 5.2.1 X-ray reference spheres

- Free use in the mouth without fixation in an X-ray or surgical stent

##### 5.2.2 Guide sleeves

- Use of drills larger or smaller than  $\varnothing$  2.0 mm

##### 5.2.3 Paralleling pin

- Use as restoration

##### 5.2.4 tioLogic® ST implant hex key AF 2.5 and implant insertion key

- Use as an insertion key for other abutments and components
- Use of the manual insertion key without securing with the counter screw
- Use with incorrect torque
- Use without torque limit
- Use without the insertion key properly engaging with the implant (PentaStop)

##### 5.2.5 Hex key ISO shank hex AF 1.3

- Use as an insertion key for implants or other abutments and components
- Use with incorrect torque
- Use without torque limit

##### 5.2.6 Insertion key, ratchet

- Use as an insertion key for implants or other abutments and components
- Use with incorrect torque
- Use without torque limit

##### 5.2.7 Adapter – ISO shank / ratchet

- Use as insertion key or hex key
- Use with incorrect torque
- Use without torque limit

##### 5.2.8 Drill extension

- Use as insertion key or hex key
- Use with incorrect torque
- Use without torque limit

##### 5.2.9 Sure-grip wheel

- Use with articles without manual pick-up

##### 5.1.10 Holding key for placement aid

- Screwing in the dental implant with the placement aid

##### 5.2.11 Silicone O-rings

- Multiple use

##### 5.2.12 ProstheticSet

- Used as cleaning container
- Sterilization and storage without approved sterilization film

### 6. Patient target group

All skeletally mature patients undergoing implant treatment. Please consider the indications / contraindications.

### 7. Information on use, processing, fitting and insertion

Prior to use, the dental professional ensures that they have carefully read and understood the tioLogic® ST Instructions for use.

We strongly recommend that, prior to using Dentaaurum products, dental professionals attend a Dentaaurum training course to familiarize themselves with the Dentaaurum products.

Carry out thorough radiological, psychological and physical examinations before implantation and inform patients of the process in detail. Dentaaurum recommends full clinical, radiological, photographic and statistical documentation.

The bar, bridge and AngleFix abutments and components of the tioLogic® ST implant system can be documented, e.g. in the patient file or PatientPass, using the additional labels.

#### Precautions:

Before surgical implantation, it is important to carry out a careful patient anamnesis and, if required, a third-party anamnesis by the patient's doctor (GP/PCP) in order to clarify if 1) implantation may be difficult due to an anatomical situation, 2) there is the risk of a serious surgical problem or a general risk, 3) wound healing/osseointegration may

be difficult, or 4) good hygiene and care of the implant, abutment and restoration may worsen. Do not fit dental implants if any of the above cases or a combination hereof is severe or uncontrollable.


Prosthetic restoration of the implant must not be performed until the soft tissue has healed without irritation. Furthermore, a radiographic check-up should be performed before the prosthetic restoration.

**Important:**  
The implantologist should ensure the products cannot be aspirated during intra-oral use.

CAD/CAM data sets for 3shape™, dental wings and exocad can be downloaded free of charge from www.dentaurum.de/cadcam. The data sets include all the information needed for all types of restorations as a complete package, i.e. for printed implant models, hybrid abutments, one-piece abutments, as well as bridges and bars.


**7.1 X-ray reference spheres**

X-ray reference sphere To produce the X-ray reference foil, the X-ray reference spheres (∅ 5.0 mm) are polymerized into the plastic template and used as X-ray reference points. For this purpose, they are positioned directly on the mucosa to determine its thickness.

**7.2. Guide sleeves**

Guide sleeves are polymerized into the plastic template in the ideal implant position and alignment direction for the prosthetic restoration to produce the surgical stent. The guide sleeves are available in a lengths of 6.0 mm. If the drill that corresponds to the outer diameter of the guide sleeves is used, the guide sleeves can be pressed directly into the planned position and direction.

Model analysis for measuring the ridge height and width after initial examination can also be used for integrating the guide sleeves in the plastic template. When planning, it is important to take into account that a surface cutter, that removes bone material, may be used. During model analysis the relationship to the adjacent teeth and opposing jaw is assessed and transferred to a special sectioned model. The surgical stent is placed on the sectioned model and the implant alignment checked. If the checks on the sectioned model are positive, an orthopantomogram can be taken with the X-ray reference foil. The position, diameter, length of the implants and their alignment in relation to the adjacent teeth can be checked two-dimensionally and the position of the guide sleeves can be corrected if necessary.

**7.3 Paralleling pin**

The paralleling pin can be used following depth drilling. This is used as orientation for subsequent depth drilling and as a depth gauge of the depth drilling performed. It is available in a diameter of: 1.9 mm (for depth drills).

In addition, a paralleling pin with a short thread can be fixed in an implant that is already in place. In addition to providing orientation for further implants, the mucosal thickness can also be determined. Secure the paralleling pin with a cord.


**7.4 tioLogic® ST implant hex key AF 2.5 and implant insertion key**

The implant holder and the placement aid attached to the implant are designed to ensure contact-free insertion with all indications. After removing the implant holder from the glass ampoule, pick up the implant using the hex key AF 2.5 via the placement aid. The implant hex key is available in various lengths for handpiece and manual pick-up and insertion.

**7.4.1 Handpiece insertion:**

Insert the hex key, ISO shank AF 2.5 (in two different lengths) into the placement aid on the implant, loosen the implant from the implant holder by turning the insertion key by a ¼ turn, and insert it with a handpiece into the prepared implant site. The insertion key can be extended using a drill extension or converted into a manual instrument using the ISO shank/ratchet adapter.

**7.4.2 Manual insertion:**

Insert the hex key, ratchet AF 2.5 (in three different lengths) into the placement aid on the implant, loosen the implant from the implant holder by turning the insertion key by a ¼ turn, and insert it manually into the prepared implant site.

The implant insertion key, ratchet can be used for final positioning (insertion height) and alignment (PentaStop) of the implants. To do this, insert the implant insertion key, ratchet into the implant and secure it in the implant using the supplied matching counter screw.


A torque of 40 Ncm should not be exceeded with any insertion procedure. Maintain a maximum motor speed of 10 rpm during handpiece insertion. Excessive torque or motor speed can damage the implant site.

The five markings on the implant hex key and the implant insertion keys correspond to the five antirotational mechanisms on the implant and allow alignment of these with regard to the subsequent prosthetic restoration. Insert the implant into the bone as far as the lower edge of the polished neck area, i.e. its final position is slightly transcrestal (0.3 mm).

The AF 2.5 implant hex keys can be removed from the placement aid after implant insertion (in cases of reduced horizontal bone conditions, remove the hex key with the necessary caution). Before removing the implant insertion key, ratchet, first loosen the counter screw.

Tightening torque	
Implant in the mouth depending on bone quality	max. 40 Ncm

**7.5 Hex key AF 1.3**

The hex keys (AF 1.3) are provided in various lengths for different intra-oral and extra-oral situations and are intended for handpiece or manual tightening or unscrewing of the corresponding screws (AF 1.3). The manual hex keys can be used

with the torque ratchet or the sure-grip wheel taking into account the respective torque values of the different screws (AF 1.3). They are intended for multiple use and reprocessing.

The conically tapered working end of the hex keys (AF 1.3) creates a clamping effect when engaging the screws, which, when used correctly, ensures secure positioning of the screw.

The hex ball point key ISO shank is specially designed for tightening or unscrewing in angulated screw apertures. The screw aperture can be inclined at an angle of 20° to the implant axis, allowing discrete emergence in the palatal area, particularly in the anterior region. For angulated screw apertures, use the AnOTite screw for angulated screw apertures L 8.5 mm and the hex ball point key AF 1.3.

The ball-point and hex keys with ISO shank connection are either handpiece-operated or can be converted into a manual hex key using the ISO shank/ratchet adapter. In both cases, observe the appropriate torques for the different screws (AF 1.3).


**7.6 Insertion key, ratchet**

Insertion key – ratchet for manually inserting, tightening, or unscrewing various abutment variants.

The insertion keys can be used with the torque ratchet or the sure-grip wheel taking into account the respective torque values.

They are intended for multiple use and reprocessing.

7.6.1 The bar/bridge/AngleFix insertion key, ratchet are specifically designed for tightening or unscrewing tioLogic® ST bar, bridge and AngleFix abutments and are available in two lengths.

Insert the bar/bridge/AngleFix insertion key on the bar, bridge or AngleFix abutment and secure it to the respective abutment using the screw provided. Longitudinal grooves in the bar/bridge/AngleFix insertion key enable the abutments to be screwed in. The bar/bridge/AngleFix abutments are inserted taking into account the corresponding torques. The bar/bridge/AngleFix insertion key is designed for repeated use and reprocessing.

#### Tightening torque

Ball abutment on the model / in the mouth	35 Ncm
---	--------

Insertion key ball abutment

7.6.2 The ball abutment insertion key, ratchet is designed specifically for tightening or unscrewing tioLogic® ST ball abutments.

The ball abutment insertion key is inserted on the ball abutment. The hexagon inside the ball abutment insertion key allows the ball abutments to be screwed in. The ball abutments are inserted taking into account the corresponding torques. The ball abutment insertion key is designed for repeated use and reprocessing.

#### Tightening torque

Ball abutment on the model: / in the mouth	35 Ncm
--	--------

Insertion key tioLOC abutment - ratchet


7.6.3 tioLOC The insertion key, ratchet is designed specifically for tightening or unscrewing tioLOC abutments.

The triangular section at the working end of the insertion key fits precisely into the recess provided in the tioLOC abutments. The tioLOC abutments are inserted taking into account the corresponding torques. The tioLOC insertion key is designed for repeated use and reprocessing.

#### Tightening torque

tioLOC Abutment on the model / in the mouth	30 Ncm
---	--------

## 7.7 Adapter ISO shank / ratchet

The adapters – ISO shank / ratchet – are used for manual handling of handpieces and components (ISO shank) with the torque ratchet or the sure-grip wheel. Depending on the available space and the situation, two different adapter lengths can be selected. Instruments with a handpiece connection are coupled via the ISO shank to the ISO shank/ratchet adapter and locked in place by the spring on the adapter. They can then be used with the torque ratchet or sure-grip wheel, taking into account the respective torque values for the different screws (AF 1.3). They are intended for multiple use and reprocessing.

## 7.8 Drill extension

To better access hard-to-reach areas in the mouth with handpiece instruments, they can be connected to the drill extension via the ISO shank. It is intended for multiple use and reprocessing.

## 7.9 Sure-grip wheel

If there is sufficient intra-oral space, the sure-grip wheel can be used instead of the torque ratchet for the manual handling of manual instruments and components. The sure-grip wheel is fitted over the manual connection of the instrument. The internal snap ring locks the respective tool in place. The sure-grip wheel is intended for multiple use and reprocessing.

## 7.10 Holding key for placement aid

The screw in the placement aid is loosened with the hex key AF 1.3 and the insertion aid removed. If the implant turns when the screw is loosened (e.g. with reduced horizontal bone conditions), the holding key for the placement aid is to be used to provide rotational safety.

## 7.11 Silicone O-rings

A silicone O-ring is pre-mounted on the screws for impression posts, open and prevents the screw from dropping out when positioning the impression post in the implant. The silicone O-rings are intended for single use.

## 7.12 ProstheticSet

The ProstheticSet is a tray equipped with the torque ratchet and two hex keys with lengths 16.0 mm and 26.0 mm. Cleaning and disinfection were performed in disassembled condition.

The ProstheticSet is used for the safe storage of the instruments and can be sterilized together with these items.

## 8. Technical specification

### Composition

Product	Composition
X-ray reference spheres	1.4301
Guide sleeves	3.7065
Paralleling pin	3.7165
tioLogic® ST implant hex key, AF 2.5 handpiece	1.4305
tioLogic® ST implant hex key, AF 2.5 manual	1.4305
tioLogic® ST implant, insertion key, manual	1.4197 / 1.4305
Hex key 1.3	1.4197
Hex ball point key for angulated screw aperture AF 1.3, ISO shank	1.4108
Adapter, ISO shank, hexagon / ratchet	1.4197 / 1.4310
Insertion key tioLOC abutment	1.4197
Drill extension	1.4197 / 1.4310
Sure-grip wheel	1.4305 / 1.4310
Silicone O-rings	Silicone
ProstheticSet	PPSU

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see [www.dentaurum.com/ifu](http://www.dentaurum.com/ifu). The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaurum.com](http://www.dentaurum.com).

## 9. Adverse reactions and residual risks

Possible complications and risks during prosthetic treatment (temporary adverse reactions) are:

Pain, swellings, difficulties speaking, inflammation of the gums.

Possible complications and risks during prosthetic treatment (permanent adverse reactions) are:

fractured dentures, aesthetic problems.

## 10. Safety information and precautions

- Do not use the product if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. Therefore, secure all articles used intra-orally against swallowing and/or aspiration.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- All tioLogic® ST prosthetic components may only be used in combination with tioLogic® ST implants! Use of third-party components can negatively impact the function and safety. In such cases, Dentaurum GmbH & Co. KG accepts no liability for replacement products or guarantees.
- Observe the indications for the individual products. Consider all indications that are not listed explicitly to be contraindications.
- Disinfection, cleaning and sterilizing should not be done in the original packaging.
- Products must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.
- Technical advice (oral and written) is based on the scientific and technical knowledge available at the time the product is put on the market. It does not release the dental professional from the responsibility of personally checking the suitability of the products for the intended indications and applications. Advice is only given as a non-binding recommendation and cannot be assumed to provide any form of assurance, liability or guarantee.

**Severe incidents:**

Report all severe incidents arising from the use of the product to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.

Severe incidents are identified as malfunctions and a drop in performance or deterioration of the properties of a medical device which resulted or may have resulted in a direct or indirect temporary or permanent deterioration in a patient's state of health, in the death of a patient, a dental professional or other persons or as posing a serious danger to public health.

**11. Storage and shelf life**

Store the tioLogic® ST implants, abutments and components always in a dry, dark place at room temperature in the original packaging. If the packaging shows damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.

**12. Disposal**

Adhere strictly to the rules that apply for your region and the details outlined in the safety data sheets, if available. Products to be disposed of should be treated as surgical waste and be decontaminated in compliance with the relevant regulations for medical waste disposal.

**13. Information for single use products**

Guide sleeves and silicone-O-rings are supplied non-sterile and are intended for single use on a single patient only. The product must be disinfected, cleaned and sterilized before first use.

Reprocessing of guide sleeves and silicone-O-rings (recycling) as well as reuse on patients are not permitted as safe reconditioning and/or functional safety cannot be guaranteed and may lead to a loss in mechanical, chemical and/or biological properties. Reuse may cause infection. Do not use implants after the expiry date stated on the label.

**14. Information for multiple use products**

X-ray reference spheres, paralleling pins, implant hex key AF 2.5, implant insertion key, hex key AF 1.3, insertion key (ratchet), ISO shank/ratchet adapter, drill extension, sure-grip wheel, and the ProstheticSet are supplied non-sterile and are intended for multiple use. The respective product must be disinfected, cleaned and sterilized before being used initially and before being used again. All reprocessed items must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels.

**15. Information on disinfection / cleaning / sterilizing**

All additional components are supplied non-sterile. They are for single use on one patient only, unless explicitly intended for multiple use. They must be cleaned and disinfected before and after being used on a patient (e.g. transfer to dental laboratory). Dentaurem generally recommends that the product be sterilized in addition, as the use of non-sterile components may cause infections. Only sterilized prosthetic components may be used if bleeding occurs.

**16. Product overview and scope of delivery**

Please refer to the tioLogic® ST product catalog for the tioLogic® ST product range.

**17. Quality information**

Dentaurem ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurem has no influence on the processing on site.

**18. Explanation of symbols**

Please refer to the label. See chapter "Additional information".

**19. Information on MRT safety**

The products within these Instructions for use are MRT safe.

**20. Additional information**

Additional information can be found at [www.dentaurem.com/ifu](http://www.dentaurem.com/ifu)

Instructions for use	REF
Screws	989-804-98
Torque ratchet	989-719-00
tioLogic® ST bar abutments and components	989-802-38
tioLogic® ST bridge abutments and components	989-801-91
tioLogic® ST AngleFix abutments and components	989-801-89
Preparation instructions for the implant system components	989-801-07
Explanation of symbols used on the label	989-313-00
MRI safety of Dentaurem products	989-804-43

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