remanium[®]
LMF140 powder 

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15–45 µm

Powder

CoCrMo dental alloy based on cobalt for firing and partial denture technique, Type 5 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 / DIN EN ISO 5832-4 for processing with laser melting machines

1. Composition (% by mass)

Co	Cr	Mo
65.5	28.5	6.0

Other elements <1%: Mn, Fe, W, Si, N. This product contains cobalt. It contains no nickel, beryllium, lead, and cadmium. This alloy meets the requirements for the chemical composition as specified in DIN EN ISO 5832-4 and ASTM F75 as well as ASTM F3213.

2. Technical data

Yield strength*	$R_{p0.2}$	800 MPa
Tensile strength*	R_m	1.250 MPa
Hardness*	H	395 HV10
Elongation at rupture*	A_5	12 %
Modulus of elasticity*	E	235 GPa
Density	ρ	8.4 g/cm ³
Solidus temperature	T_s	1.390 °C
Liquidus temperature	T_L	1.430 °C
CTE (25 °C - 500 °C)	α	14.8 x 10 ⁻⁶ K ⁻¹

The technical data represent typical values that are achieved when using remanium® LMF140 powder with the heat treatment recommended here. A laser melting machine of the type Zonelab 2Create was used for the production of test specimens. The layer thickness was 30 µm.

* Variations of ± 10 % are possible depending on the batch.

3. Intended purpose

Non-precious alloy for the manufacture of medical products by means of an additive powder bed fusion process, e.g. selective laser melting.

4. Indication

remanium® powders for additive manufacturing are cobalt-based dental alloys. remanium® LMF140 powder is particularly suitable for the manufacture of fixed and removable dental restorations, appliances and metal-ceramic frameworks.

5. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction must be in line with dental regulations in order to fulfil mechanical requirements. The finished product must, however, have a material thickness of at least 0.3 mm.

6. Processing in the powder bed system

Before processing in the powder bed system begins, it must be clarified that the used system is suitable for processing this material. Product-optimized and verified process parameters are available for the machine manufacturers listed under "Compatibility" at <https://www.dentaurum.de/lasermelting>.

You can find more information in the download section at <https://www.dentaurum.de/lasermelting>.

If you would like to work with a different supplier, please contact your machine manufacturer in order to clarify this or contact us. We will be pleased to help you check the system requirements.

7. Post treatment

Powder residues on the components manufactured by the additive manufacturing process must be removed with due care, taking safety precautions into consideration.

8. Stress relief heat treatment

The process outlined below describes stress relief heat treatment for products manufactured by additive manufacturing technology and ensures that the material properties are as indicated and guarantee that the components are of high quality and lasting durability. If the heat treatment program is altered, results may differ. This, in turn, may void the guarantee. Stress relief heat treatment optimizes the material structure of the CoCrMo alloy. Inner stress is relieved and the accuracy of fit of the dental component is improved.

We recommend that heat treatment of the components is carried out under a protective atmosphere on a build plate in a suitable furnace. When using a gas box, observe the gas flow rates recommended by the manufacturer. Argon is recommended as shielding gas.

Heat treatment:

1. Preheat the furnace to 800 °C and then place the build plate with the supported components into the furnace.
2. Hold the temperature at 800 °C for 40 minutes.
3. Then, remove the build plate and allow to cool down at room temperature.
4. Processing of the components can continue once they have cooled down.

9. Post treatment

The components can be removed from the build plate using a separating disc or band saw. Reaction layers as well as support residues can be removed and surfaces smoothened using a finely-toothed hard metal bur. The surface can be conditioned either by polishing manually or mechanically. When using a polishing unit with a dry suspension, the surfaces are processed in preparation using blasting abrasive (REF 128-016-00) and, if required, using finely-toothed hard metal burs.

10. Powder processing

Surplus powder can, for example, be passed through a vibratory sieve shaker from the AS 200 series (Retsch) with a mesh of approx. 63 µm. For compact screening machines with ultrasound stimulation and outlet tube, a screening mesh width of 63 µm is recommended, e.g. CSM 500 (Assonic). After a visual inspection confirms its good condition, the separated powder can be further processed in the powder melting system or stored in its original container, kept dry and tightly sealed.

11. Welding

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Suitable welding wires are available as filler materials for all remanium® alloys, for example CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

12. Soldering

Soldering should be avoided if possible to reduce the number of materials. If, however, soldering is necessary, please use a suitable solder and flux which are appropriate for the composition and melting range of the alloy used, for example CoCrMo-Sold 1 (REF 102-306-00) in connection with rema® Flux 1 (REF 102-304-00).

13. Ceramic veneering

remanium® bonding alloys can be veneered with all suitable ceramic masses, e.g. ceraMotion® Me (Dentaurum), as long as the CTE is suitable for the alloy. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

14. Contraindications and adverse reactions

Signs of intolerance to non-precious remanium® alloys are extremely rare if the manufacturer's instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reports of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur. Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California, elements of this alloy are considered to be carcinogenic.

15. Warnings and precautions

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. It is possible that dust may be generated when the package is opened or the remanium® powder is transferred. Metal dust is combustible. Any type of ignition source should be kept at a distance. Metal fire powder extinguishers should always be positioned in close proximity and be in working order. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn.

Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes.

**Please observe the relevant safety data sheet
(see www.dentaurum.com)!**

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence available on the safety or efficacy of treatment of pregnant women, nursing mothers or children.

16. Reuse

The alloy may be changed in its composition and properties if reused. In principle, we recommend using new material to manufacture high-quality dental restorations.

17. Storage conditions

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the remanium® alloys. remanium® powder should be stored in a dry place in its original packaging with the lid closed tightly.

18. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

19. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, inform the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident.

The SSCP is available at <https://ec.europa.eu/tools/eudamed> and www.dentaurum.com.

These instructions for use contain the main steps for processing the non-precious remanium® alloys. You can find more information in the download section at <https://www.dentaurum.de/lasermelting>.

For more information on our products and services, visit www.dentaurum.com. You will find these Instructions for use and information on the current status at www.dentaurum.com/ifu.

Should further questions arise on laser-welding and processing titanium, please contact your local representative.

20. Explanation of symbols used on the label

Please refer to the label. Additional information can be found on the internet at www.dentaurum.com/ifu (Explanation of symbols used on the label REF 989-313-00).

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