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***rematitan***® Finishing kit

**Coffret de finition *rematitan***®

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**EN Warranty**

Dentaurum guarantees faultless quality of the products manufactured by us. The information in these instructions for use is based upon our own experiences. The user himself is responsible for the correct processing or use of the products. We cannot be held responsible for failures as we have no influence on the processing. However, any claim for compensation only applies to the commercial value of our products.

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REF 135-500-00

### 1. Safety Notes

- always wear goggles when preparing or finishing titanium restorations
- always switch on the exhaust system
- do not exceed the maximum rpm of the finishing instruments

### 2. Description

Pure titanium is a soft, tough material that requires a special finishing and polishing process. The finishing kit contains all the materials necessary for the efficient finishing and polishing of titanium for dental purposes. The sequence of use and essential features of the instruments are described.

### 3. General remarks

- finishing tools used for titanium should not be used for other alloys
- excess material on grinding tools should be avoided. This means that, for example, diamond grinders are not suitable for this purpose
- localized overheating of the workpiece must be avoided. For this reason, great care is required when using rubber polishing materials
- always work with a low pressure and low rotation speed

### 4. Preparation requirements

As titanium castings are always covered with a skin, this must be carefully removed before the actual finishing operation can start. This is done by using a sand-blasting unit with an appropriate size grain of sand.

**Important:** Sparks are always created when cast titanium is being finished.

**For rematitan® Plus investment material:**

Model casting: sand-blasting material  $\text{Al}_2\text{O}_3$  (< 250  $\mu\text{m}$ )

Crowns and bridges: sand-blasting material  $\text{Al}_2\text{O}_3$  (< 125  $\mu\text{m}$ )

**Important: Do not damage the edges of the crown – use low pressure!**

**For Trinell investment material:**

Sand-blasting material  $\text{Al}_2\text{O}_3$  (30–125  $\mu\text{m}$ )

For very fine parts (inlays) careful blasting with polishing beads is sufficient.

## 5. Procedure for finishing titanium

Observe the notes given on page 9 in the correct sequence. This assures that an excellent polish can be obtained at very little expense.

**Preparation of ceramic-coated surfaces:**

- use only hard-alloy cutters
- blast carefully using  $\text{Al}_2\text{O}_3$  (125–250  $\mu\text{m}$ ) at a pressure of 2–3 bars

**Important:** In order to allow the titanium oxide layer to build up, wait 10 minutes before placing the finished piece in the Ultra Sonic or Steam Cleaner.

In addition to the component numbers, remember to state the order numbers for spare parts as given on page 10.

## Finishing

Operations	Finishing material	Remarks
<b>Removal of sprue runners</b>	ST cutting disk (large) (1) STM cutting disk (thick) (2) TX cutting disk (thin) (3)	– thick sprue runners (quick cutter) or STM (hand tool) – thin runners TX cutting disks – avoid tilting, cool workpieces with water!
<b>Finishing (rough)</b>	Hard-alloy cutting, mini (4) Hard-alloy cutting, midi (5) Hard-alloy cutting, maxi (6) Hard-alloy cutting, maxi plus (7)	– maximum rotation speed 10 000 min <sup>-1</sup> – low pressure
<b>Finishing (fine)</b>	Aloxin grinder blue, small, B (8) Aloxin grinder blue, large, C (9)	– circular grind on workpiece – Aloxin grinder not suitable for surfaces to be coated with ceramics!
<b>Grinding (fine)</b>	Emery cloth -500- (10)	– recommended before rubber polishing
<b>Rubber polishing</b>	Rubber disks (grey) (11) Rubber disks (red) (12) Rubber bobs (red) (13) Rubber rollers (red) (14)	– use sparingly and without pressure; low rotation speed (-> danger of overheating) – grey= rough polish, red = fine polish (align rollers before first use)
<b>Polishing</b>	Large polishing brush (polishing motor) (15) Black polishing brushes (16) Small brushes (17) Tiger brilliant, polishing paste (18)	– change direction of polishing frequently – for high lustre add no more polishing paste

Components supplied with finishing kit			Replacement parts supplied	
(1)	ST cutting disk, $\varnothing$ 40 mm	1 piece	10 pieces	REF 130-100-00
(2)	STM cutting disk, mini, $\varnothing$ 25 mm	2 pieces	20 pieces	REF 130-110-00
(3)	TX cutting disk, $\varnothing$ 22 mm	5 pieces	100 pieces	REF 130-512-00
(4)	Ti hard-alloy cutter, mini	1 piece	5 pieces	REF 123-610-00
(5)	Ti hard-alloy cutter, midi	1 piece	1 piece	REF 123-611-00
(6)	Ti hard-alloy cutter, maxi	1 piece	1 piece	REF 123-612-00
(7)	Ti hard-alloy cutter, maxi plus	1 piece	1 piece	REF 123-613-00
(8)	Aloxin grinder, B, blue	1 piece	12 pieces	REF 135-852-00
(9)	Aloxin grinder, C, blue	1 piece	12 pieces	REF 135-853-00
(10)	Emery cloth -500-	70 cm	25 m	REF 140-016-00
(11)	Rubber polisher grey, disk	4 pieces	100 pieces	REF 138-102-00
(12)	Rubber polisher red, disk	4 pieces	100 pieces	REF 138-601-00
(13)	Rubber polisher red, bob	4 pieces	100 pieces	REF 138-603-00
(14)	Rubber polisher, roller	4 pieces	100 pieces	REF 138-602-00
(15)	Polishing brushes, $\varnothing$ 20 mm, black	4 pieces	10 pieces	REF 141-800-00
(16)	Polishing brushes, $\varnothing$ 60 mm, Chunkin bristle, 2 rows	1 piece	10 pieces	REF 141-711-00
(17)	Small brushes, black	1 piece	10 pieces	REF 141-810-00
(18)	Tiger brilliant, polishing paste	60 g	400 g	REF 190-350-00
(19)	Mandrels for rollers	1 piece	10 pieces	REF 139-100-00
(20)	Mandrels for disks and bobs	6 pieces	10 pieces	REF 139-300-00
(21)	Mandrels for emery paper	1 piece	10 pieces	REF 139-500-00

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