AIRVAC 2002

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo

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Contents

1.	Declaration of conformity	15
2.	Safety instructions	16
2.1	Approved use	16
2.2	Symbols indicating risk	16
2.3	Guarantee and liability	16
2.4	Obligations of authorised user	17
2.5	Obligations of personnel	17
2.6	General information for safe use of the unit	17
2.7	High risks	17
2.8	Disposal	17
2.9	Packing	18
2.10	Package contents	18
3.	Function and description of the unit	18
3.1	Function	18
4.	Installation	18
4.1	Positioning the unit	18
4.2	Mounting the unit on the wall	19–20
4.3	Connecting the unit to the power supply	20
5.	Operating the AIRVAC 2002	21
5.1	General information	21
5.2	Procedure for mixing plaster or investment	21–22
6.	Cleaning and maintenance	23
6.1	Cleaning and maintenance of the mixing mechanism and bowl flange	23
6.2	Cleaning the unit	23
6.3	Cleaning and changing the sinter filter	23
6.4	Changing the fuse	23
7.	Troubleshooting	24
8.	Spare parts	24
9.	Accessories	24
10.	Block circuit diagram	25
11.	Technical data	25

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hereby declares that the design and construction of the laboratory equipment described below, including the version marketed by us, comply with the basic regulations governing safety and health as stated in the EEC Guidelines. This declaration will become invalid if the laboratory equipment is modified or altered in any way without our prior consent.

Description of unit:	AIRVAC 2002 Vacuum mixing unit REF 095-072-00
Start with Unit No.:	114-001 115-001
EEC guidelines:	98/37/EC Guideline for machine 73/23/EEC Electrical equipment used within certain voltage limits 93/68/EEC Modified Guideline 89/336/EEC Electromagnetic compatibility 91/263/EEC Modified Guideline 92/31/EEC Modified Guideline 93/68/EEC Modified Guideline 93/97/EEC Modified Guideline
Applied unified standards:	EN 292-1 and -2 EN 294 EN 61010-1+2-051 EN 50081-1 EN 6100-3-2+3 EN 61000-6-2

Date and manufacturers signature: 01.01.2010
Signatory:

- i.V. Dipl.Ing. (FH) K. Merkle -
Production Manager Mechanic

Printing Date: 05.11.10

2. Safety instructions

Caution:

Read through the instructions for installing and operating the unit carefully before switching it on.

2.1 Approved use

The **AIRVAC 2002** vacuum mixer is intended solely for homogeneous mixing of investment, plaster or silicone used for fabricating models in the dental laboratory. The unit should only be used for the purpose for which it is intended. Dentaureum GmbH & Co. KG, is not liable for damages resulting from any other use of the unit. Approved use also includes adhering to the instructions for use as well as servicing and carrying out maintenance work at regular intervals.

2.2 Symbols indicating risk

The following symbols are used in the instructions to indicate risks:

Warning:

Indicates a potential danger to the life and health of the person. Disregarding this sign can critically affect a person's health and may even result in life-threatening injuries.

Caution:

Indicates a potentially dangerous situation. Disregarding this sign can result in minor injury or damage to the equipment.

Note: Indicates general information on operating the unit.

2.3 Guarantee and liability

Our general conditions of sale and delivery apply. Guarantee and liability claims for damages to persons or equipment are invalid if the damage is found to have resulted from any of the following:

- improper installation, operation, assembly or maintenance of the unit
- use of the unit other than for the purpose intended
- operation with faulty safety equipment or improperly installed or non-functioning safety and protective precautions
- disregarding the guidance in the instructions for use relating to transport, storage, assembly, operation and maintenance of the unit
- inadequate monitoring of parts liable to wear and tear
- unauthorised structural alterations to the unit
- repairs which have not been carried out properly.

2.4 Obligations of authorised user

The authorised user is obliged to ensure that all those who operate the unit

- Are familiar with safety regulations and accident prevention and have been instructed on how to operate the unit.
- Have read and understood the safety precautions and instructions for use.
- Have been instructed in the current regulations for the prevention of accidents.

2.5 Obligations of personnel

Before commencing work, all those who operate the unit undertake

- to adhere to the basic regulations regarding safety at work
- to sign that they have read and understood the safety precautions and instructions for use.

2.6 General information for safe use of the unit

- Keep the working area clean. Untidiness at the workplace increases the risk of accidents.
- For reasons of personal safety, use only accessories and materials recommended in the instructions.
- Always switch off the unit and disconnect it from the mains before undertaking any maintenance work.

2.7 Additional risks

- Do not reach into the unit when the mixing mechanism is rotating.
- Always use original Dentaurem bowls for mixing.

Caution:

The unit may excessively heat up if more than 25 consecutive mixing cycles are performed within one hour. In this case it is recommended to allow the machine to cool for 3 hours.

2.8 Disposal

Note:

This symbol indicates, products marked with it should not be disposed off together with garbage of households. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available through Dentaurem or the dental trade.

2.9 Packing

The packing provides only limited protection against wet, heat and any mechanical force. The packing should be disposed of in the proper manner.

2.10 Package contents

- | | | |
|---|--|----------------|
| 1 | AIRVAC 2002 vacuum mixer – 230 V | REF 095-072-00 |
| 1 | Mixing bowl, medium | REF 095-063-00 |
| 1 | Instructions for use and declaration of conformity | |
| 3 | Fixing screws for wall mounting | |
| 1 | Support to attach the AIRVAC 2002 to the wall | |
| 1 | Screw for the support | |
| 1 | Drilling template | |

3. Function and description of the unit

3.1 Function

The **AIRVAC 2002** is a universal vacuum mixer for use in dental technology. It can be used for mixing any standard investment, plaster or silicone. A feature of the integrated industrial motor is its very high torque. The high-performance vacuum pump incorporated in the unit is completely maintenance-free.

Other features of the **AIRVAC 2002**:

- very easy to operate
- precise positioning of the mixing shaft ensures the mixing motor has a long life-span
- compact
- automatic regulation of the vacuum provides increased reliability
- easy to service

4. Installation

4.1 Positioning the unit

The unit should be installed in a place which is clean and as free from dust as possible. The wall unit should be mounted on a straight wall. The standing unit should be placed on a solid, clean base.

The unit should be placed on a sturdy table, if it is to be used as a standing unit with the unit base (REF 095-075-00).

Caution:

Furnaces or other equipment, which emit heat, should not be placed next to the the **AIRVAC 2002**.

4.2 Mounting the unit on the wall

Use a spirit level and the drilling template supplied to ensure that the mounting holes are in the correct position:

- The mounting position on the wall should be dry and able to withstand the weight of the unit.
- Position the drilling template and mark the three given drill holes on the wall with a pencil through the template. The position can be checked with a spirit level.
- Drill holes at the three marks using a \varnothing 6 mm masonry drill. Insert the rawlplugs and screw in the two upper screws so that the heads still protrude approx. 2–3 mm (illus. 1).



Illustration 1: Screwing the three fixing screws into the wall

Caution:

Plasterboard rawlplugs should be used to mount the unit onto plaster or plasterboard walls. Check the load capacity of the rawlplug!

- Screw the support to the wall loosely so that it stays still movable.
- Set the unit onto the screws and pull it downwards carefully (see picture 2).
- Line up the support according to the drill hole in the bottom of the **AIRVAC 2002** and tighten both screws (to the wall and to the bottom of the unit).

1. set the unit onto the screws

2. and pull downwards carefully

3. screw tidely the support with the bottom of the unit

Illustration 2: Mounting the **AIRVAC 2002** onto the fixing screws

Caution:

Before operating the unit, check that the rawlplugs and screws are tight and the vacuum mixer is securely mounted on the wall with the support otherwise rotations during the operating process could remove the unit from the wall!

4.3 Connecting the unit to the power supply

The electricity specifications on the model plate of the unit should be the same as the mains voltage.

- The unit should be connected to a 230 V socket.
- The socket should be fitted with a 10 or 16 amp fuse.

Caution:

All electrical work on the socket or electricity supply should only be carried out by a qualified electrician.

5. Operating the AIRVAC 2002

5.1 General information

In the Airvac system, the bowl and mixing mechanism are connected to the base of the unit mechanically by turning them. Air is pumped out directly through lid of the bowl and through the underside of the unit.

Sinter filters (illus. 3, Pos. 2+3) on the flange of the bowl prevent dirt getting into the vacuum pump. There is an additional fine filter inside the unit.

The vacuum is permanently monitored through a specially drilled sensor hole (illus. 3, Pos. 3) by the integrated microelectronics.

There is a choice of three sizes of mixing bowl with different mixing mechanisms. These can be attached directly to the **AIRVAC 2002** without an adaptor or support. The bowls of earlier versions of the Airvac (REF 095-060-00 to REF 095-071-00) are compatible with the **AIRVAC 2002**.

Caution:

Separate bowls and mixers should be used for investment and plaster! Disregarding this warning can result in problems with the setting time of the plaster or investment.

- Pos. 2
Sinter filter
- Pos. 3
Drilled sensor hole
with Sinter filter
- Pos. 1
Multiple unit with socket,
mains switch and fuse holder


Illustration 3: Underside of the unit with flange for the bowl

5.2 Procedure for mixing plaster or investment

- Switch the unit on underneath using the ON/OFF switch (illus. 3, Pos. 1). The green POWER indicator (illus. 4, Pos. 4) lights up. The last mixing time set is displayed on the time scale indicator (illus. 4, Pos. 7)
- Select the required mixing time using the time setting keys (illus. 4, Pos. 9).
- Then fill the bowl with liquid and powder.

i Note: add the liquid, then the powder.

max. filling level in grams of powder incl. liquid	max. amount of investment e.g. Rema® Exakt	max. amount of plaster or fine-grained investment e.g. Castorit®-super
large bowl (REF 095-062-00)	900 g powder	700 g powder
medium bowl (REF 095-063-00)	500 g powder	400 g powder
small bowl (REF 095-064-00)	120 g powder	100 g powder

The maximum filling levels are marked on each of the bowls.

- Spatulate the mixture briefly until it is completely moist.
- Wipe and moisten the marginal seal of the bowl.
- Put on the appropriate mixer and turn it slightly to connect it to the underside of the unit.


Illustration 4: AIRVAC 2002 control panel

- When the start key (illus. 4, Pos. 8) is pressed, the unit starts extracting the air from the bowl. The yellow START indicator (illus. 4, Pos. 6) lights up. The mixer starts automatically after a few seconds and mixes the contents, while the vacuum pump continues to extract air from the bowl.
- An integrated pressure sensor automatically monitors the vacuum during mixing. If the vacuum is too low, the red VACLOW indicator (illus. 4, Pos.5) lights up.

i Note: If the mixer does not switch on after approx. 5 seconds, press the bowl more firmly to the unit.

- The vacuum pump and motor automatically switch off at the end of the programmed mixing time and the vacuum in the mixing bowl is released. The yellow START indicator goes out when mixing is complete.
- Release the bowl and mixer from the unit by turning them to the left (anti-clockwise). The mixture is now ready to use.
- Clean the mixer and bowl under running water immediately after use.

6. Cleaning and maintenance

Warning:

Always switch off the unit and disconnect it from the mains during any servicing and maintenance work. Never work alone. There should always be two people with electrical experience present when servicing and repairs are being carried out.

6.1 Cleaning and maintenance of the mixing mechanism and bowl flange

The following parts should be kept clean to achieve a good vacuum:

- vacuum passageway in the mixer lid
- sinter filter on the flange of the bowl
- seal on the mixer lid and connection for the bowl on the unit. The rubber seals can be removed for cleaning under running water
- all contact surfaces of the O-rings and bowl
- the notches on the edge of the bowl. These are used to lock the bowl in position in the support and must be kept clean
- Though the bowls are made from unbreakable plastic, extreme loading of the mixing bowls should be avoided.

6.2 Cleaning the unit

Wipe the unit occasionally with a damp cloth. Do not use abrasive cleaning agents or tools.

Caution:

Never use compressed air, water or solvent to clean the unit. Never steam clean the unit.

6.3 Cleaning and changing the sinter filters

The sinter filters in the flange of the bowl (illus. 3, Pos. 2) should be checked once a week for dirt. A visual check can be made without dismantling the filter.

If the filter is very dirty, it should be unscrewed with a slotted screwdriver and cleaned in an ultrasonic cleaner.

The sinter filters (REF 908-587-00) should be renewed every 1–2 years depending on how much the unit is used. This is the only way to ensure that the fine filter inside the unit does not get dirty.

Note: Never use the **AIRVAC 2002** without the sinter filters.

6.4 Changing the fuse

To change the fuse, first switch off the unit and disconnect it from the mains.

Remove the fuse holder on the multiple unit (illus. 3, Pos. 1) and insert a new fuse T 3, 15 A.

7. Troubleshooting

The following errors can cause faults.

Fault / Fault indicator:	Cause:	Remedy:
- Unit does not mix or build up a vacuum - The POWER indicator on the control panel does not light up	Not plugged in	1. Insert plug
	Unit not switched on	2. Switch the unit on at the ON/OFF switch
	Fuse blown	3. Change the fuse in the multiple unit (illus. 3, Pos. 1). Δ Refer to section 6.4!
- Vacuum pump builds up vacuum, but the mixer does not switch on	Mixing bowl is not properly connected to the mixing mechanism	4. Connect the mixing bowl and mixing mechanism more firmly at the flange of the bowl
	Dirty seals	5. Check seals, clean and change if necessary
- Red LED VACLOW indicator lights up - There are bubbles in the mixture after mixing	O-rings, seal surfaces, hoses and connections dirty, faulty or leaking	6. Check seals, clean and change if necessary
	Sinter filter (illus. 3, Pos. 2) in the flange of the bowl is blocked	7. Clean or change the sinter filter
	Fine filter inside the unit is blocked	8. Δ If there is no improvement after carrying out procedures 6 and 7, have the fine filter changed by customer services
	Blocked vacuum or sensor hole	9. Push a wire through the hole and clean

If faults cannot be remedied using this fault-finding chart, call your service technician or call Dentaurum customer services directly.


8. Spare parts

Mains cable	REF 907-027-10
Fuse T 3, 15 A	REF 907-037-00
Sinter filter	REF 908-587-00
Fine filter	REF 908-589-00
O-ring, large 100 x 5 for bowl flange	REF 908-405-00
O-ring, small 32 x 5 for bowl flange	REF 908-528-00
Set of seals for vacuum pump	REF 908-612-10
Pressure switch 300 mbar	REF 908-509-01
Pressure switch 720 mbar	REF 908-509-02

9. Accessories

Unit base	REF 095-075-00
Mixing bowl, large	REF 095-062-00
Mixing bowl, medium (supplied)	REF 095-063-00
Mixing bowl, small	REF 095-064-00

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