Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo

Paramil 3

Fräsgerät | Milling machine | Fraiseuse | Microfresadora

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Paramil 3
Fräsgerät

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Paramil 3
Milling machine

Contents

1 Declaration of Conformity	20
2 Safety instructions	21
2.1 Approved use	21
2.2 Symbols indicating risk	21
2.3 Guarantee and liability	21
2.4 Obligations of authorised user	22
2.5 Obligations of personnel	22
2.6 General information for safe use of the machine	22
2.7 High risks	22
2.8 Packing	23
2.9 Delivery	23
2.10 Package contents	23
3 Function and description of the machine	24
3.1 Function	24
4 Installation	24
4.1 Positioning the machine	24
4.2 Connecting the machine to the power supply	25
4.3 Connecting the units to the machine	25
4.4 Connecting the compressed air supply	25
5 Setting up and operating the machine	27
5.1 Inserting the cutter and drill into the chuck	27
5.2 Adjusting the height of the support bar	27
5.3 Unlocking the swivel arm	28
5.4 Using the milling head	28
5.5 Setting the motor speed and changing the rotational direction	28
5.6 Using the model table	28
5.6.1 Locking and moving the model table pneumatically	28–29
5.6.2 Clamping the model on the model table	29
5.6.3 Using the extension adapter for the model table	29
5.7 Using the compressed air cooling	29–30
5.8 Adjusting the height of the milling motor	30
5.9 Using the halogen lamp	30
5.10 Changing the chuck	30
5.11 Mount for attachment – parallel holder	30
6 Maintenance	31
6.1 Cleaning the milling machine	31
6.2 Changing the halogen bulb	31
6.3 Changing the fuse	31
7 Troubleshooting	32
8 Spare parts	32
9 Circuit diagram	33
10 Technical data	33

1 Declaration of Conformity

EC-Declaration of Conformity

DENTAURUM GmbH & Co. KG
Turnstr. 31
75228 Ispringen

hereby declares that the design and construction of the laboratory equipment described below, including the version marketed by us, comply with the basic regulations governing safety and health as stated in the EC Guidelines. This declaration will become invalid if the laboratory equipment is modified or altered in any way without our prior consent.

Unit:	Paramil 3 (REF 094-220-00)	
Description of unit:	Dental Milling machine	
Start with Unit No.:	107-01PAR 0010	
EC directives:	98/37/EC	Machinery
	73/23/EEC	Electrical equipment used within certain voltage limits
	89/336/EEC	Electromagnetic compatibility
Applied harmonized standards:	EN 292-1	
	EN 292-2	

Date and manufacturers signature: 01.01.2010
Signatory:

- i.V. Dipl. Ing. (FH) K. Merkle -
Production Manager Mechanic

2 Safety instructions

Caution: Read through the instructions for installing and operating the machine carefully before switching it on!

2.1 Approved use

The Paramil 3 milling machine has been developed solely for milling and drilling wax or metal in the dental laboratory. The machine should only be used for the purpose for which it is intended. Dentaureum GmbH & Co. KG is not liable for damages resulting from any other use of the machine. Approved use also includes adhering to the instructions for use as well as servicing and carrying out maintenance work at regular intervals.

When a Dentaureum product is finally no longer in use the disposal regulations of that particular country apply. Dentaureum or the dental trade is available to answer questions regarding the correct disposal of any specific product.

2.2 Symbols indicating risk

The following symbols are used in the instructions to indicate risks:

Warning: Indicates a potential danger to the life or health of the person. Disregarding this sign can critically affect a person's health and may even result in life-threatening injuries.

Caution: Indicates a potentially dangerous situation. Disregarding this sign can result in minor injury or damage to the equipment.

2.3 Guarantee and liability

Our general conditions of sale and delivery apply. Guarantee and liability claims for damages to persons or equipment are invalid if the damage is found to have resulted from any of the following:

- improper installation, operation, assembly or maintenance of the machine
- use of the machine other than for the purpose intended
- operation of the machine with faulty safety equipment or improperly installed or non-functioning safety and protective precautions
- disregarding the guidance in the instructions for use relating to transport, storage, assembly, operation and maintenance of the machine
- inadequate monitoring of parts liable to wear and tear
- unauthorised structural alterations to the machine
- repairs which have not been carried out properly

2.4 Obligations of authorised user

The authorised user is obliged to ensure that all those who operate the machine

- are familiar with safety regulations and accident prevention and have been instructed on how to operate the machine
- have signed that they have read and understood the safety precautions and instructions for use
- have been instructed in the current regulations for the prevention of accidents.

2.5 Obligations of personnel

Before commencing work, all those who operate the machine undertake

- to adhere to the basic regulations regarding safety at work
- to sign that they have read and understood the safety instructions and instructions for use.

2.6 General information for safe use of the machine

- Do not alter the machine in any way
- The machine should not be used if it is in any way faulty
- Keep the working area clean. Untidiness at the workplace increases the risk of accidents
- For reasons of personal safety, use only accessories and materials recommended in the instructions. Using materials other than those recommended in the instructions constitutes an accident risk for the operator and cannot therefore be regarded as approved use
- Always disconnect the machine from the mains before undertaking any maintenance work

2.7 High risks

- Do not touch the cutter or drill when it is rotating
- Instruments should only be changed when the machine is switched off
- The milling handpiece should only be used for milling work on the machine
- Only light pressure should be used when milling

Caution: Always wear protective glasses when milling or drilling.

2.8 Packing

The packing provides only limited protection against wet, heat and any mechanical force.

2.9 Delivery

- Note** This symbol shows that the marked product must not be disposed of as normal domestic waste. The law prohibits the disposal of electronic waste in communal collection areas. Further information can be obtained through Dentarurum or your dental distributors.

2.10 Package contents

- 1 Paramil 3 milling machine
- 1 Control unit for the micromotor
- 1 Pneumatic model table
- 1 Model plate for the model table
- 1 Foot control
- 1 Mains cable
- 1 Connecting cable to the micromotor
- 1 Compressed air hose \varnothing 6 mm (6 x 4 x 1 DIN 73378) with quick-release coupling KS-PK4
- 1 Chuck \varnothing 3.0 mm
- 1 Chuck \varnothing 2.35 mm
- 1 Chuck adaptor \varnothing 1.6 mm
- 1 Chuck key
- 5 Allen keys 5 mm, 3 mm, 2.5 mm, 2 mm, 1.5 mm
- 1 Instructions for use and EC Declaration of Conformity
- 1 bottle of oil for lubricating the moveable parts
- 1 bottle of milling oil
- 1 extension adapter for the model table

3 Function and description of the machine

3.1 Function

The Paramil 3 is a universal milling machine for use in dental technology. It can be used both for milling and drilling in wax and metal.

The built-in motor is exceptionally quiet and has a high torque. It can be operated clockwise or anti-clockwise and has a variable speed of 1500 min⁻¹ to 27000 min⁻¹.


Other features of the Paramil 3 are:

- variable motor speed control with overload cut-out
- well-positioned drive shaft allowing smooth movement without any free play
- double-hinged arm with a locking mechanism
- movable milling table with a pneumatic lock
- working area well lit by a swivel halogen lamp
- cooled milling with oil misting option
- easy operation
- compactness
- easy maintenance

4 Installation

4.1 Positioning the machine

The machine should be placed on a solid base, e.g. on a sturdy table, and set level. It should be installed in a place which is clean and as free from dust as possible.

Caution: Furnaces or other equipment, which emit heat, should not be placed next to the milling machine.

Illustration 1: Main view of the milling machine and control unit

4.2 Connecting the machine to the power supply

The electricity specifications on the model plate of the machine should be the same as the mains voltage.

- The machine should be connected to a 230 V socket
- The socket should be fitted with a 10 or 16 amp fuse

Caution: All electrical work on the socket or electricity supply and fuses should only be carried out by a qualified electrician.

4.3 Connecting the units to the machine

Position the control unit (illus. 1, 3) under the milling table (illus. 1, 1) and connect the socket on the control unit (illus. 3, 3.7) with the socket on the back of the milling machine (illus. 2, 2.2) using the cable supplied. The two 5-pin plugs should be inserted and secured with the swivel nut.

Plug the mains cable supplied into the back of the control unit (illus. 3, 3.2) and connect to the power supply.

Insert the 3-pin plug of the foot control cable into the socket at the back of the control unit (illus. 3, 3.8) and secure with the swivel nut.

The light connector (illus. 3, 3.9) is intended for external lighting.

4.4 Connecting the compressed air supply

A compressed air supply is required for operating the movable milling table with the pneumatic lock as well as for air-cooling.

The machine is connected to the existing compressed air supply in the laboratory or to an external compressor using the following method.

Caution: Check the pressure of the compressed air supply in the laboratory before connecting the machine to it. Ensure that the compressed air pressure is at least 3 bar but does not exceed 4 bar. If the pressure is too high, a pressure reduction valve should be fitted.

- The foot control also has two compressed air hoses to turn the air cooling for milling on or off. To install the hoses, hold the blue socket (illus. 2, 2.8) on the machine firmly and insert the blue compressed air hose of the foot control. The socket automatically clamps the hose tightly, so no additional retention is required. Then hold the black socket (illus. 2, 2.7) on the machine firmly and insert the blue compressed air hose with the black marking. The socket again automatically clamps the hose tightly.
- Insert the blue compressed air hose supplied (ø 6 mm) into the socket on the side of the milling machine (illus. 2, 2.1). Hold the socket on the side of the machine firmly and insert the compressed air hose. The socket automatically clamps the hose tightly, so no additional retention is required. To remove the compressed air hose, hold the blue ring on the socket firmly and pull out the hose.
- The compressed air hose (ø 6 mm) is then connected to the compressed air supply in the laboratory using the coupling. The coupling (KS-PK4) supplied is generally for universal use. It is possible that it may not be suitable for use in your laboratory. In this case contact your local supplier.


Illustration 2: Rear view of the milling machine



Illustration 3: Rear view of the control unit

5 Setting up and operating the machine

After the machine has been properly connected, it can be set up for milling.

5.1 Inserting the cutter and drill into the chuck

- Turn the chuck (illus. 4, 5.8) clockwise as far as it will go.
- Insert the cutter or drill into the chuck of the milling machine handpiece.
- Turn the chuck (illus. 4, 5.8) anti-clockwise to lock the cutter or drill in place. This automatically secures the instrument in the chuck.

⚠ Caution: Instruments with a cutting diameter of more than 3 mm should not be used.

⚠ Caution: To prevent accidents, remove the instrument after the work has been completed.

i Please note: Do not tighten the chuck without an instrument in it. Ensure that the instrument has the correct shank diameter.

i Please note: In addition to the \varnothing 2.35 mm chuck, a \varnothing 3 mm chuck and a \varnothing 1.60 mm chuck adaptor are available. How to change the chuck and use the chuck adaptor is described in section 5.10.


Illustration 4: Milling motor and chuck

5.2 Adjusting the height of the support bar

Loosen the adjustment screw on the support bar (illus. 1, 4.1) and set the bar to the required height. Then lock the support bar in place again by tightening the adjustment screw.

⚠ Caution: If the adjustment screw is not tightened, the bar should be supported to prevent it accidentally falling down.

i Please note: The friction screw on the left column (illus. 1, 4.2) should be tight enough to prevent the bar falling down, but loose enough to allow the support bar to be moved when the adjustment screw is loosened.

5.3 Unlocking the swivel arm

The milling arm lock (illus. 5, 5) can be released using the locking lever (illus. 5, 5.1). The milling arm then moves freely and can be positioned as required. Once it is in the correct position, the arm is locked again by pressing the lever.


Illustration 5: Detail view

5.4 Using the milling head

The machine is switched on at the ON/OFF switch (illus. 3, 3.1) on the back of the control unit. The motor is started by pressing the left foot control. If continuous operation is required please push the start-button (illus. 3, 3.5). The appropriate cutter is selected for the work to be carried out and inserted in the chuck before starting the motor.

5.5 Setting the motor speed and changing the rotational direction

The motor speed can be set between 1500 min^{-1} and 27000 min^{-1} using the motor speed controllers (illus. 3, 3.4). The motor speed is shown on the display screen (illus. 3, 3.3).

The rotational direction is set using the switch (illus. 3, 3.6) on the front of the control unit.

Caution: The rotational direction switch should only be used when the motor has stopped.

5.6 Using the model table

The model table (illus. 6, 6.1) is supplied with air through the compressed air connection (illus. 6, 6.2) enabling it to be moved easily over the base plate or locked in place.

5.6.1 Locking and moving the model table pneumatically

A secure connection to the compressed air supply (illus. 6, 6.1) is provided by a quick-release coupling which allows the hose to be easily inserted into the model table. To disconnect the supply, push the ring back and at the same time pull the hose out of the connection.

The pressure control (illus. 2, 2.3) on the model table is used for locking or moving the model table with the model. If the adjustment screw on the pressure control is turned anti-clockwise, the model table (illus. 6.1, 6.1.2) is locked onto the base plate by suction.

If the adjustment screw is turned clockwise, the air in the model table base is directed downwards forming an air cushion between the model table and the base plate. The model table can then be moved easily.

The basic pressure of the air supply to the model table is set using the red compressed air control (illus. 2, 2.3) on the column support.

Turning the control clockwise reduces or stops the flow of compressed air and turning it anti-clockwise increases the flow.

Illustration 6: Pneumatic model table


Illustration 6.1: Extension adapter

5.6.2 Clamping the model on the model table

The basic position of the clamp is adjusted to fit differently sized models by pressing down the jaws of the clamp and simultaneously pulling out the knurled screw (illus. 6, 6.4). The setting is finely adjusted by turning the knurled screw (illus. 6, 6.4).

The front clamp jaws (illus. 6, 6.3) can be moved eccentrically, if required, using an Allen key.

The model table can be swivelled by loosening the lever (illus. 6, 6.5). The model table is set at the required angle by tightening the lever and, if necessary, by using the additional locking mechanism.

5.6.3 Using the extension adapter for the model table

In order to raise the height of the clamping device, screw the extension adapter (illus. 6.1) between the lower part and the clamping device (illus. 6.1, 6.1.1).

5.7 Using the compressed air cooling

Instruments can be cooled using compressed air. There is also a small reservoir above the air nozzle, which can be filled with oil. The preset air jet is operated using the right foot control.

The force of the air jet cooling is regulated using the blue compressed air control (illus. 2, 2.4).

The air and oil nozzles (illus. 4, 5.10) are clamped to the handle of the micromotor (illus. 4, 5.7) and adjusted to the correct height.

Milling oil from the reservoir can be added as required, providing a film of oil which makes it easier to smooth the metal when milling.

Oil misting (illus. 8) is turned on or off using the sliding mechanism (illus. 4, 5.11).


Illustration 8: Diagram of the air coolant nozzle

5.8 Adjusting the height of the milling motor

The height adjustment mechanism of the milling motor, which has spring suspension, can be locked using the adjustment screw (illus. 5, 5.2).

The guide bar (illus. 5, 5.3) is inserted into the socket next to the micromotor and has several functions.

- It allows the height to be adjusted. The height is set using the lateral adjustment screw. The eccentric head is supported on the swivel arm.
- It provides an aid for milling channels. The adjustable height setting (illus. 5, 5.5) of the guide bar enables channels to be milled precisely when the milling arm is locked.
- It provides an aid for drilling and milling interlocks. The adjustable socket (illus. 5, 5.4) of the guide bar (illus. 5, 5.3) enables interlocks to be precisely drilled or milled when the milling arm is locked.

! Please note: The milling motor must be moved to the left side of the swivel arm before the guide bar can be inserted into the socket (illus. 5, 5.4).

5.9 Using the halogen lamp

The halogen lamp is attached to a flexible arm which allows the angle of the lighting to be adjusted. The lamp is switched on and off by turning the ring on the head of the lamp.

5.10 Changing the chuck

- Turn the chuck (illus. 4, 5.8) clockwise as far as it will go.
- Remove the drill or cutter.
- Put the triangular chuck key supplied onto the chuck.
- A screw in the shank of the chuck locks the chuck into the drive shaft. The 1.5 mm Allen key is inserted into the chuck to loosen the screw.
- Hold the chuck firmly with the chuck key and turn the 1.5 mm Allen key anti-clockwise until the chuck can be loosened manually.
- Turn the chuck clockwise to remove it from the shaft.
- The new chuck with an instrument already inserted is now screwed manually into the drive shaft as far as it will go.
- Unscrew the chuck by a half turn in a clockwise direction. This allows the instruments to be inserted more easily by creating a little more space between the chuck and the instrument. Then remove the instrument.
- Place the triangular chuck key supplied onto the chuck.
- The screw in the shank of the chuck is now tightened again using the 1.5 mm Allen key.

5.11 Mount for attachment – parallel holder

An attachment holder can be fixed in the mount at the milling arm for the parallel alignment of prefabricated attachments (illus. 5, 5.12).

6 Maintenance

Warning: Always disconnect the machine from the mains during any servicing and maintenance work. Never work alone! There should always be two people with electrical experience present when servicing and repairs are being carried out.

6.1 Cleaning the milling machine

Wipe the milling machine occasionally with a damp cloth. Do not use abrasive cleaning agents or tools.

After each milling, the machine should always be cleaned of any milling dust and splinters with a brush and checked to ensure that it moves freely. It is important to ensure that no metal splinters get into the guiding mechanism of the columns.

Caution: Never use compressed air, water or solvents to clean the milling machine. Do not steam clean the machine.

6.2 Changing the halogen bulb

Turn the ring on the lamp head right towards the arrow .

The cover can now be taken off. Remove the halogen bulb from the socket using tweezers or pliers.


To insert the light bulb and assemble the lamp, simply reverse the procedure.

Caution: Do not handle the glass of the halogen bulb, as fingerprints can destroy the light bulb when it is switched on.

6.3 Changing the fuse

In order to replace the fuses all connections must be disconnected from the control unit. Pull forward the control unit under the milling table, then unscrew the 4 upper fixing screws of the cover and remove the latter. The fuse holders can now be removed upwards and the fuses can be replaced.

Fuse for the halogen lamp
Glass fuse 5X20 – T 630mA

Fuse for the motor
Glass fuse 5X20 – F4A

Illustration 7: interior view of the control unit

7 Troubleshooting

The following errors can cause faults.

Fault	Cause	Remedy
The milling machine does not start when the foot control is operated	Not connected to the mains	Insert the plug on the control unit (illus. 3, 3.2)
	Control unit is not switched on	Switch the control unit on at the switched on ON/OFF switch (illus. 3, 3.1)
	Foot control is not connected	Insert the connecting cable of the connected foot control into the socket on the control unit (illus. 3, 3.8) and screw on tightly
	Micromotor is not connected	Insert the connecting cable of the connected micromotor into the socket on the spindle column and control unit (illus. 3, 3.7) and screw on tightly
	Faulty fuse	Change the fuse in the control unit (illus. 7).  Refer to section 6
	Faulty milling head	Inform Dentaaurum customer services
Halogen light does not function	Faulty halogen light bulb	Change halogen light bulb  Refer to section 6
	Faulty fuse	Change the fuse in the control unit (illus. 7).  Refer to section 6
Vertical movement not smooth	Dirt in the columns	Clean the columns with a brush

If faults cannot be remedied using this fault finding chart, call your service technician or call the Dentaaurum customer services directly.

Tel. +49 72 31/803-211

E-Mail: info@dentaaurum.de

8 Spare parts

REF

Mains cable	907-027-10
Fuse T 630 mA	907-032-00
Fuse F 4 A	
Halogen light 12V– 5W	907-999-00
Chuck ø 2,35 mm	
Chuck ø 3,00 mm	907-647-00
Milling oil	094-222-20

Instructions for use and EC Declaration of Conformity

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