rematitan® Stege

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EN

rematitan® bars

Bars for the welding process of titanium materials according to DIN EN ISO 5832-2

Composition (% by mass)

Ti grade 4 - 99.0 % Ti

Additional elements < 1 ‰: Fe, O

Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG I Turnstraße 31 I 75228 Ispringen I Germany

2. General product description

The rematitan® bars are used in conjunction with laser welding technology. Prefabricated implant sleeves in titanium as well as abutments in titanium are used for the abutment teeth.

Due to the same cross-section as in the case of the known precious metal Dolder® bars, the known auxiliary parts such as placeholders, paralleling bars and activation keys can be used. The sleeves can also be interchanged.

3. Area of application

rematitan® bars are suitable for edentulous jaws. They are used to fabricate implant-supported bar bases (patrices and matrices) produced by laser welding and serve as anchorage for full plastic dental prostheses.

4. Processing

During welding, the point of contact between the objects to be welded must always be clean. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Welding may only take place under a shielding gas atmosphere (preferably argon) and is an unconditional prerequisite for good welding results.

5. Contraindications and adverse reactions

Signs of intolerance to titanium and titanium alloys are extremely rare if the manufacturer’s Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur. Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Warnings and safety data sheets must be observed.

6. Warnings and precautions to be taken

Metal chips and metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. Due to the oxygen affinity of titanium, there is a risk of combustion and deflagration, particularly in the hot state. Any type of ignition source shall be kept at a distance. Metal fire powder extinguishers should always be positioned in close proximity to extinguish a fire if necessary. They should be in working order. The workpiece must be cooled during mechanical machining in the presence of oxygen. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The material should not be ground or polished in the patient’s mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

7. Reuse

The material may be changed in its chemical composition and material properties if reused. In principle, we recommend using new material to manufacture high-quality prosthetic restorations so that clear bath traceability is also possible.

8. Storage conditions

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of titanium and titanium alloys.

9. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

10. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

11. Additional tips

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

12. Explanation of symbols used on the label

 Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germania