Smile on Line

Orthodontie

N° 18 - Novembre 2008

SOMMAIRE

Edito

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Page 2

rematitan® sl

Page 3

Publi-communicé discovery® sl

Pages 4 à 7

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Pages 8 à 11

Mini-vis «TAD insertion :
Angles for success»

Pages 12 à 13

Nouveautés

Pages 14 à 15

Innovations

Page 16

Formations 2008 / 2009

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Introduction

By **Dr. Sebastian BAUMGAERTEL**

Sebastian Baumgaertel, DMD, MSD, FRCD(C) received his Orthodontic education at Case Western Reserve University where he now holds the position of assistant clinical professor and is co-director of the subspecialty clinic for skeletal anchorage. In addition, Dr. Baumgaertel maintains an active private practice in the Cleveland area. He is a Diplomate of the American Board of Orthodontics, a Fellow of the Royal College of Dentists of Canada, and a certified orthodontic specialist in Germany.

Dear Colleagues,

Without doubt, we have the privilege to practice in very exciting times. Our specialty has made incredible progress over the past years, largely due to the introduction of new technologies and treatment techniques.

One of the most interesting new concepts in orthodontics is the use of absolute anchorage through implementation of temporary anchorage devices (TADs) in our treatment mechanics. These mini-implants or mini-screws allow us to move teeth in a never before possible fashion. Their use results in more efficient mechanics and thus reduced treatment time, a reduction in required patient compliance, and often times a reduction in the level of invasiveness of the treatment.

I have been using mini-implants since their infancy stage here in the United States. While initially I used up to four different systems simultaneously, over time, Dentaureum's tomas® System proved to work best in my hands which prompted me to make it the exclusive system in my practice. The tomas® Pin has evolved quite a bit since the early days. Today clinicians can choose between the original non-drill free thread design and a drill-free thread design (self-drilling) which does not require any pre-drilling (my favorite) and three different lengths. What I particularly value is the versatility of the implant head which allows me to use direct and indirect anchorage solutions as needed in the individual patient. To me, a great proof for how well the Pin's head is designed is the fact that it has changed very little over time. While other systems constantly redesign the heads to compensate for shortcomings, the universal tomas® Pin head is time tested and proved successfully.

Placing a mini-implant is not difficult and the learning curve is steep. Following my recommended steps and respecting the information from Dr. Petrey's study will improve your clinical success rates. Nevertheless, even if one masters all the theory, implementation into clinical practice can be a challenge. That is why Dentaureum developed the tomas® Auxiliary Kit which, after placement of the TAD, allows for quick and easy immediate loading using simple mechanics and prefabricated auxiliaries. Biomechanics are the cornerstone of successful clinical TAD use. Therefore I did not have to think twice when Dentaureum and Orthodontic Products Magazine approached me to compile the seven most commonly used TAD indications for you and explain my biomechanical designs that lead to successful and predictable attainment of the treatment goals. You can keep this little guide book as a reference when installing your mini-implant biomechanics.

If you haven't already started using TADs, today would be a great day to get your fingers wet.

Yours truly,

Sebastian Baumgaertel, D.M.D., M.S.D., FRCD(C)

TAD insertion : angles for success

By Joseph S. PETREY

Temporary Anchorage Devices are becoming the preferred method of skeletal anchorage and a significant component to the modern orthodontist's armamentarium. While not useful in every case, TADs offer the clinical orthodontist the opportunity to use many mechanics not possible using tooth-borne or extra-oral anchorage systems. As with most orthodontics, the key to proper utilization of temporary anchorage devices is in diagnosis. Once it is determined that TADs would benefit treatment, the next critical component is the utilization of the TAD in a manner which will lead to the anchor's success throughout treatment.

Numerous components of TAD utilization contribute to the anchor's success or failure. Components which have been referenced include : local inflammatory factors, force levels utilized, loading techniques, the use of a pilot hole, implant design characteristics, and the method in which the TADs are inserted.^{1,2} Many of the components of TAD success are outside of the orthodontist's control. One factor for success which is well within the orthodontist's control is the manner in which the anchor is inserted and loaded.

As outlined, TADs should be placed in areas of attached gingiva with local periodontal factors controlled, but beyond these basic rules, the angle the implant is placed and loaded is critical for success. Previous reports have suggested TADs be placed at a vertical angle to the cortical plate to reduce the risk of root contact.³ This orientation, however, does not provide for the best retention of the TAD or its successful use. A more appropriate approach is to select the proper TAD location and then utilize traditional mechanics to position roots prior to TAD placement, allowing for ample room for the implant.

Once an appropriate position has been selected for the placement of the TAD, it is important to ensure the TAD is properly placed in this location. In research conducted at the University of Kentucky and independent studies from China Medical University, it has been shown there is a dramatic drop in the force levels TADs will withstand when placed in an orientation other than perpendicular to the cortical plate.^{4,5} Contradicting reports have suggested placement of TADs in a Tweed or tent-peg fashion in the cortical plate to the force is the most retentive.⁶ Unfortunately, these studies have defined failure as complete removal of the TAD from the bone, and in many cases fracture of the surrounding bone. The force levels of these experiments are well above those reached in traditional orthodontic mechanics, and the definition of failure is well beyond the usefulness of a TAD. For clinical practice, failure of a TAD is when anchorage is no longer absolute.

Defining TAD failure as the initial loss of anchorage, clinical tests have demonstrated that implants placed at extreme angles 45° and 135° to the cortical plate fail at a markedly reduced force levels than those placed 90° to the cortical plate. Implants placed in a tent-peg fashion have been shown to fail at 33% less force than those placed 90 degrees to the cortical plate.⁶ These failure force levels were found to be well within traditional orthodontic loads.

Dr. Petrey completed his orthodontic training at the University of Kentucky Orthodontic Graduate Program. He maintains an active private practice in Somerset, KY.

Not only is it important to consider the angle of the implant with respect to the cortical plate when placing implants, but also the direction in which orthodontic force is placed on the TAD. TADs should be placed so they may be loaded in the same plane as the force acted upon them, with the TAD perpendicular to the direction of force placed upon them. Implants placed in an angular fashion and loaded tangentially, have been shown to fail at 45% less force than implants placed 90° to the cortical plate and perpendicular to the force acted upon them.⁶ Although, previous studies have recommended inserting implants at an extreme angle to the cortical plate to reduce the risk of root contact on insertion³, this recommendation significantly reduces the amount of force the implant will withstand before failure. A more appropriate technique is to complete pre-implant alignment preparing the implant site so that the TAD may be placed perpendicular to the cortical plate without risk of root contact.

The following images demonstrate the insertion of TADs for the purpose of protracting maxillary molars in a case without maxillary second premolars. The implants were placed perpendicular to the cortical plate, and the load placed upon the implants was designed to be in the same plane as the TAD, and perpendicular to the direction of the force.

TAD being placed perpendicular to the cortical plate

TAD placed perpendicular to the cortical plate


TAD placed perpendicular to the cortical plate and loaded in preparation of molar protraction

Note not only is the implant perpendicular to the cortical plate, but it is perpendicular to the force acted upon it.

TAD success is dependent on many factors for optimal retention. Whenever possible, clinicians should develop absolute anchorage systems with the greatest opportunity for success. These attempts include: proper placement protocol following clinical guidelines, and when possible implants should be placed 90° to the cortical plate, in the same plane as the force placed upon them, and loaded with the TAD perpendicular to the direction of the force acting upon it.

1. properly placed 2. occlusally placed 3. gingivally placed

1. A properly inserted TAD with the implant perpendicular to the cortical plate.
2. TAD improperly placed with an extreme occlusal angle to the cortical plate. TADs placed in this manner will fail at significantly lower force levels and placement will be more difficult; many times requiring a pilot hole or notching of the bone to secure a purchase point for insertion.
3. TAD improperly placed with an extreme gingival angle to the cortical plate. TADs placed in this manner will fail at significantly lower force levels and placement will be more difficult; many times requiring a pilot hole or notching of the bone to secure a purchase point for insertion.

1. Park HS, Jeong SH, Kwon OW. Factors affecting the clinical success of screw implants used as orthodontic anchorage. Am J Orthod Dentofacial Orthop. 2006 Jul;130(1):18-25.
 2. Miyawaki S, Koyama I, Inoue M, Mishima K, Sugahara T, Takano-Yamamoto T. Factors associated with the stability of titanium screws placed in the posterior region for orthodontic anchorage. Am J Orthod Dentofacial Orthop. 2003 Oct;124(4):373-8.
 3. Kyung HM, Park HS, Bae SM, Sung JH, Kim IB. Development of orthodontic micro-implants for intraoral anchorage. J Clin Orthod. 2003 Jun;37(6):321-8.
 4. Zang Y, Zang D, Cuijuan F, Peng P, Hu H, Kwakami T, Takagi T, Nagai N. A three dimensional finite element analysis for biomechanical characteristics of orthodontic anchorage micro-implant. J Hard Tissue Biology. 2006; 15(2): 69-72.
 5. Petrey JS, Kluemper GT, Beeman CS, Cunningham LL, Saunders MM. Orthodontic Mini-Implant Insertion Variables: Effects on Retention. 2008 Article in press.
 6. Pickard MB, Dechow P, Rossouw PE, Buschang PE. Effects of Mini-Screw Orientation on Implant Stability and Resistance to Failure. Am J Orthod Dentofacial Orthop.

Overview of auxiliaries

tomas®-coil springs

Light, Medium, or Heavy tension
Custom-designed NiTi spring with
a larger eyelet to fit the tomas®-pin.
10 pieces, [REF 302-012-\(00/10/20\)](#)


tomas®-Nikodem springs

7mm, 11mm, or 14mm lengths
Superelastic NiTi spring designed by
Dr. Nikodem. Custom-designed to
easily fit over the tomas®-pin head.
10 pieces, [REF 302-016-\(07/11/14\)](#)

tomas®-cross tube

Cross tube for connecting the **tomas®**
pin to the archwire. Ideal for indirect ancho-
rage. Available for the 18 and 22 technique.
10 pieces each,
[REF 302-014-18](#) / [REF 302-014-22](#)

tomas®-Uprighting spring

Provides 3-D tooth control for
simultaneous intrusion or
extrusion while uprighting
molars. 10 pieces,
[REF 302-009-00](#)

tomas®-power arms

square or round available
Crimpable square power arm that
easily connects the **tomas®**-pin
to the existing arch. Crimpable round
power arm for custom bending hooks
at desired vector level. 10 pieces,
[REF 302-015-00](#) / [302-019-00](#)


tomas®-hook

Additional option for securing
elastic components to the head
of the **tomas®**-pin.
10 pieces, [REF 302-009-10](#)

tomas®-crimp hook 4L / 4R

Crimpable tube with 4 steps for atta-
ching springs and elastics at various
vector points. Compatible with 18 or 22
slot. 10 pieces each,
left side: [REF 302-020-00](#)
right side: [REF 302-020-10](#)

tomas®-compression spring

Superelastic compression spring
can be custom-cut to desired length.
Maximum elasticity with constant,
medium recovery. 10 pieces,
[REF 302-022-00](#)

tomas®-power chain

Plastic chain with high elasticity
and constant tension force.
1 piece, [REF 302-018-00](#)

tomas®-T-wire

21 x 25 SS T-wire used for
connecting and locking in
segments of teeth to a palatal
TAD 1 piece,
[REF 302-024-00](#)

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