



Multi-indicative vestibular screen for Face former therapy.

by Berndsen





Neurophysiological – cause-oriented.

Wide range of indications.

- Swallowing disorders (dysphagia)
- Malfunctions of the tongue
- Impaired mastication
- Disorders of speech, speaking and fluency
- Avoiding thumb or pacifier sucking, lip biting behavior etc.
- Muscle weakness or paralysis of the mouth and/or face
- Before or after face, jaw and tongue surgery
- Malpositioned teeth, jaw anomalies
- Teeth grinding (bruxism)
- CMD (craniomandibular dysfunction)
- Dental restorations and implants
- Sleep disorder, snoring, sleep apnea
- Hearing loss, tinnitus
- Dizziness, Menière's disease
- Problems compensating pressure in the middle ear
- Switching from breathing through the mouth to breathing through the nose
- Migraine
- Activation of cranial nerve functions
- Problems with head balance and body statics
- Impaired salivation
- Influence on metabolism, diabetes
- Head, face, neck and back pain
- Cosmetic: slack skin and muscles of neck, face and chin
- Asthma
- Reflux

Face former training.

Regular Face former training has an effective impact on cranial nerve function and brain performance, neurophysiological patterning, coordination, tension of tissues and muscles, functions in the mouth, throat and facial regions, body posture, body movement and sensory abilities. Causal therapies, such as the Face former therapy, are not successful from one day to the next. The objective is rather a more gradual approach to the desired goal and then stabilization of achievement, as well as gaining independence from harmful medication or aids.

Excellent functionality.

Face former anatomical

Face former anatomical is delivered in a standard size and usually does not need adjusting.

It is the standard therapy device.

by Berndsen



Direct effect.

First signs of success can already be seen after three weeks. A positive development can be expected after 12 – 16 weeks. Compared with other techniques for myofunctional therapy, lip strength can be improved to a maximum. Visible success after a short time encourages the patient to continue therapy.

High biocompatibility.

The material contains no latex, bisphenol-A or harmful plasticizers. The careful manufacturing of the Face former is carried out in a complex procedure. It is made from the silicone rubber Silopren®. Silopren® is one of the highest quality materials produced today. It is biocompatible and was especially developed for use in medical products.

Excellent hygiene.

Face former is cleaned by rinsing it under hot water. A toothbrush can be used if necessary.

Dentaurum

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE

	www.dentaurum.com		
			

Date of information: 11/19
Subject to modifications

