Estimado cliente

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 · 75228 Ispringen · Germany

2. General product description

Herbst IV is a further development of the standard Herbst I hinge. It has a ball and socket hinge that allows the patient greater freedom of move- ment, increasing comfort and tolerance during treatment. Retention using lock washers (C-clips) facilitates fitting by the orthodontist and replaces the relatively awkward insertion of screws.

3. Intended purpose

Correction of distocclusions.

4. Intended users

All products are intended for dental professionals only.

5. Indications

Distocclusion/class II (unilateral/bilateral), molar distalization, gap closure (aplasia), temporomandibular dysfunction, sleep apnea therapy, dentoalveolar compensation of occlusion (elastics substitution).

6. Contraindications

The Herbst appliance made of stainless steel must not be used for patients with a nickel allergy. Warnings are found on the respective product labels. Also see 9. Composition.

7. Patient target group

All patients for whom a fixed appliance can be placed as part of ortho- dontic treatment.

8. Fabrication of the Herbst IV appliance

Cast the impressions, which include the waxed out bands fitted on the first premolars and molars, with class IV dental stone, e.g. Rapidur® (REF 164-450-00). Mount the models in the correct position in the orthodontic FKO-split-post fixator (REF 072-004-00) with plaster using the construction bite recorded by the dentist. Bend remaloy® wire, half-round, e.g. ☐ 1.50 mm x 0.75 mm (REF 528-155-00), or round, e.g. ø 1.00 mm (REF 528-100-00) to form a lingual bow and then weld the wire to the bands. A palatal arch can be placed in the upper jaw to provide stabilization, e.g. remanium® Goshgarian palatal bar or one made from round ø 1.0 mm remaloy® wire (REF 528-100-00). Weld the base of the retention hinges to the bands, e.g. using the laser welding unit desktop Compact (REF 090-578-50).

Alternatively, splints can be cast instead of using bands. The bases of the retention hinges should not be cast on, but welded afterwards, for accuracy and biocompatibility reasons. Ensure that the screws are at the correct angle to the axis between the upper and lower jaw. It is important that the retention hinges operate without friction to ensure smooth functioning of the telescope.

Fitting the telescope tubes and rods:

First, shorten the telescope tubes according to the construction bite and then the telescope rods which are approx. 2 mm longer. The correct length should be checked again intra-orally and adjusted if necessary.

Note: "Double bands" should be used to minimize the risk of the bands being damaged by shear forces. This involves welding a second band over the first band, e. g. using the spot welding unit Assistent 3000 (REF 078-000-00) or desktop Compact laser (REF 090-578-50).

Generally the second band is 2 to 3 sizes larger than the first band fitted. Ready-welded "double bands" can be ordered from Dentaurum. The order form is available on request.

Fitting intra-orally:

The bands are fitted intra-orally using standard band cement. After the cement has set, the prepared retention hinges are fitted on the bases and secured with a lock washer, C-clip. When the hinges are fitted correctly on the base, it should be possible to insert the C-clip into the two slots of the ball unit so that the ends are clipped in at both sides.

The first recall appointment should be one week following fitting of the appliance!

If necessary, the Herbst IV appliance can be activated using spacer rings, (please refer to our product range).

9. Composition

The Herbst IV appliance is made of the following materials: DIN 1.4301, DIN 1.4305 and DIN 1.4310. The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSPC is available at https://ec.europa.eu/ tools/eudamed and www.dentaurum.com.

10. Safety information

Ensure that the C-clip is fitted parallel to the slot of the ball hinge! If the C-clip is fitted at an angle, it can be bent and become loose in the mouth. Always check that the C-clips are properly seated after fitting the Herbst IV appliance and at each recall appointment. If the C-clips are repeatedly fitted and removed, they lose their spring tension and should be replaced!

The product should not be used if there is a known allergic reaction to one or more of the material components. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accord- ingly.

For reasons of biocompatibility, welding (e.g. with a laser welding unit) should be chosen over soldering when preparing the appliance.

13. Notes on disposal

When disposing of the product, please adhere to the national regulati- ons that are valid for your region and - if available - the details outlined in the safety data sheets.

14. Product range

For the complete range of Herbst products, please see the orthodontic catalog or visit www.dentaurum.com.

15. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental pro- fessional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

16. Explanation of symbols used on the label

 **Please refer to the label. The Instructions for use and additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).**

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