

Ortho-Cast™ mini buccal tubes



CE 0483

Perfect in combination with equilibrium® mini brackets

- Biocompatible 1-piece construction
- 3D-base for excellent fit
- Extremely easy to handle – due to large chamfered archwire entry slots and “large footprint” base

Teeth	Torque	Distal Offset	Width	System	REF ¹⁾
17-16	- 14°	0°	2,5 mm	Roth 18 /	724-250-51
26-27	- 14°	0°	2,5 mm	MBT* 18	724-251-51
47-46	- 30°	0°	2,5 mm	Roth 18	724-252-51
36-37	- 30°	0°	2,5 mm		724-253-51
17-16/36-37	0°	0°	2,5 mm	Standard	724-262-51
26-27/47-46	0°	0°	2,5 mm	Edgewise 18	724-263-51

Teeth	Torque	Distal Offset	Width	System	REF ¹⁾
17-16	- 14°	0°	2,5 mm	Roth 22 /	724-254-51
26-27	- 14°	0°	2,5 mm	MBT* 22	724-255-51
47	- 10°	0°	2,5 mm	MBT* 22	724-260-51
37	- 10°	0°	2,5 mm		724-261-51
46	- 20°	0°	2,5 mm		724-258-51
36	- 20°	0°	2,5 mm		724-259-51
47-46	- 30°	0°	2,5 mm	Roth 22	724-256-51
36-37	- 30°	0°	2,5 mm		724-257-51
17-16/36-37	0°	0°	2,5 mm	Standard	724-264-51
26-27/47-46	0°	0°	2,5 mm	Edgewise 22	724-265-51

¹⁾ each 10 pieces

Our orthodontic technicians are available to answer any questions you may have regarding our products:

Fixed Appliance Technique Hotline +49 72 31/803-550
Removable Appliance Technique Hotline +49 72 31/803-555

Order Department +49 72 31/803-210

▶ Additional information regarding Dentaureum products can be found on the internet.

www.dentaureum.de


▶ Also note our special course offers:
 Further information can be obtained by contacting our hotline: +49 72 31/803-470

Date of information: 07/07


Subject to modifications

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

equilibrium®



The modern
1-piece bracket line
aesthetic – biocompatible – economical



989-525-20 Printed by Dentaureum Germany 07/07/CR1


D
DENTAUREUM

D
DENTAUREUM

equilibrium®

bracket line

Still easier bracket identification by FDI-marking on the base



equilibrium® mini

The smallest bracket in the world

- Aesthetics as yet unknown with metal brackets along with maximum functionality
- Biocompatible 1-piece construction
- Highest precision

CE 0483



equilibrium® mini

	1 case		10 cases		50 cases	
	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5
Roth 18	718-401-01	718-441-01	718-401-10	718-441-10	718-401-50	718-441-50
Roth 22	722-401-01	722-441-01	722-401-10	722-441-10	722-401-50	722-441-50
MBT* 18						
MBT* 22	722-402-01	722-442-01	722-402-10	722-442-10	722-402-50	722-442-50

equilibrium® 2

The proven standard bracket

- Biocompatible 1-piece construction
- Excellent functionality by innovative design
- Perfect fit – easy handling
- Also available for the systems Ricketts® Universal 18, Standard Edgewise 18/22 and Tweed 22

CE 0483

equilibrium® 2

	1 case		10 cases		50 cases	
	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5
718-301-01	718-341-01	718-301-10	718-341-10	718-301-50	718-341-50	
722-301-01	722-341-01	722-301-10	722-341-10	722-301-50	722-341-50	
718-302-01	718-342-01	718-302-10	718-342-10	718-302-50	718-342-50	
722-302-01	722-342-01	722-302-10	722-342-10	722-302-50	722-342-50	

equilibrium® ti

The perfectly biocompatible titanium bracket

- The only 1-piece construction made of pure titanium
- Highest precision
- Unbeatable price performance ratio

CE 0483

equilibrium® ti

	1 case		10 cases		50 cases	
	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5	Hooks on 3	Hooks on 3-5
718-501-01		718-501-10		718-501-50		
722-501-01		722-501-10		722-501-50		
718-502-01		718-502-10		718-502-50		
722-502-01		722-502-10		722-502-50		

