

	Description		REF	Quantity
	SUS²-Starter-Set 2 x Telescopic tube (with Allen screw); 2 x Telescopic rod, assembled; 2 x Ball retainer; 2 x Spacer spring, length 10 mm; 2 x Compressive spring, stainless steel, length 15 mm; 1 x Activation key (Allen); 1 x Instructions for use	1 Set	607-140-00	<input type="text"/>
	SUS²-Set 2 x Telescopic tube (with Allen screw); 2 x Telescopic rod, assembled; 2 x Ball retainer; 2 x Spacer spring, length 10 mm; 2 x Compressive spring, stainless steel, length 15 mm; 1 x Instructions for use	1 Set	607-145-00	<input type="text"/>
	SUS²-Telescopic rod, assembled Comprising: telescopic rod, retention screw, arch adapter	2 pieces	607-148-00	<input type="text"/>
	SUS²-Telescopic rod	2 pieces	607-149-00	<input type="text"/>
	SUS²-Arch adapter/ screw	2 pieces	607-143-00	<input type="text"/>
	SUS²-Telescopic tube (with Allen screw)	1 piece	607-147-00	<input type="text"/>
	SUS²-Ball anchor	10 pieces	607-136-00	<input type="text"/>
	SUS²-Compression spring (stainless steel) Length 15 mm	6 pieces	607-141-00	<input type="text"/>
	SUS²-Compression spring (nickel-titanium) Length 15 mm	6 pieces	607-150-00	<input type="text"/>
	SUS²-Spacer spring Length 10 mm	6 pieces	607-142-00	<input type="text"/>
	SUS²-Activation key	1 piece	607-144-00	<input type="text"/>

Customer:

Customer number	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>
Name	<input type="text"/>
Street	No. <input type="text"/>
Postcode	Town <input type="text"/>

Date / Signature _____

Dentaurum's terms of delivery and payment and other conditions apply when placing an order.