


Herbst TS



Telescoping Herbst Bite Jumping Hinges

In 3 different lengths

Herbst TS

With F-bases and hexagon socket screws, complete

- Contents:**
- 1 Telescopic attachment for right and left each
 - 4 F-bases
 - 4 Hexagon socket screws
 - 1 Instruction for use

F-Base

Hexagon Socket Screws

Key for Hexagon Socket Screw

For Spacer Rings, slotted, Herbst TS / SUS²

Clamping Plier for Herbst TS / SUS²

For clamping on slotted spacer rings for Herbst TS and SUS².

Premium-Line · 10 Year Guarantee

Length: 20 mm	607-126-20	1 Set	CE 0483
Length: 22 mm	607-126-22	1 Set	CE 0483
Length: 24 mm	607-126-24	1 Set	CE 0483

607-127-30	10 pieces	CE 0483
------------	-----------	---------

607-122-30	10 pieces	CE 0483
------------	-----------	---------

607-144-00	1 piece	Sterilizable	CE
607-123-10	1 piece	For laboratory use only!	—

Length: 1 mm	607-103-11	10 pieces	CE 0483
Length: 2 mm	607-104-11	10 pieces	CE 0483
Length: 3 mm	607-105-11	10 pieces	CE 0483
Length: 4 mm	607-107-11	10 pieces	CE 0483

003-710-00	1 piece	Sterilizable	CE
------------	---------	--------------	----