

Adhesives.

For perfect bonding results.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Our adhesives for your dental practice.

The range of adhesives from Dentaureum gives users a wide selection for all bonding tasks in the orthodontic practice.

The light-curing adhesives CONTEC Ic and CONTEC Icr have been especially designed for bonding fixed orthodontic appliances – CONTEC Ic is for bonding metal, ceramic and plastic brackets and buccal tubes, CONTEC Icr is for bonding orthodontic retainers using adhesive. Both adhesives have a very high bond strength. An alternative to the light-curing adhesive CONTEC Ic is the chemically-curing bracket adhesive CONTEC go. CONTEC Ic contains no BisGMA (bisphenol A glycerolate dimethacrylate), TEGDMA (triethylene glycol dimethacrylate) or HEMA (hydroxyethyl methacrylate). CONTEC Icr and CONTEC go are free of TEGDMA and HEMA. The adhesives are therefore also suitable for patients with an allergy to methacrylates.

CONTEC

Best compatibility, no TEGDMA or HEMA.



Adhesives – areas of application.

	CONTEC lc	CONTEC lcr	CONTEC go
discovery® sl	+++	+	+++
discovery® sl 2.0	+++	+	+++
discovery®	+++	+	+++
discovery® smart	+++	+	+++
discovery® pearl	+++	+	+++
discovery® delight	+++	++	+
topic	+++	+	+++
equilibrium® 2 / ti	+++	+	+++
equilibrium® mini	+++	+	+++
Retainer	++	+++	—
Buccal tubes	++	+	++

- +++ highly suitable
- ++ very suitable
- + suitable with restrictions
- not suitable

Adhesives and accessories



CONTEC Ic set

Light-curing.

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

For approx. 800 brackets (40 treatment cases).

A plastic primer must be used to bond plastic brackets.

Content:

- 3 x CONTEC Ic bracket adhesive, each 3 g (REF 163-111-00)
- 1 x CONTEC Ic / lcr primer, 5 ml (REF 163-113-10)
- 2 x CONTEC etching gel, each 2.5 ml (REF 163-115-00)
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

When used with metal and ceramic brackets, CONTEC Ic is free of BisGMA, TEGDMA and HEMA.

CONTEC Ic, used with plastic brackets, is free of TEGDMA and HEMA.

REF 163-110-00	1 set
-----------------------	-------



CONTEC Ic bracket adhesive

Light-curing.

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

A plastic primer must be used to bond plastic brackets.

When used with metal and ceramic brackets, CONTEC Ic is free of BisGMA, TEGDMA and HEMA.

CONTEC Ic, used with plastic brackets, is free of TEGDMA and HEMA.

REF 163-111-00	2 pieces (3 g each)
-----------------------	---------------------



CONTEC lcr set

Light-curing.

For bonding of retainers.

For approx. 80 retainers, braided or approx. 250 retainer 3-3.

Content:

- 2 x CONTEC lcr retainer adhesive, each 3.4 g (REF 163-131-00)
- 1 x CONTEC Ic / lcr primer, 5 ml (REF 163-113-10)
- 2 x CONTEC etching gel, each 2,5 ml (REF 163-115-00)
- 6 x CONTEC lcr disposable tube
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

CONTEC lcr is free of TEGDMA and HEMA.

REF 163-130-00	1 set
-----------------------	-------

CONTEC Icr retainer adhesive

Light-curing.

For bonding retainers.

Without disposable tube.

REF 163-131-00

3.4 g



CONTEC Icr is free of TEGDMA and HEMA.

CONTEC Ic / Icr primer

Light-curing.

For CONTEC Ic bracket adhesive (REF 163-111-00)

and CONTEC Icr retainer adhesive (REF 163-131-00).

REF 163-113-10

5 ml



CONTEC etching gel

35% phosphoric acid

For etching dental enamel.

Suitable for all CONTEC adhesives.

Without disposable tubes.

REF 163-115-00

2 pieces (2.5 ml each)



Adhesives and accessories



CONTEC go set

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

For approx. 600 brackets (30 treatment cases).

Content:

- 3 x CONTEC go bracket adhesive, each 2.5 g (REF 163-121-00)
- 1 x CONTEC go primer, 5 ml (REF 163-114-00)
- 2 x CONTEC etching gel, je 2,5 ml (REF 163-115-00)
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

CONTEC go is free of TEGDMA and HEMA.

Only compatible with REF 163-114-00 CONTEC go primer and

REF 163-121-00 CONTEC go bracket adhesive. No other

Dentaurum adhesive is compatible and may, therefore, not be

used in combination.

REF 163-107-00

1 set

NEW



CONTEC go bracket adhesive

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

CONTEC go is free of TEGDMA and HEMA.

Only compatible with REF 163-107-00 CONTEC go set and

REF 163-114-00 CONTEC go primer. No other Dentaurum

adhesive is compatible and may, therefore, not be used in

combination.

REF 163-121-00

2 pieces (each 2.5 g)

NEW



CONTEC go primer

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets

Only compatible with REF 163-107-00 CONTEC go set and

REF 163-121-00 CONTEC go bracket adhesive. No other

Dentaurum adhesive is compatible and may, therefore, not be

used in combination.

REF 163-114-00

5 ml

NEW



CONTEC etching gel

35% phosphoric acid

For etching dental enamel.

Suitable for all CONTEC adhesives.

Without disposable tubes.

REF 163-115-00

2 pieces (2.5 ml each)

Adhesive residue remover

For gentle removal of residual adhesive and cement without damaging the enamel.

Use water cooling!

Recommended speed: 5000 min⁻¹ - max. 10000 min⁻¹

Form	REF	Quantity
conical	123-700-00	10 pieces
radial	123-701-00	10 pieces
pear, large	123-702-00	10 pieces
pear, small	123-703-00	10 pieces

**Carbide bur****sterilizable**

For the removal of adhesive residue.

Shaft ø: 2.35 mm

Recommended speed: 40000 min⁻¹

Form	REF	Quantity	REF	Quantity
fine	123-604-00	1 piece	123-604-30	10 pieces
coarse	123-603-00	1 piece	123-603-30	10 pieces

