

Bracket accessories



Bite-Bumper™

Bonding technique
Material: Plastic
Palatal bite raiser for raising a deep overbite in Class II treatment.

REF 705-000-00	10 pieces
----------------	-----------



Titanium button with chain acc. to Watted

Bonding technique
With laser-structured titanium base and titanium chain.
For aligning retracted anteriors and cuspids.
Length: 5 cm

REF 750-003-51	1 piece
----------------	---------



Lingual buttons

Bonding technique
With laser-structured base, flat base, square. For incisors.
Dimensions (H x W): 3.5 mm x 3.5 mm

REF 750-002-51	10 pieces
----------------	-----------



Lingual buttons

Bonding technique
With laser-structured base, contoured base, round.
For cuspids and bicuspid.
ø 3.0 mm

REF 750-001-51	10 pieces
----------------	-----------



Lingual cleats

Bonding technique
With mesh base for bonding, flat.
For additional elastics or traction springs.
Dimensions (H x W): 3.8 mm x 3.8 mm

REF 750-602-50	10 pieces
----------------	-----------

Mesh base

Bonding technique
Anatomically contoured.



Indication	Base	Dimensions (H x W)	REF	Quantity
for incisors	flat	3.8 mm x 3.8 mm	799-202-00	10 pieces
for buccal tubes and 2nd molars	small	3.0 mm x 6.5 mm	799-206-00	10 pieces

Protection wax

Protects mucosa and lips in first days after bracket bonding until patient has become used to the appliance. Also for relief from irritation during treatment.

Form of delivery:

50 small boxes with 5 bars of protection wax each.

Color: white

REF 124-305-01

50 packs



Bracket tray

With self sticking foil to presort brackets.

REF 891-002-00

50 pieces



Bracket container

For storage of previously sorted brackets.

REF 891-001-00

1 piece



Storage box

For sorting brackets individually.

REF 766-003-10

1 piece

