



Quality Certificate for remanium® **GM 800+**



1 Indication

remanium® GM 800+ is a Cobalt Chromium alloy (CoCr alloy) for use in dental laboratories for manufacturing rigid frameworks with a very high yield strength.

2 Ingredients in %

Co 58.3, Cr 32.0, Mo 6.5, W 1.5, Si 1.0 / others < 1: (Mn, N, C)

3 Biocompatibility

Tests run by independent testing institutes, show excellent resistance to corrosion and no cytotoxicity.

4 Safety

remanium® alloys are certified by CE 0483 and are in compliance with the regulations of the Medical Device Directive EG 93/42.

5 Quality

In-house production ensures excellent quality due to continuous control during the manufacturing process (certified DIN EN ISO 22674, DIN EN ISO 13485, DIN EN ISO 9001).

Dr. Thomas Wiest
Technical Director Chemistry