



CAD/Vest und CAD/Press

Wachs fräsen leicht gemacht.

CAD/Vest and CAD/Press

Wax milling made easy.

Wachsblanks mit hervorragenden Fräseigenschaften zur Herstellung von Kronen und Brücken oder Modellguss im CAD/Vest Verfahren.

Verarbeitung in allen gängigen Fräsanlagen für Standard-Abmessungen von Fräsblanks.

StarWax Blank grey

ø 98,5 mm mit Bund

Ideal für CAD/Vest Kronen- und Brückenarbeiten zum Pressen oder Gießen:

- Rückstandsloses Ausbrennverhalten
- Hohe Kantenstabilität
- Leichte Elastizität



Farbe	Stärke	REF	Menge
grau	20 mm	120-230-00	1 Stück

StarWax Blank blue

ø 98,5 mm mit Bund

Ideal für CAD/Vest Modellgussarbeiten zum Gießen:

- Spannungsfrei
- Hohe Bruchfestigkeit
- Hohe Passgenauigkeit



Farbe	Stärke	REF	Menge
blau	25 mm	120-235-00	1 Stück

Wax blanks with excellent milling properties for the manufacture of crowns and bridges or cast partial dentures using the CAD/Vest process.

Can be processed on all milling machines suitable for milling blanks with standard dimensions.

StarWax Blank grey

ø 98.5 mm with collar

Ideal for CAD/Vest crowns and bridges to be pressed or cast:

- Residue-free burnout
- High edge strength
- Slightly elastic



Color	Thickness	REF	Quantity
grey	20 mm	120-230-00	1 piece

StarWax Blank blue

ø 98.5 mm with collar

Ideal for CAD/Vest cast partial dentures to be cast:

- Tension-free
- High fracture strength
- Accurate fit



Color	Thickness	REF	Quantity
blue	25 mm	120-235-00	1 piece