

Minimal effort – *enhanced aesthetics*



ceraMotion®
PZr

Photo: © Christian Ferrario

ceraMotion® PZr

Press-on ceramic for zirconia frameworks



D
DENTAURUM

ceraMotion® PZr




Press-on ceramic for zirconia frameworks

ceraMotion® PZr is a press-on ceramic for use with zirconia frameworks. It is compatible with the ceraMotion® Zr veneering ceramic.

- ▷ high cost-effectiveness due to the full anatomical reproduction of the modelled pattern
- ▷ simplified working steps for large span bridges
- ▷ functional design of the occlusal surfaces, shoulder region or contact points
- ▷ fully compatible with ceraMotion® Zr veneering ceramic

**The common properties shared by all product lines combine to produce the special product advantages of the ceraMotion® ceramic system:
The highest standards of aesthetics, reliability and cost-effectiveness**

The innovative Touch Up concept enhances the aesthetics with minimum effort. ceraMotion® press ceramics harmonise perfectly with the veneering materials. Maximum shade stability is attained due to the thermocoloration process. The special production process guarantees outstanding shape and edge stability as well as a constant CTE value. The ceraMotion® ceramic concept also demonstrates that a clear overview can be maintained with a very comprehensive product range. The unparalleled product variety of the ceraMotion® ceramic system is based on the unique spectrum of expertise of the Dentaaurum Group: development and manufacture of ceramic, alloys and investment materials from a single source coupled with an exemplary range of services.



Minimal effort –
enhanced aesthetics



ceraMotion® PZr product overview

Press Liner	PZr L	A-D
Ingot	PZr I	A-D

Fully compatible with all materials within the ceraMotion® Zr line.

Date of information: 01/14

Subject to modifications

ceraMotion® – Development and manufacture by:

