




Distalization appliance.




DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





advanced
molar
distalization
appliance.



advanced molar distalization appliance.

Together with the tomas® system, amda® offers a safe system for the distalization of upper molars. As a skeletal anchorage device it can be used for the bi- or unilateral distalization or the mesialization of upper molars and for the retraction of the upper anterior teeth depending on how it is fitted.

The appliance serves as a skeletally anchored transpalatal arch during the retention phase. amda® is more than merely a distalization appliance!

The amda® system was invented by Prof. Dr. med. dent. Moschos A. Papadopoulos. Dentaureum developed the system in close cooperation with him. Dr. Papadopoulos is a Professor and Chairman of the Department of Orthodontics at the School of Dentistry, Aristotle University in Thessaloniki, Greece.



Indications.

- Distalization of molars.
- Mesialization of molars.
- Retraction of anterior teeth.
- Combination of distalization and mesialization.
- Retention of molars.


Advantages:

- No patient cooperation required.
- In situ – 24 hours.
- Hygienic and comfortable for the patient.
- No contamination of the spring.
- Easy to fit/Easy to remove.



Distalization appliance

amda® / Starter set



amda® starter set (assembled)

Skeletally anchored appliance for the distalization of molars with the tomas®-pin EP.



REF 607-180-00	1 set
----------------	-------

CE 0483

Contents:

2 x amda® telescope (2 pieces)	REF 607-181-02
1 x amda® connector, pair 7 mm	REF 607-184-07
1 x amda® connector, pair 10 mm	REF 607-184-10
1 x amda® connector, pair 13 mm	REF 607-184-13
1 x amda® palatal arch	REF 607-185-01
1 x Hexagon socket key 0.9	REF 607-129-00

amda® telescope

Encapsulated stainless steel spring with two stop screws
(REF 302-013-01).
Spring force: max. 4,7 N

REF 607-181-02	2 pieces
REF 607-181-10	10 pieces

CE 0483



amda® connector, left

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-182-07	10 pieces
10 mm	607-182-10	10 pieces
13 mm	607-182-13	10 pieces

CE 0483



amda® connector, right

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-183-07	10 piece
10 mm	607-183-10	10 pieces
13 mm	607-183-13	10 pieces

CE 0483




amda® connector, pair

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-184-07	1 pair
10 mm	607-184-10	1 pair
13 mm	607-184-13	1 pair

CE 0483



Distalization appliance

amda® / Appliance and accessories



amda® palatal arch

(incl. 2 x tomas®-abutments EP, REF 302-027-00).

Palatal arch for amda® for fixing the amda® telescopes.

REF 607-185-01	1 piece
REF 607-185-10	10 pieces

CE 0483



Hexagon socket key 0.9

sterilizable

Key width: 0.9 mm

REF 607-129-00	1 piece
----------------	---------

CE



Orthodontic demonstration model tomas®/amda®

Fitted with:

- tomas®-pin EP
 - amda®
 - lingual / palatal sheaths "easy in", welded to standard molar bands
- Maxilla made of plastic.

REF 892-000-26	1 piece
----------------	---------

NEW



