

Premium-Line Pliers / Fixed Appliance Technique



004-350-00	1 piece
------------	---------

CE

Special Bracket Debonding Pliers, 45°, Angled

Sterilizable

For the removal of brackets and bonding adhesive. 45° angled cutting surface. Hard metal inserts. Sharpening possible by the manufacturer.

Premium-Line · 10 Year Guarantee



003-349-00	1 piece
------------	---------

CE

Bracket Removing Pliers

Sterilizable

For removing brackets between bracket base and adhesive. Sharpening possible by the manufacturer.

Replacement plastic insert see Accessories REF 038-356-00.

Premium-Line · 10 Year Guarantee



004-349-00	1 piece
------------	---------

CE

Bracket Removing Pliers

Sterilizable

For removing brackets and adhesive residue. Sharpening possible by the manufacturer. With rubber insert.

Premium-Line · 10 Year Guarantee



037-355-00	2 pieces
------------	----------

CE

Rubber Inserts

Sterilizable

For Bracket-Removing Plier REF 004-349-00.



026-347-00	1 piece
------------	---------

CE

Oliver-Jones Band Removing Pliers

Sterilizable

With replaceable plastic tip.

With replaceable hard metal tip. Sharpening possible by the manufacturer.

Plastic tip see Accessories REF 039-357-00.

Premium-Line · 10 Year Guarantee



003-114-00	1 piece
------------	---------



Johnson Contouring Pliers

Sterilizable

With curvature of the beaks only on one side. For contouring and adapting bands and crowns.

Premium-Line · 10 Year Guarantee

Premium-Line Pliers / Fixed Appliance Technique

Double Beak Pliers


Sterilizable

With concave beaks. For forming bicuspid and molar bands in the mouth. Can be used on both sides.

Premium-Line · 10 Year Guarantee

015-157-00	1 piece
------------	---------

CE



Separating Pliers

Sterilizable

For expanding Dentalastics® separating rings before positioning.

Premium-Line · 10 Year Guarantee

015-160-00	1 piece
------------	---------

CE



How Pliers, Straight

Sterilizable

Universal pliers for the multi-band technique. Due to the slender working tips perfect for placing bands with prewelded brackets, especially in the lower anterior area.

Premium-Line · 10 Year Guarantee

003-112-00	1 piece
------------	---------

CE



How Pliers, Curved

Sterilizable

These curved pliers are especially suited for the posterior areas. The curved lip-safe tips adapt especially well to the shape of the tooth.

Premium-Line · 10 Year Guarantee

003-111-00	1 piece
------------	---------

CE



Side Cutter, Maxi

Sterilizable

For wires up to \varnothing 1.0 mm / 40, hard. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

Premium-Line · 10 Year Guarantee

004-266-00	1 piece
------------	---------

CE



Distal End Cutter, Mini

Sterilizable

This distal end cutter has a cut and hold design.

Delicate style. Inserts made of hard metal. Sharpening of inserts possible by the manufacturer.

For wires up to \varnothing 0.56 x 0.70 mm / 22 x 28.

Premium-Line · 10 Year Guarantee

Length: 12 cm / 4 3/4"

003-701-00	1 piece
------------	---------

CE



Premium-Line Pliers / Fixed Appliance Technique



003-702-00	1 piece
------------	---------

CE

Distal End Cutter, Mini

Sterilizable

This small size distal end cutter has longer handles plus a cut and hold design. Inserts made of hard metal.

Sharpening of inserts possible by the manufacturer.

For wires up to \square 0.56 x 0.70 mm / 22 x 28.

Premium-Line · 10 Year Guarantee

Length: 14 cm / 5 1/2"



003-700-00	1 piece
------------	---------

CE

Distal End Cutter, Maxi

Sterilizable

This sturdy distal end cutter has a cut and hold design. Inserts made of hard metal.

Sharpening of inserts possible by the manufacturer.

For wires up to \square 0.56 x 0.70 mm / 22 x 28.

Premium-Line · 10 Year Guarantee

Length: 12 cm / 4 3/4"



014-151-00	1 piece
------------	---------

CE

Ligature Cutter, Mini

Sterilizable

Accuracy right to the tip for smooth cutting of soft ligature wires. With hard metal tips replaceable by the manufacturer.

Sharpening of inserts possible by the manufacturer.

For wires up to \varnothing 0.4 mm / 16, soft.

Premium-Line · 10 Year Guarantee



014-152-00	1 piece
------------	---------

CE

Ligature Cutter 45°

Sterilizable

45° angled cutting tips for easy cutting of ligatures, especially in the posterior areas. With hard metal tips replaceable by the manufacturer.

Sharpening of inserts possible by the manufacturer.

For wires up to \varnothing 0.4 mm / 16, soft.

Premium-Line · 10 Year Guarantee



016-153-00	1 piece
------------	---------

CE

Side Cutters, Medium

Sterilizable

Black anodized hard metal tips replaceable by the manufacturer. Smooth cuts guaranteed by accurate cutting edges.

Sharpening of inserts possible by the manufacturer.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee



003-130-00	1 piece
------------	---------

Light Wire Pliers

Sterilizable

For especially fine wire bending work in the light wire technique.

For wires up to \varnothing 0.5 mm / 20, spring hard.

Premium-Line · 10 Year Guarantee

Premium-Line Pliers / Fixed Appliance Technique

Begg Light Wire Pliers

Sterilizable

For especially fine wire bending work in the light wire technique. Additional cutting device. Inserts made of hard metal.

For wires up to \varnothing 0.5 mm / 20, spring hard.

Premium-Line · 10 Year Guarantee

004-130-00	1 piece
------------	---------



Angle/Tweed Ribbon Arch Pliers

Sterilizable

For applying torque to edgewise wires.

Parallel grip prevents crushing of wire.

Hard metal inserts.

For wires up to \square 0.56 x 0.70 mm / 22 x 28.

Premium-Line · 10 Year Guarantee

003-442-00	1 piece
------------	---------



Weingart Universal Pliers

Sterilizable

For positioning the arch and bending the arch wire ends. Accurately closing serrated tips. With hard metal inserts replaceable by the manufacturer.

For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee

003-120-00	1 piece
------------	---------

CE



Weingart Universal Pliers Mini

Sterilizable

Very slender design. For positioning the arch and bending the arch wire ends. Accurately closing serrated tips.

For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee

003-121-00	1 piece
------------	---------

CE



Nance Loop Bending Pliers

Sterilizable

With four-step beaks for producing loops of different sizes in round and rectangular wires.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

000-001-00	1 piece
------------	---------



Nance Loop Closing Pliers

Sterilizable

For bending of stop loops and as special pliers for Adams clasps.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

001-002-00	1 piece
------------	---------



Jarabak Light Wire Pliers

Sterilizable

For wires up to \varnothing 0.5 mm / 20, spring hard.

Premium-Line · 10 Year Guarantee

003-125-00	1 piece
------------	---------



Premium-Line Pliers / Fixed Appliance Technique



012-350-00

1 piece

Tweed Loop Forming Pliers

Sterilizable

For especially fine bending work.
For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee



012-351-00

1 piece

Tweed Loop Forming Pliers

Sterilizable

For especially fine bending work. With individually replaceable hard metal tip.

For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee



038-357-00

3 pieces

Hard Metal Tip

Sterilizable

For Tweed Loop Forming Plier
REF 012-351-00.



004-109-00

1 piece

De La Rosa Contouring Pliers

Sterilizable

With guiding grooves for forming round and square arches 16, 18, 22.

Premium-Line · 10 Year Guarantee



003-185-00

1 piece

Arch Forming Pliers

Sterilizable

For individual construction and activation of lingual/palatal bars. Using these pliers, the bent tips of the wire will fit into all Dentaurem lingual/palatal tubes \varnothing 2 x 0.9 mm / 36.

For wires up to \varnothing 0.9 mm / 36, hard.

Premium-Line · 10 Year Guarantee



003-139-00

1 piece

Angle Wire Bending Pliers, Standard

Sterilizable

Sturdy constructed plier with a cone and pyramid shaped beak.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee



004-139-00

1 piece

Angle Wire Bending Pliers

Sterilizable

Sturdy constructed plier with a cone and pyramid shaped beak. Inserts made of hard metal.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

Premium-Line Pliers / Fixed Appliance Technique

Loop Forming Pliers, Extra Mini

Sterilizable

Ideally suited for forming different loops and adjusting preshaped arches for the fixed appliance technique.

For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee

003-350-00	1 piece
------------	---------



Loop Forming Pliers, Mini

Sterilizable

For loop forming and bending of thin round and edgewise wires.

For wires up to \varnothing 0.5 mm / 20, hard.

Premium-Line · 10 Year Guarantee

003-351-00	1 piece
------------	---------



Ligature Forming Pliers

Sterilizable

For forming ligatures.

For wires up to \varnothing 0.3 mm / 12, soft.

Premium-Line · 10 Year Guarantee

015-158-00	1 piece
------------	---------



Aderer 3-Prong Pliers, Extra Mini

Sterilizable

For bending wires / arches for the fixed appliance technique. Specially suited for bending NiTi wires / arches such as rematitan® "LITE" and Tensic®.

For wires up to max. \varnothing 0.5 mm / 20, hard or \square 0,48 x 0,64 mm / 19 x 25.

Premium-Line · 10 Year Guarantee

003-205-00	1 piece
------------	---------

CE



Face Bow Bending / Cutting Pliers, Magnum

Sterilizable

Ideal for bending and cutting the outer bows of a face bow.

For wires up to \varnothing 1.8 mm / 71.

Premium-Line · 10 Year Guarantee

003-180-00	1 piece
------------	---------



Premium-Line Pliers / Removable Appliance Technique



003-053-00

1 piece

Loop Forming Pliers, Medium

Sterilizable

For bending small springs and loops.
For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee



003-233-00

1 piece

Loop Forming Pliers, Maxi

Sterilizable

These sturdy pliers have one round and one concave beak for heavier wires, such as space maintainers, etc.

For wires up to \varnothing 0.8 mm / 32, hard.

Premium-Line · 10 Year Guarantee



001-650-00

1 piece

Ruhland Loop Forming Pliers, Magnum

Sterilizable

To bend coffin springs, labial arches, loops, etc.

For wires up to \varnothing 1.1 mm / 44, hard.

Premium-Line · 10 Year Guarantee



003-515-00

1 piece

Adams Pliers, Medium

Sterilizable

With two smooth, rectangular beaks.
For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

Length: 12 cm / 4 3/4"



004-515-00

1 piece

Adams Pliers, Medium

Sterilizable

With two smooth, rectangular beaks. Inserts made of hard metal.

For wires up to \varnothing 0.7 mm / 28, spring hard.

Premium-Line · 10 Year Guarantee

Length: 12 cm / 4 3/4"



003-094-00

1 piece

Adams Pliers, Maxi

Sterilizable

With two smooth, rectangular beaks.
For wires up to \varnothing 0.8 mm / 32, hard.

Premium-Line · 10 Year Guarantee

Length: 14.5 cm / 5 3/4"



003-516-00

1 piece

Adams Pliers, Medium

Sterilizable

With one smooth, rectangular and one smooth, rounded beak.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

Length: 14 cm / 5 1/2"

Premium-Line Pliers / Removable Appliance Technique

Adams Pliers, Maxi

Sterilizable

With one smooth, rectangular and one smooth, rounded beak. Inserts made of hard metal.

For wires up to \varnothing 0.8 mm / 32, hard.

Premium-Line · 10 Year Guarantee

Length: 16 cm / 6 1/4"

004-095-00	1 piece
------------	---------



Wire Bending Pliers, Wide

Sterilizable

Serrated tips.

Suitable for activating finger springs.

For wires up to \varnothing 0.8 mm / 32, hard.

Premium-Line · 10 Year Guarantee

033-452-00	1 piece
------------	---------



Wire Bending Pliers, Maxi

Sterilizable

For bending wires up to \varnothing 0.9 mm / 36, hard.

For cutting wires up to \varnothing 0.7 mm / 28, hard.

With hard metal inserts replaceable by manufacturer.

Premium-Line · 10 Year Guarantee

013-522-00	1 piece
------------	---------



Universal Wire Bending Pliers

Sterilizable

With serrated tips and guiding grooves.

For wires up to \varnothing 0.8 mm / 32, hard.

Premium-Line · 10 Year Guarantee

003-290-00	1 piece
------------	---------



Pointed Beak Pliers

Sterilizable

Smooth faces with guiding grooves. Inserts made of hard metal.

For wires up to \varnothing 0.7 mm / 28, spring hard.

Premium-Line · 10 Year Guarantee

004-455-00	1 piece
------------	---------



Pointed Beak Pliers

Sterilizable

Serrated faces with guiding grooves.

Fine pointed beaks.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

035-455-00	1 piece
------------	---------



Young Loop Bending Pliers

Sterilizable

Universal pliers for different size loops.

Ideal for labial arches.

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

012-074-00	1 piece
------------	---------



Premium-Line Pliers / Removable Appliance Technique



003-201-00	1 piece
------------	---------

Aderer 3-Prong Pliers, Mini

Sterilizable

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee

Length: 12.5 cm / 5"



003-200-00	1 piece
------------	---------

Aderer 3-Prong Pliers, Medium

Sterilizable

For wires up to \varnothing 0.9 mm / 36, hard.

Premium-Line · 10 Year Guarantee

Length: 13 cm / 5 1/4"



003-202-00	1 piece
------------	---------

Aderer 3-Prong Pliers, Maxi

Sterilizable

For wires up to \varnothing 1.5 mm / 59, hard.

Premium-Line · 10 Year Guarantee

Length: 15.5 cm / 6 1/4"



012-520-00	1 piece
------------	---------

Arrow Clasp Former

Sterilizable

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee




012-521-00	1 piece
------------	---------

Arrow Clasp Bending Pliers

Sterilizable

For wires up to \varnothing 0.7 mm / 28, hard.

Premium-Line · 10 Year Guarantee



Pliers / Removable Appliance Technique


Adams Clasp Bending Pliers

Can be disinfected but not sterilized.

Special pliers for easy forming of Adams Clasps in one step.

For wires up to \varnothing 0.7 mm / 28, hard.

012-523-00	1 piece
------------	---------



Side Cutter, Maxi

Can be disinfected but not sterilized.

Made of extra hard vanadium steel, chrome plated. Sharpening of inserts possible by the manufacturer.

For wires up to \varnothing 1.0 mm / 40, hard.

037-266-00	1 piece
------------	---------



Heavy Gauge Side Cutter, Magnum

Can be disinfected but not sterilized.

Little force required for cutting through large size wires due to special design.

For wires up to \varnothing 1.3 mm / 51, hard.

044-177-00	1 piece
------------	---------



Instruments



075-601-00	1 piece
------------	---------

CE

Lip Retractor, small size

Can be disinfected but not sterilized.
For efficient work without assistance.
With flexible plastic connector.



075-600-00	1 piece
------------	---------

CE

Lip Retractor, standard size

Can be disinfected but not sterilized.
For efficient work without assistance.
With flexible plastic connector.



025-277-00	1 piece
------------	---------

CE

Clamping Tweezers

Sterilizable
For bracket placement. The blade at the end of the tweezers can be used for angulating brackets.
Delicate design.



025-276-00	1 piece
------------	---------

CE

Clamping Tweezers

Sterilizable
With long placement and angulation aid. For bracket placement. The blade at the end of the tweezers can be used for angulating brackets.



025-279-00	1 piece
------------	---------

CE

Clamping Tweezers

For posterior teeth.
Sterilizable
Specially suited for placing brackets and buccal tubes in the posterior region.



027-349-00	1 piece
------------	---------

CE

Scaler and Band Seater, double ended

Sterilizable
Sharp scaling tip on one end for removing adhesive residue.
Serrated band seater on the other end.



018-000-00	1 Set
------------	-------

CE

Separating Strip Placer

Sterilizable
Holding instrument for separating strips. For smoothing protruding edges of fillings and separating. Complete with 6 separating strips.

Separating Strips, coated on both sides

Sterilizable

For Separating Strip Placer REF 018-000-00.

Width: 4 mm

018-001-00

12 pieces



Separating Strips, coated on one side

Sterilizable

For Separating Strip Placer REF 018-000-00.

Width: 4 mm

018-002-00

12 pieces



Scissors, Straight

Sterilizable

For cutting thin metal bands.

007-024-00

1 piece



Scissors, Curved

Sterilizable

For cutting thin metal bands.

008-025-00

1 piece



Eby Band Seater

Sterilizable

With curved tip for posterior region.

Stainless Steel.

Circular snap device.

Seating pressure range 0.8 - 10 N
(approx. 80 - 1,000 g).

051-023-00

1 piece

CE



Hard Metal Tip, Curved

Sterilizable

For Eby Band Seater REF 051-023-00.

051-024-00

1 piece

CE



Hard Metal Tip, Straight

Sterilizable

For Eby Band Seater REF 051-023-00.

051-025-00

1 piece

CE



Instruments



026-355-00	1 piece
------------	---------

CE

Band Seating Instrument, Nylon

Sterilizable to a maximum of 180 °C/356 °F

With hard metal, triangular serrated tip.



003-006-00	1 piece
------------	---------

CE

Band Seater

Sterilizable

Large surface tin inlay.



026-354-00	1 piece
------------	---------

CE

Molar Band Seater

Sterilizable to a maximum of 180 °C/356 °F

With replaceable plastic head and hard metal triangular tip.



026-357-00	1 piece
------------	---------

CE

Plastic Head with triangular hard metal tip

Sterilizable

For Molar Band Seater REF 026-354-00.



024-300-00	1 piece
------------	---------

CE

Mershon Band Pusher

Sterilizable

With serrated rectangular tip and anatomically formed grip for good transfer of force when positioning bands.



030-390-00	1 piece
------------	---------

CE

Bracket Positioning Gauge

Sterilizable

Simplifies seating the brackets on the teeth with exactly the right distance between bracket slot and incisal edge.

With metal tips.



030-391-00	1 piece
------------	---------

CE

Bracket Positioning Gauge

Can be disinfected but not sterilized.

Simplifies seating the brackets on the teeth with exactly the right distance between bracket slot and incisal edge.

With pencils.

Safety Caps

Can be disinfected but not sterilized.

For Bracket Positioning Gauge
REF 030-391-00.

Color: white

030-392-00	4 pieces
------------	----------

CE



“Mathieu” Needle Holder, Medium

Sterilizable

Ideal for ligating. Convenient and practical design with serrated tips.

000-030-00	1 piece
------------	---------

CE



“Bachmann” Needle Holder with Finger Ring, Medium

Sterilizable

Eases ligating with the middle finger under constant pressure. Can also be used in ligating jaw fracture splints. Serrated tips.

000-035-00	1 piece
------------	---------

CE



Needle Holder, Mini

Sterilizable

With especially fine serrated tips. For Dentalastics® plastic ligatures.

Not for wire ligatures.

000-031-00	1 piece
------------	---------

CE



Needle Holder, Mini

Sterilizable

With extremely fine serrated tips. For Dentalastics® plastic ligatures.

Not for wire ligatures.

000-032-00	1 piece
------------	---------

CE



Mosquito Forceps, without hook

Sterilizable

For placing Dentalastics® plastic ligatures.

000-730-00	1 piece
------------	---------

CE



Mosquito Forceps, with hook

Sterilizable

For placing Dentalastics® plastic ligatures.

000-731-00	1 piece
------------	---------

CE



Ligator

Sterilizable

An instrument for fast ligating of the complete arch. Y-shaped tip simplifies the twisting of the ligatures.

For ligating ligature wires up to ø 0.3 mm / 12, soft.

023-280-00	1 piece
------------	---------

CE



Instruments



023-282-00	1 piece
------------	---------

CE

Spool

Sterilizable

For Ligator REF 023-280-00.



023-281-00	1 piece
------------	---------

CE

Hard Metal Tip

Sterilizable

For Ligator REF 023-280-00.



022-274-00	1 piece
------------	---------

CE

Ligating Instrument

Sterilizable

Ligature adjuster and tucker for long preformed ligatures.



021-273-00	1 piece
------------	---------

CE

Ligating Instrument, double ended

Sterilizable

Zum Umbiegen und Adaptieren von Drahtligaturen.



019-000-00	1 piece
------------	---------

CE

Distal End Bender (Langer), double ended

Sterilizable

Sturdy instrument for bending arch wire ends.
Handle \varnothing 8 mm



004-180-00	1 piece
------------	---------

Arch Turret without Torque

Sterilizable

For forming flat, plain arches in edgewise wires.

For edgewise wires:

- 0.41 mm / 16
- 0.43 mm / 17
- 0.46 mm / 18
- 0.51 mm / 20
- 0.56 mm / 22

Blue tinted.



004-186-00	1 piece
------------	---------

Arch Turret with Torque / 16

Sterilizable

For forming flat or angled (torque) rectangular arches:

- 0.41 x 0.41 mm / 16 x 16
- in 0°, 7°, 10°, 13° and 16° Torque

Red tinted.

Arch Turret with Torque / 16 x 16

Sterilizable

For forming flat or angled (torque) rectangular arches:

☑ 0.41 x 0.41 mm / 16 x 16

in 0°, 10°, 16° and 23° Torque

Specially designed for Ricketts® technique.

Green tinted.

004-187-00	1 piece
------------	---------



Arch Turret with Torque / 18

Sterilizable

For forming flat or angled (torque) rectangular arches:

☑ 0.46 x 0.64 mm / 18 x 25

in 0°, 7°, 10°, 13° and 16° Torque

Gold tinted.

004-182-00	1 piece
------------	---------



Arch Turret with Torque / 22

Sterilizable

For forming flat or angled (torque) rectangular arches:

☑ 0.56 x 0.70 mm / 22 x 28

in 0°, 7°, 10°, 13° and 16° Torque

Silver tinted.

004-184-00	1 piece
------------	---------



Tension and Compression Gauge

Can be disinfected but not sterilized.

Handy device for measuring intraoral tensile and compressive forces. Pressure side with guide slot and fork. Tension side with hook. Measuring scales on pressure and tension sides.

Measuring range 28 - 450 g

006-013-00	1 piece
------------	---------

CE



Measuring Gauge

Can be disinfected but not sterilized.

For measuring tensile and compressive forces.

Measuring range 25 - 250 g

040-711-00	1 piece
------------	---------

CE



Measuring range 50 - 500 g

040-712-00	1 piece
------------	---------

CE



Measuring range 100 - 1,000 g

040-713-00	1 piece
------------	---------

CE



Diagnostic



042-750-00	1 piece
------------	---------

CE

Dental Vernier (Beerendonk)

Sterilizable

With fine, long tips for intraoral measurements.



042-751-00	1 piece
------------	---------

CE

Dental Vernier (Muenchner Design)

Sterilizable

For intra and extraoral measurements, with extra graceful tips.



044-730-01	1 piece
------------	---------

CE

Orthodontic Ruler (Bernklau)

Sterilizable

Made of Stainless Steel. Perfect for all diagnostic measurements, both intraoral and on the model.



766-000-00	1 piece
------------	---------

Arch Gauge

For measuring labial arches, Adams clasps, palatal bars, etc.

White plastic



044-731-00	1 piece
------------	---------

CE

Ruler (Muenchner Design)

Sterilizable

Stainless Steel, with 1/2 mm and mm increments.

Length: 15 cm / 6"



003-366-00	1 piece
------------	---------

Measuring Grid (Schmuth)

Can be disinfected but not sterilized.

This screen is printed red with a 1 x 1 cm scale on one side and the other side in black with a 1 cm x 2 mm scale. When the red and black scales are in alignment, the parallax correction has been achieved for precise measurement of distances.



029-360-01	1 Set
------------	-------

Orthometer Kit (Korkhaus)

Can be disinfected but not sterilized.

Set with orthometer, mm disk, cm disk and bow divider.



029-361-01	1 piece
------------	---------

Bow Divider

Can be disinfected but not sterilized.

For Korkhaus Orthometer Kit REF 029-360-01.

Three Dimensional Bow Divider (Korkhaus)

Can be disinfected but not sterilized.

Allows the direction and extent of tooth movement to be checked by a single intraoral measurement.

028-353-00	1 piece
------------	---------

CE



Bow Divider

Can be disinfected but not sterilized.

Standard design.

To measure tooth width, bracket distance to the incisal edge, etc.

030-395-00	1 piece
------------	---------



Cephalometric Tracing Kit

Contents:

- 4 lead pencils with lead 0.5 mm
- 1 package (= 12 pieces) red refills
- 1 package (= 12 pieces) green refills
- 1 package (= 12 pieces) blue refills
- 1 protractor (ruler)

075-400-01	1 Set
------------	-------




Diagnostic Ruler-Protractor

For evaluation of radiographs.

Dimensions: 153 x 38 mm

Plastic transparent, clear

004-367-00	1 piece
------------	---------



Tracing-Protractor Template

Transparent template with protractor and tooth symbols. For evaluation of radiographs.

Dimensions: 207 x 132 mm

Plastic transparent, clear

004-368-00	1 piece
------------	---------



Transparent Template

Ricketts® method

For analyzing radiographs. With mm scale, protractor and cut out tooth symbols for tracing.

Dimensions: 197 x 74 mm

Plastic transparent, clear

004-369-00	1 piece
------------	---------



Diagnostic



075-100-00

100 pieces

Acetate Tracing Foil

For radiograph tracings. One side matte, other side shiny.
Tear-off pad.

Dimensions: 17.5 x 23.5 cm



075-101-00

100 pieces

Dimensions: 24 x 30 cm



075-102-00

50 pieces

Acetate Tracing Foil Cover

For tracing and filing of radiographs.
One side matte, other side shiny.

Dimensions: 18 x 24 cm

Corrosion Protection Oil for Pliers and Instruments

For pliers and instrument care.
To use in particular before sterilization.

055-201-00	50 ml
------------	-------



Plastic Inserts

Sterilizable

For Bracket-Removing Pliers REF 003-349-00 and 045-005-00.

038-356-00	4 pieces
------------	----------

CE



Plastic Tips

Sterilizable

For Oliver-Jones Band Removing Pliers REF 026-347-00 and 045-004-00.

039-357-00	10 pieces
------------	-----------

CE



Stand for Wires and Pliers, combination

Made of clear acrylic.

With 6 tubes for straight wires.
Do not use disinfectants containing ethanol, phenol or alcohol.

Height: Approx. 13 cm / 5 1/8"

Width: Approx. 18 cm / 7"

055-101-01	1 piece
------------	---------



Stand for Pliers

Made of clear acrylic.

Do not use disinfectants containing ethanol, phenol or alcohol.

Length: 13 cm / 5 1/8", for 7 pliers

055-000-00	1 piece
------------	---------



Length: 25 cm / 9 7/8", for 14 pliers

055-100-00	1 piece
------------	---------



Stand for Pliers

Made of white acrylic.

For 15 pliers.
Do not use disinfectants containing ethanol, phenol or alcohol.

Width: 30 cm / 12"

Depth: 20 cm / 8"

Height: 6.5 cm / 2 1/2"

055-105-00	1 piece
------------	---------



Accessories



089-430-00

1 piece

Bur Stand Magnetic Holder

Partly plastic.

Do not use disinfectants containing ethanol, phenol or alcohol.



089-431-00

1 piece

Bur Stand

Made of clear acrylic.

Do not use disinfectants containing ethanol, phenol or alcohol.



043-700-00

6 pieces

Marking Pencils

White Wax Marker



043-800-00

1 piece

Marking Pencils

Wipe resistant

Black Marker



043-900-00

1 piece

Red Marker

