





Frequency of drill application

Batch / LOT

Product image	Description	ø	REF	Frequency of application ¹																																							
	Marking drill	ø 1.5 mm	382-015-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
	Stepped drill for ø 1.8 mm	ø 1.5 mm	302-118-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
	Stepped drill for ø 2.2 mm	ø 1.9 mm	302-122-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
	Stepped drill for ø 2.5 mm	ø 2.1 mm	302-125-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40

¹ **Note:** With careful and correct maintenance, the drills and thread taps can be used for 30 – 40 applications thanks to the high-quality materials and workmanship (see also Surgery + Prosthetics Manual 989-360-20).
Further information is available at www.dentaurum.com.