

# Overview of individual components.


ceraMotion®  
One Touch


Included in the ceraMotion® Set: ● One Touch (REF 250-800-70)

● No Limits (REF 250-801-70)


● Pink (REF 250-803-70)


## Glaze pastes


	<b>Paste Glaze bright</b> bright, fluorescent glaze	250-272-03	●
---	--	------------	---


	<b>Paste Glaze transpa</b> fluorescent glaze	250-271-03	●
---	---	------------	---


## 2D pastes


	<b>Paste 2D Body A</b> Body paste for Vita classical® A shades	250-306-03	●
---	---	------------	---


	<b>Paste 2D Body B</b> Body paste for Vita classical® B shades	250-307-03	●
---	---	------------	---


	<b>Paste 2D Body C</b> Body paste for Vita classical® C shades	250-308-03	●
---	---	------------	---


	<b>Paste 2D Body D</b> Body paste for Vita classical® D shades	250-309-03	●
---	---	------------	---


	<b>Paste 2D Body R</b> Body paste to shift the shade toward red (Vita 3D-Master®)	250-317-03	●
---	--	------------	---


	<b>Paste 2D Body L</b> Body paste to shift the shade toward yellowish (Vita 3D-Master®)	250-318-03	●
---	--	------------	---


	<b>Paste 2D Body R+</b> Body paste to shift the shade toward reddish (dark) (Vita 3D-Master®)	250-319-03	●
---	--	------------	---


	<b>Paste 2D Body orange</b> Body paste to shift the shade toward yellowish (dark) (Vita 3D-Master®)	250-320-03	●
---	--	------------	---


	<b>Paste 2D Incisal I2</b> Incisal for medium shades	250-315-03	●
---	---	------------	---


	<b>Paste 2D Incisal I3</b> Incisal for dark shades	250-316-03	●
---	---	------------	---

	<b>Paste 2D honey</b> Characterization: orange	250-312-03	●
---	---	------------	---


	<b>Paste 2D white</b> Characterization: white	250-310-03	●
---	--	------------	---


	<b>Paste 2D grey</b> Characterization: grey	250-311-03	●
---	--	------------	---

	<b>Paste 2D blue</b> Characterization: blue	250-313-03	●
---	--	------------	---

	<b>Paste 2D violet</b> Characterization: violet	250-314-03	●
---	--	------------	---

## Liquids

	<b>Diluting Liquid</b> Adjusting the viscosity of the pastes	254-012-01	● ● ●
---	---	------------	-------

	<b>Refreshing Liquid</b> Adjusting the viscosity and maintaining the stability of 3D pastes	254-014-01	● ● ●
---	--	------------	-------

\*VITA classical® and VITA 3D-Master® are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co., Bad Säckingen, Germany

# Overview of individual components.














ceraMotion®  
One Touch

Included in the ceraMotion® Set: ● One Touch (REF 250-800-70)








● No Limits (REF 250-801-70)

● Pink (REF 250-803-70)




## 3D pastes

	Paste 3D Incisal opal bright very bright incisal (bleach)	250-298-03	●
	Paste 3D lumin bright incisal for bright shades	250-301-03	●
	Paste 3D I2 medium incisal for medium shades	250-300-03	●
	Paste 3D I3 dark incisal for dark shades	250-299-03	●
	Paste 3D transpa Increasing translucency	250-297-03	●
	Paste 3D opal grey Modifier – grey opalescent	250-304-03	●
	Paste 3D opal honey Modifier – orange opalescent	250-305-03	●
	Paste 3D opal blue Modifier – blue opalescent	250-303-03	●
	Paste 3D neutral reduces the intensity of opal pastes	250-302-03	●
	Paste 3D Dentin bleach for bleach shades	250-293-03	●
	Paste 3D Dentin light (bright) for bright shades	250-294-03	●
	Paste 3D Dentin medium for medium shades	250-295-03	●
	Paste 3D Dentin dark for dark shades	250-296-03	●

## Gingival pastes

	Paste 3D Gingival 1 bright pink	250-330-03	●
	Paste 3D Gingival 2 bright pink/orange	250-331-03	●
	Paste 3D Gingival 3 medium pink	250-332-03	●
	Paste 3D Gingival 4 dark pink/red	250-333-03	●
	Paste 3D Gingival white Individualizing 3D Gingival 1-4/Modifier	250-334-03	●
	Paste 3D Gingival transpa Individualizing – Increasing translucency	250-335-03	●
	Paste 3D Gingival Modifier dark red, slightly opaque	250-336-03	●

## Stains

	Stains ST white Powder stains for strong effects	254-201-02	●
	Stains ST red brown Powder stains for strong effects	254-212-02	●
	Stains Universal Set (16 stains) Compatible with all 2D and 3D pastes	254-800-30	