



my first
ceramic
in paste



ceraMotion®
One Touch



ceraMotion® One Touch is the first veneering ceramic in paste form. The 2D and 3D pastes have been specially developed for the quick and aesthetic finish of all-ceramic restorations made of lithium disilicate and zirconium oxide.



smart.
fast.
aesthetic.

Smart. What you see is what you get.
Fast. All in one firing.
Aesthetic. Natural translucency, opalescence and fluorescence.

The complete ceraMotion® One Touch system.



Paste 2D

Paste 3D



The **translucency**
of the 2D pastes maintains the light dynamics.

Paste 2D Incisal

Paste Glaze

Colorless glazes in transpa and bright.


Paste 3D Incisal

Incisal pastes for layering with natural in-depth effect.



Paste 3D Dentin

Thanks to the combination of larger ceramic particles with a thixotropic paste it is possible to create an individual **morphology** in the dentin and incisal area.



Paste 2D Modifier

The particle size of the 2D pastes is ideally suited to deliver a **homogeneous** and natural finish to the surface.



Paste 2D Body

Simple and precise adaptation of the desired tooth shade (VITA classical®* A-D, VITA 3D-Master®** and individual shades).



Paste 3D Opal

Natural **opalescence**.
Shade-coordinated 3D paste with an opalescent effect for individualization.



Paste 3D Gingival

For a perfect copy of **nature**.



All ceraMotion® One Touch pastes are fluorescent and have been adapted to reflect the **fluorescence** of natural teeth.



"When the fine tip of our brush comes in contact with the pastes, it creates a **symphony of forms and fissures**. It also produces the subtle effects".
Christian Ferrari®

* VITA classical® and VITA 3D-Master® are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co., Bad Säckingen, Germany

Dentaurum

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE



Our thanks go to Christian Ferrari for kindly providing the photographs used in this brochure.
Realisation and Photo: © Christian Ferrari

Date of information: 2024-01
Subject to modifications



Turnstr. 31 | 75228 Ispringen | Germany | Phone +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaurum.com | info@dentaurum.com