

# Chroma Concept

Color chart to reproduce the VITA System 3D Master\* colors.

| Color VITA System 3D-Master* | Dentin mixing table ceraMotion® Me Ti Zr |                 |     | Opaker or Liner | Incisal     |
|------------------------------|--|-----------------|-----|-----------------|-------------|
| M1                           | 100%                                     | CC D 1          |     | CC1 bleach      | CC11 bleach |
| M2                           | 50%                                      | CC D 1 + CC D 2 | 50% | CC1 bleach      | CC11 bleach |
| M3                           | 100%                                     | CC D 2          |     | CC1 bleach      | CC11 bleach |
| 1M1                          | 50%                                      | D A1 + CC D 2   | 50% | B1              | I1          |
| 1M2                          | 60%                                      | CC D 1 + CC D 5 | 40% | B2              | I1          |
| 2M1                          | 75%                                      | CC D 3 + CC D 1 | 25% | B2              | I2          |
| 2L1.5                        | 67%                                      | D A2 + D B1     | 33% | B2              | I2          |
| 2R1.5                        | 50%                                      | CC D 2 + CC D 6 | 50% | D2              | I2          |
| 2M2                          | 90%                                      | D A2 + D B3     | 10% | B2              | I2          |
| 2L2.5                        | 80%                                      | D B3 + CC D 1   | 20% | B2              | I2          |
| 2R2.5                        | 55%                                      | CC D 2 + CC D 9 | 45% | B2              | I2          |
| 2M3                          | 75%                                      | CC D 5 + CC D 2 | 25% | B2              | I2          |
| 3M1                          | 90%                                      | CC D 3 + CC D 5 | 10% | CC2             | I3          |
| 3L1.5                        | 50%                                      | CC D 3 + CC D 4 | 50% | CC2             | I3          |
| 3R1.5                        | 50%                                      | CC D 3 + CC D 6 | 50% | CC2             | I3          |
| 3M2                          | 50%                                      | CC D 3 + CC D 5 | 50% | CC2             | I3          |

| Color VITA System 3D-Master* | Dentin mixing table ceraMotion® Me Ti Zr |                   |     | Opaker or Liner | Incisal |
|------------------------------|--|-------------------|-----|-----------------|---------|
| 3L2.5                        | 50%                                      | CC D 4 + CC D 5   | 50% | CC2             | I3      |
| 3R2.5                        | 50%                                      | CC D 5 + CC D 6   | 50% | CC2             | I3      |
| 3M3                          | 90%                                      | CC D 5 + CC D 3   | 10% | CC2             | I3      |
| 4M1                          | 50%                                      | CC D 7 + CC D 10  | 50% | CC3             | I3      |
| 4L1.5                        | 75%                                      | CC D 7 + CC D 9   | 25% | CC3             | I3      |
| 4R1.5                        | 75%                                      | CC D 10 + CC D 8  | 25% | CC3             | I3      |
| 4M2                          | 50%                                      | CC D 7 + CC D 9   | 50% | CC3             | I3      |
| 4L2.5                        | 75%                                      | CC D 8 + CC D 10  | 25% | CC3             | I3      |
| 4R2.5                        | 75%                                      | CC D 9 + CC D 7   | 25% | CC3             | I3      |
| 4M3                          | 50%                                      | CC D 8 + CC D 9   | 50% | CC3             | I3      |
| 5M1                          | 50%                                      | CC D 10 + CC D 11 | 50% | CC4             | I3      |
| 5M2                          | 50%                                      | CC D 10 + CC D 12 | 50% | CC4             | I3      |
| 5M3                          | 50%                                      | CC D 9 + CC D 12  | 50% | CC4             | I3      |

\* VITA System 3D-Master is a registered brand of Vita Zahnfabrik H. Rauter GmbH & Co.KG, 79704 Bad Säckingen, Germany