

ceraMotion® Z Blanks

Kunden-Nummer

| | | | | | |
|--|--|--|--|--|--|
| | | | | | |
|--|--|--|--|--|--|



The symbol for smarter zirconia

ceraMotion® Z Hybrid

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge |
|----------|-------------|-------|-------------|-------|-------------|-------|
| A1 | 272-411-14 | | 272-411-18 | | 272-411-22 | |
| A2 | 272-412-14 | | 272-412-18 | | 272-412-22 | |
| A3 | 272-413-14 | | 272-413-18 | | 272-413-22 | |
| A3,5 | 272-415-14 | | 272-415-18 | | 272-415-22 | |
| A4 | 272-414-14 | | 272-414-18 | | 272-414-22 | |
| B1 | 272-421-14 | | 272-421-18 | | 272-421-22 | |
| B2 | 272-422-14 | | 272-422-18 | | 272-422-22 | |
| C2 | 272-432-14 | | 272-432-18 | | 272-432-22 | |
| Bleach 1 | 272-300-14 | | 272-300-18 | | 272-300-22 | |
| Bleach 2 | 272-400-14 | | 272-400-18 | | 272-400-22 | |



A1 A2 A3 A3.5 A4 B1 B2 C2 D2

- Material
- Durchmesser
- Biegefestigkeit (Biaxial)

ZrO₂ 3Y-TZP → 4Y-PSZ
 98,3 mm
 1300 MPa → 1.020 MPa



ceraMotion® Z Cubic Multishade

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge |
|-------|-------------|-------|-------------|-------|-------------|-------|
| A1 | 272-311-14 | | 272-311-18 | | 272-311-22 | |
| A2 | 272-312-14 | | 272-312-18 | | 272-312-22 | |
| A3 | 272-313-14 | | 272-313-18 | | 272-313-22 | |
| A3,5 | 272-315-14 | | 272-315-18 | | 272-315-22 | |
| A4 | 272-314-14 | | 272-314-18 | | 272-314-22 | |
| B1 | 272-321-14 | | 272-321-18 | | 272-321-22 | |
| B2 | 272-322-14 | | 272-322-18 | | 272-322-22 | |
| C2 | 272-332-14 | | 272-332-18 | | 272-332-22 | |
| D2 | 272-333-14 | | 272-333-18 | | 272-333-22 | |



A1 A2 A3 A3.5 A4 B1 B2 C2 D2

- Material
- Durchmesser
- Biegefestigkeit (Biaxial)

ZrO₂ 5Y-PSZ
 98,3 mm
 600 MPa



ceraMotion® Z HT Multishade

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge | REF (25 mm) | Menge |
|-------|-------------|-------|-------------|-------|-------------|-------|-------------|-------|
| A1 | 272-211-14 | | 272-211-18 | | 272-211-22 | | 272-211-25 | |
| A2 | 272-212-14 | | 272-212-18 | | 272-212-22 | | 272-212-25 | |
| A3 | 272-213-14 | | 272-213-18 | | 272-213-22 | | 272-213-25 | |
| A3,5 | 272-215-14 | | 272-215-18 | | 272-215-22 | | 272-215-25 | |
| A4 | 272-214-14 | | 272-214-18 | | 272-214-22 | | 272-214-25 | |
| B1 | 272-221-14 | | 272-221-18 | | 272-221-22 | | 272-221-25 | |
| B2 | 272-222-14 | | 272-222-18 | | 272-222-22 | | 272-222-25 | |
| C2 | 272-232-14 | | 272-232-18 | | 272-232-22 | | 272-232-25 | |
| D2 | 272-233-14 | | 272-233-18 | | 272-233-22 | | 272-233-25 | |



A1 A2 A3 A3.5 A4 B1 B2 C2 D2

- Material
- Durchmesser
- Biegefestigkeit (Biaxial)

ZrO₂ 3Y-TZP
 98,3 mm
 1.100 MPa



ceraMotion® Z Blanks

Kunden-Nummer

| | | | | | |
|--|--|--|--|--|--|
| | | | | | |
|--|--|--|--|--|--|

The symbol for smarter zirconia

ceraMotion® Z HT Shade

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge | REF (25 mm) | Menge |
|-------|-------------|-------|-------------|-------|-------------|-------|-------------|-------|
| A1 | 272-111-14 | | 272-111-18 | | 272-111-22 | | 272-111-25 | |
| A2 | 272-112-14 | | 272-112-18 | | 272-112-22 | | 272-112-25 | |
| A3 | 272-113-14 | | 272-113-18 | | 272-113-22 | | 272-113-25 | |
| A3,5 | 272-115-14 | | 272-115-18 | | 272-115-22 | | 272-115-25 | |
| A4 | 272-114-14 | | 272-114-18 | | 272-114-22 | | 272-114-25 | |
| B1 | 272-121-14 | | 272-121-18 | | 272-121-22 | | 272-121-25 | |
| B2 | 272-122-14 | | 272-122-18 | | 272-122-22 | | 272-122-25 | |
| B3 | 272-123-14 | | 272-123-18 | | 272-123-22 | | 272-123-25 | |
| B4 | 272-124-14 | | 272-124-18 | | 272-124-22 | | 272-124-25 | |
| C1 | 272-131-14 | | 272-131-18 | | 272-131-22 | | 272-131-25 | |
| C2 | 272-132-14 | | 272-132-18 | | 272-132-22 | | 272-132-25 | |
| C3 | 272-133-14 | | 272-133-18 | | 272-133-22 | | 272-133-25 | |
| C4 | 272-134-14 | | 272-134-18 | | 272-134-22 | | 272-134-25 | |
| D2 | 272-142-14 | | 272-142-18 | | 272-142-22 | | 272-142-25 | |
| D3 | 272-143-14 | | 272-143-18 | | 272-143-22 | | 272-143-25 | |
| D4 | 272-144-14 | | 272-144-18 | | 272-144-22 | | 272-144-25 | |



- Material
- Durchmesser
- Biegefestigkeit (Biaxial)

ZrO₂ 3Y-TZP
98,3 mm
1.100 MPa



ceraMotion® Z HT White

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge | REF (25 mm) | Menge |
|-------|-------------|-------|-------------|-------|-------------|-------|-------------|-------|
| White | 272-094-14 | | 272-094-18 | | 272-094-22 | | 272-094-25 | |

- Material
 - Durchmesser
 - Biegefestigkeit (Biaxial)
- ZrO₂ 3Y-TZP
98,3 mm
1.200 MPa



ceraMotion® Z White

| Farbe | REF (14 mm) | Menge | REF (18 mm) | Menge | REF (22 mm) | Menge | REF (25 mm) | Menge |
|-------|-------------|-------|-------------|-------|-------------|-------|-------------|-------|
| White | 272-098-14 | | 272-098-18 | | 272-098-22 | | 272-098-25 | |

- Material
 - Durchmesser
 - Biegefestigkeit (Biaxial)
- ZrO₂ 3Y-TZP-A
98,3 mm
1.200 MPa

