

IT'S EASY
#CLICK&
GO. SO
FAST &
RELIABLE

discovery[®]
SMART sl



PROGRAMMA DI FORNITURA
discovery[®] smart sl

GENTLE, STABLE AND EFFICIENT.

#CLINICAL & GO



discovery[®]
SMART sl

Assortimenti

Bracket e tubi:

discovery[®]
SMART sl

(15-11/21-25, 45-41/31-35)

ORTHOCAST[®] SL

(16/26, 46/36)

ORTHOCAST[®]
M-SERIES MINI

(17/27, 47/37)

McLaughlin-Bennett-Trevisi* 22

Slot: 0,56 x 0,76 mm / 22 x 28



Torque	- 14°	- 14°	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°	- 14°	- 14°
Angolazione			0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°		
In/Out			1,0 mm	0,7 mm	0,7 mm	1,1 mm	0,7 mm	0,7 mm	1,1 mm	0,7 mm	0,7 mm	1,0 mm		
Larghezza	2,5 mm	3,6 mm	2,8 mm	2,8 mm	2,8 mm	2,6 mm	2,8 mm	2,8 mm	2,6 mm	2,8 mm	2,8 mm	2,8 mm	3,6 mm	2,5 mm
Distal Offset	0°	10°											10°	0°
REF senza gancio			722-002-15	722-002-14		722-002-12	722-002-11	722-002-21	722-002-22		722-002-24	722-002-25		
REF con gancio	822-412-17	822-912-16	722-042-15	722-042-14	722-042-13					722-042-23	722-042-24	722-042-25	822-912-26	822-412-27
Quantità	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi

SUP dx.

7	6	5	4	3	2	1	1	2	3	4	5	6	7
7	6	5	4	3	2	1	1	2	3	4	5	6	7

SUP sx.

INF dx

INF sx

Torque	- 10°	- 20°	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°	- 20°	- 10°
Angolazione			+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°		
In/Out			0,7 mm	0,7 mm	0,7 mm	1,1 mm	1,1 mm	1,1 mm	1,1 mm	0,7 mm	0,7 mm	0,7 mm		
Larghezza	2,5 mm	3,6 mm	2,8 mm	2,8 mm	2,8 mm	2,6 mm	2,6 mm	2,6 mm	2,6 mm	2,8 mm	2,8 mm	2,8 mm	3,6 mm	2,5 mm
Distal Offset	0°	0°											0°	0°
REF senza gancio			722-002-45	722-002-44		722-002-42	722-002-41	722-002-31	722-002-32		722-002-34	722-002-35		
REF con gancio	822-412-47	822-912-46	722-042-45	722-042-44	722-042-43					722-042-33	722-042-34	722-042-35	822-912-36	822-412-37
Menge	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi	10 pezzi

Assortimenti:

1 caso	SUP/INF 5-5	gancio sui 3	20 barcket	in blister	REF 722-002-01
1 caso	SUP/INF 5-5	gancio sui 3, 4, 5	20 barcket	in blister	REF 722-042-01
1 caso	SUP/INF 7-7	gancio sui 3	20 barcket + 8 tubi	in blister	REF 722-002-07
1 caso	SUP/INF 7-7	gancio sui 3, 4, 5	20 barcket + 8 tubi	in blister	REF 722-042-07
10 casi	SUP/INF 5-5	gancio sui 3	200 barcket	in contenitore	REF 722-002-10
10 casi	SUP/INF 5-5	gancio sui 3, 4, 5	200 barcket	in contenitore	REF 722-042-10
10 casi	SUP/INF 6-6	gancio sui 3	200 barcket + 40 tubi	in contenitore	REF 722-002-17
10 casi	SUP/INF 6-6	gancio sui 3, 4, 5	200 barcket + 40 tubi	in contenitore	REF 722-042-17

Accessori:

Strumento sl per discovery [®] smart sl	REF 023-274-00
Modello demo discovery [®] smart sl	REF 892-000-36



Strumento di apertura

* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin-Bennett and Trevisi.