

BENDER II



Retainer production
made easy



Retainer production made easy!

The Bender II is the second generation of a computer-controlled wire bending machine from the company YOAT that was especially developed for the production of individual lingual retainers. The system, unique worldwide, consists of a compact desktop unit and an easy-to-follow software (FixR)*. Complete your digital workflow and increase your profit by saving time and costs!

BENDER II

Lingual fixed retainers in 5 minutes

Highlights at a glance:

01.

Service and support

Use our "one-stop" principle: receive full service and support for the Bender II, FixR* and our Dentaurem wires – all from one single point of contact.

02.

Efficient work routine

With the Bender II, it is no longer necessary to create and send impressions or master models. Production and fitting of the retainer can therefore take place in only one sitting instead of two.

03.

Time-saving production

It is possible to produce a retainer in only 5 minutes using the Bender II. It can be produced in real time, even during the debonding process.

*We can send you the system requirements on request.

** Bender machine outer dimensions (w/d/h): 240 mm / 400 mm / 285 mm Weight: 15 kg

*** One wire in the specified quality can be processed per machine.

Advantages:

◆ Compact design

It is designed to save space** and can be set up even where space is tight!

◆ Ergonomic design

The finished retainer can be easily removed from the machine thanks to an ergonomically-formed output tray. No instruments, such as tweezers, are required.

◆ User-friendly software

Qualified personnel with little previous knowledge can design retainers thanks to the intuitive software (FixR)*.



Request you individual offer now.

For design ***, the offer includes:

1 software (FixR)*, installation and instruction

Selection Design	REF Bender II	Wire	Form	REF Wire
<input type="checkbox"/>	450-000-00	for remanium® straight wire round 0.45 mm / 18, extra spring hard	●	535-045-00
<input type="checkbox"/>	450-001-00	for remanium® straight wire, round 0.50 mm / 20, spring hard	●	535-002-00
<input type="checkbox"/>	450-002-00	for dentaflex® straight wire, 3-strand-twisted, round 0.45 mm / 18, super spring hard	●●●	545-745-00
<input type="checkbox"/>	450-003-00	for dentaflex® straight wire, 3-strand-twisted, round 0.50 mm / 20, super spring hard	●●●	545-750-00
<input type="checkbox"/>	450-004-00	for remanium® straight wire, rectangular 0.41 mm x 0.41 mm / 16 x 16, spring hard	■	537-340-00
<input type="checkbox"/>	450-005-00	for remanium® straight wire, rectangular 0.41 mm x 0.56 mm / 16 x 22, spring hard	■	537-310-00

Please get in touch with your local contact for individual offers.

Customer:

Customer number

--	--	--	--	--	--

Street

ZIP code

City

Name

Date

Signature customer



CONTACT
DENTAURUM

Dentaurum GmbH & Co. KG



Turnstr. 31
75228 Ispringen · Germany



info@dentaurum.com
www.dentaurum.com



+49 72 31 / 803 - 0



ONLINE SHOP
SHOP.DENTAURUM.COM